

Choice Framework for local Policy and Procedures 01-01 – Management and decontamination of surgical instruments (medical devices) used in acute care. Part D: Washer-disinfectors



Choice Framework for local Policy and Procedures (CFPP) 01-01

Management and decontamination of surgical instruments (medical devices) used in acute care

Part D: Washer-disinfectors

© Crown copyright 2013

Terms of use for this guidance can be found at <http://www.nationalarchives.gov.uk/doc/open-government-licence/>

Preface

Introduction

The Choice Framework for local Policy and Procedures (CFPP) is an initiative being piloted by the Department of Health.

It forms a suite of evidence-based policy and guidance documents on the management and decontamination of reusable medical devices.

Purpose

The purpose of CFPP is to enable local choices to be made regarding the management, use and decontamination of reusable medical devices at controlled costs using risk control.

CFPP is designed to reflect the need to continuously improve outcomes in terms of:

- patient safety;
- clinical effectiveness; and
- patient experience.

Essential Quality Requirements and Best Practice

The Health Act Code of Practice recommends that healthcare organisations comply with guidance establishing Essential Quality Requirements and demonstrate that a plan is in place for progression to Best Practice.

Essential Quality Requirements (EQR), for the purposes of this best practice guidance, is a term that encompasses all existing statutory and regulatory requirements. EQRs incorporate requirements of the current Medical Devices

Directive and Approved Codes of Practice as well as relevant applicable Standards. They will help to demonstrate that an acute provider operates safely with respect to its decontamination services.

Local policy should define how a provider achieves risk control and what plan is in place to work towards Best Practice.

Best Practice is additional to EQR. Best Practice as defined in this guidance covers non-mandatory policies and procedures that aim to further minimise risks to patients; deliver better patient outcomes; promote and encourage innovation and choice; and achieve cost efficiencies.

Best Practice should be considered when developing local policies and procedures based on the risk of surgical procedures and available evidence. Best Practice encompasses guidance on the whole of the decontamination cycle, including, for example, improved instrument management, where there is evidence that these procedures will contribute to improved clinical outcomes.

The CFPP suite is listed below.

- Choice Framework for local Policy and Procedures 01-01: Management and decontamination of surgical instruments (medical devices) used in acute care
- Choice Framework for local Policy and Procedures 01-04: Decontamination of linen for health and social care
- Choice Framework for local Policy and Procedures 01-06: Decontamination of flexible endoscopes

Acknowledgements

Andrew Thomas Royal College of Surgeons

Bill Keevil Southampton University

Brian Kirk IHEEM Decontamination Technology Platform

David Perrett Barts & the London School of Medicine & Dentistry, Queen Mary University of London

Gavin Maxwell Royal Society of Medicine Patients Support Group

Geoff Sjogren Institute of Decontamination Sciences

Geoffrey L. Ridgway, OBE, MD Clinical Microbiologist

Graham Stanton NHS Wales Shared Services Partnership – Facilities Services

Helen Baxter University of Edinburgh

Jackie Duggan Health Protection Agency

James Ironside University of Edinburgh

Jim Gray Birmingham Children's Hospital NHS Foundation Trust

Jimmy Walker Health Protection Agency

Ken Toal Health Estates, Northern Ireland

Margaret Hollis Great Ormond Street Hospital

Mike Painter Public Health Physician

Mike Simmons Public Health Wales

Peter Bringham Newcastle upon Tyne Hospitals NHS Foundation Trust

Peter Hoffman Health Protection Agency

Robert Baxter University of Edinburgh

Robert Kingston IHEEM Decontamination Technology Platform

Sylvia Martin University College Hospital London

Terry Durack Great Ormond Street Hospital

Tracy Coates Association for Perioperative Practice

Abbreviations

ACDP-TSE RM [subgroup]: Advisory Committee on Dangerous Pathogens – Transmissible Spongiform Encephalopathies Risk Management [subgroup] (formerly the TSE Working Group)

ACDST: Advisory Committee on Decontamination Science and Technology

AE(D): Authorising Engineer (Decontamination)

AP(D): Authorised Person (Decontamination)

BCH: Birmingham Children’s Hospital

BS: British Standard

BSE: Bovine Spongiform Encephalopathy

CFPP: Choice Framework for local Policy and Procedures

CJD: Creutzfeldt-Jakob disease

CMO: Chief Medical Officer

CP(D): Competent Person (Decontamination)

CQC: Care Quality Commission

DH: Department of Health

DIPC: Director of Infection Prevention and Control

EDIC: episcopic differential interference contrast

EDIC/EF: episcopic differential interference contrast/epifluorescence

EFSCAN: epifluorescent surface scanner

EN: European norm

FITC: fluorescein isothiocyanate

GOSH: Great Ormond Street Hospital

ISO: International Standards Organisation

MDD: Medical Devices Directive

MDR: Medical Devices Regulations

MHRA: Medicines and Healthcare products Regulatory Agency

NDS: National Decontamination Survey

NICE: National Institute for Health and Clinical Excellence

NICE IPG 196 (2006): NICE’s (2006) interventional procedure guidance 196 – ‘Patient safety and reduction of risk of transmission of Creutzfeldt–Jakob disease (CJD) via interventional procedures’

OPA/NAC: o-phthalaldehyde/N-acetyl-L-cysteine

PO: posterior ophthalmic

sCJD: sporadic Creutzfeldt-Jakob disease

SSD: sterile services department

TSEs: transmissible spongiform encephalopathies

UCHL: University College Hospital London

vCJD: variant Creutzfeldt-Jakob disease

Executive summary

Choice Framework for local Policy and Procedures (CFPP) 01-01 offers best practice guidance on the whole decontamination cycle including the management and decontamination of surgical instruments used in acute care.

Part A covers the policy, management approach and choices available in the formulation of a locally developed, risk-controlled operational environment. The technical concepts are based on European (EN), International (ISO) and British (BS) Standards used alongside policy and broad guidance. In addition to the prevention of transmission of conventional pathogens, precautionary policies in respect of human prion diseases including variant Creutzfeldt-Jakob disease (vCJD) are clearly stated. Advice is also given on surgical instrument management related to surgical care efficiencies and contingency against perioperative non-availability of instruments.

Part B covers common elements that apply to all methods of surgical instrument reprocessing such as:

- test equipment and materials

- design and pre-purchase considerations
- validation and verification.

Part C covers standards and guidance on steam sterilization.

Part D covers standards and guidance on washer-disinfectors.

Part E covers low temperature (non-steam) sterilization processes (such as the use of vapourised hydrogen peroxide gas plasmas and ethylene oxide exposure).

CFPP 01-01 Part A supersedes Health Technical Memorandum 01-01 Part A (2006).

CFPP 01-01 Parts B–E supersede Health Technical Memoranda 2010, 2031 and Health Building Note 13 Supplement 1, and partially supersede Health Technical Memorandum 2030. Only washer-disinfectors used for processing surgical instruments (and not those used in laboratories or for endoscopes) are covered in CFPP 01-01.

Contents

Preface

Introduction
Purpose
Essential Quality Requirements and Best Practice

Acknowledgements

Abbreviations

Executive summary

1 Design and pre-purchase considerations

1

Purpose of washer-disinfectors
Choice of washer-disinfector
 Categorisation of washer-disinfectors by nature of load to be processed
 Categorisation of washer-disinfectors by configuration/load handling type
 When is a washer-disinfector required?
 Choice of washer-disinfector
 Assessment of workload and throughput
 Rigid endoscopes
 Powered devices
Specification and contract
 Introduction
 CE marking
 Preparing a specification
 Construction materials
 Integral air compressors
 Integral calorifiers and tanks
 Dosing systems
 Door controls
 Loading systems
 Cleaning the washer-disinfector
 Steam
 Integral steam generators
 Condensate recovery
 Compressed air
 Drainage
 Hazardous effluents
 Machine ventilation
 Chemical additives: storage
Chemical additives
 Introduction
 Compatibility with the materials of construction of the washer-disinfector
 Compatibility with the process
 Compatibility with the items to be processed
 Compatibility with the quality of water
 Compatibility with other chemical additives
 Compatibility with subsequent decontamination processes

- Disinfectants
- Washer-disinfectors for surgical instruments and associated equipment
 - Single-chamber machines
 - Multi-chamber machines
 - Standard specifications
 - Load-handling equipment
 - Test connections
- Ultrasonic cleaners
 - Applications
 - Standard specifications
 - Additional specifications
 - Wash cycle
 - Type 1 and Type 2 ultrasonicators

2 Validation and verification

15

- Testing of washer-disinfectors
 - Interrelation of test programmes 13828
- Schedule of type tests and works tests
- Schedule of installation tests
 - Checks on ancillary equipment
 - Checks on washer-disinfectors
- Schedule of operational tests
- Performance Qualification tests
 - Surrogate devices
 - Cleaning efficacy tests
 - Disinfection
 - Load dryness tests
- Schedule of periodic tests
 - Introduction
 - Weekly safety checks
 - Yearly safety checks
- Test methods
 - Drainage
 - Venting system
 - Doors and door interlocks
 - Chemical dosing
 - Water vapour emissions
 - Instrumentation fitted to a washer-disinfector
 - Load carriers
 - Thermometric tests
 - Load dryness
 - Residual chemical additives
 - Air quality
 - Sound pressure test
 - Electromagnetic compatibility
 - Cleaning efficacy tests
 - Periodic tests
 - Disinfection efficacy tests
 - Automatic control test
- Validation and periodic tests for ultrasonic cleaners
 - Introduction
 - Test for ultrasonic activity
 - Reference test loads
 - PQ tests

3 Water supply**40**

Introduction	
Water hardness	
Ionic contaminants	
Microbial population	
Concentration of bacterial endotoxins	
Water treatment	
Chemical purity	
Microbial purity	
Pipework	
Water supply by-laws	
Water system tests	
Water samples	
Water quality tests	
Choice of method	
Water supply temperature	
Apparatus	
Method	
Alternative method (for periodic testing)	
Results	
Water supply pressure	
Apparatus	
Method	
Results	
Appearance	
Apparatus	
Method	
Results	
pH	
Apparatus	
Method	
Results	
Electrical conductivity	
Equipment and materials	
Method	
Results	
Total dissolved solids	
Conductivity method	
Results	
Evaporative residue method	
Results	
Hardness (as CaCO_3)	
Ion-selective electrodes (ISE) method	
Titrimetric method	
Results	
Chloride	
Ion selective electrode (ISE) method	
Silver nitrate titration method	
Results	
Heavy metals	
Method	
Results	
Iron	


Apparatus	
Method	
Results	
Phosphate	
Apparatus	
Method	
Results	
Silicate	
Apparatus	
Method	
Results	
Bacterial endotoxins	
Apparatus	
Method	
Analytical procedure	
Sources of error	
Safety information	
Results	
Total viable count	
Apparatus	
Method: sample collection	
Method: testing	
Results	
4 Operational management	57
Introduction	
Equipment damage	
Maintenance	
Introduction	
Planned maintenance programme	
Pressure Systems Safety Regulations	
Features requiring special attention	
Returning a washer-disinfector to service	
Troubleshooting	
Appendix A: Particular specification for washer-disinfectors used for processing surgical instruments	63
References	64

1 Design and pre-purchase considerations

Purpose of washer-disinfectors

- 1.1 Washer-disinfectors are used to clean and disinfect items intended for re-use. They may be used in relation to both medical devices and medicinal products as well as other items
- 1.2 The cleaning and disinfection process is intended to:
 - make the item safe for staff to handle;
 - make the item safe for use on a patient (after any necessary additional processing such as a terminal sterilization process), including, if appropriate, ensuring freedom from contamination that could lead to an erroneous diagnosis.
- 1.3 When items being cleaned and disinfected by a washer-disinfectors are intended to be used again without further treatment (such as a terminal sterilization process), before being reused the washer-disinfectors should produce an item that is microbiologically safe for its intended use.
- 1.4 When the items being cleaned/disinfected by a washer-disinfectors are intended to be subjected to further processing (such as a terminal sterilization process) before being re-used the disinfection stage in the washer-disinfectors should produce an item that is microbiologically safe to be handled during preparation for subsequent processing.
- 1.5 Washer-disinfectors processing involves two distinct stages: cleaning and microbial inactivation. The latter can be achieved by disinfection, sterilization or both (disinfection followed by sterilization). Washer-disinfectors are used to decontaminate items intended for re-use by subjecting the items to an automated process of cleaning and disinfection. This is shown in Figure 1.

Figure 1 Elements of the decontamination process



1.6 The efficacy of the cleaning stage of the process is of crucial importance to the successful outcome of the disinfection stage.

Choice of washer-disinfector

Categorisation of washer-disinfectors by nature of load to be processed

- 1.7 Washer-disinfectors may be categorised according to their intended use, i.e. the nature of the load they are intended to process:
- 1.8 Washer-disinfectors may be further defined both by their configuration/load handling system and by the nature of the operational process.
- 1.9 Washer-disinfectors may be used for processing a wide range of products used in clinical practice. Loads will include surgical instruments, anaesthetic equipment, bowls, dishes, receivers, utensils and glassware. Further information, including validation protocols can be found in BS EN ISO 15883.

Categorisation of washer-disinfectors by configuration/load handling type

- 1.10 Washer-disinfectors may also be categorized by their construction and the manner in which the load is processed through the machine.

Single-chamber machines

- 1.11 Single-chamber machines have one chamber in which the full range of process stages are carried out. They are machines in which all stages of the cycle are completed on one chamber load before another load can be processed in that chamber.
- 1.12 These might have either a single door through which both loading and unloading takes place or double doors with one door being used for loading and the other for unloading.

Multi-chamber machines

- 1.13 Multi-chamber machines have more than one chamber where separate stages of the cycle are performed in each chamber. The full range of process stages is only completed when the load is delivered from the final chamber. The segregation of cycle stages may be determined by specific User requirements but, typically, the chambers will be dedicated to cleaning, disinfection and drying. Different loads may be processed concurrently. An extra chamber for ultrasonic cleaning may also be

included. These machines will have doors at either end and at intermediate positions between chambers.

Conveyor machines

- 1.14 Washer-disinfectors of this type may find it difficult to meet the requirements of some validation and periodic tests and are thus not recommended. Current machines of this type should be replaced.

Thermal disinfection

- 1.15 In this process disinfection is achieved by the action of moist heat maintained on the surface to be disinfected at a particular temperature for a particular time. The combination of time and temperature should satisfy the A0 requirement of BS EN ISO 15883-1 for disinfection to be achieved.
- 1.16 The A_0 concept is defined in BS EN ISO 15883-1. It is defined as the disinfection effect resulting from an exposure to 80°C for a period of 1 s for an organism with a z-value of 10°C. Time/temperature bands meeting the requirements of an acceptable A_0 of 600 are as follows:

Disinfection temperature (°C)	Maximum allowable temperature (°C)	Exposure time period
70	75	100 minutes
80	85	10 minutes
90	95	1 minute

Drying

- 1.17 Drying can be an integral part of the cycle, usually by the circulation of hot air over the product, or provided as a separate drying cabinet. For some products, for example corrugated anaesthetic tubing, prolonged drying times are recommended and a separate drying cabinet can improve the productivity of the washer-disinfector.
- 1.18 A dry product can also be obtained by the flash evaporation of residual moisture from product items that are hot following a high temperature thermal disinfection stage.

When is a washer-disinfector required?

- 1.19 For many products used in healthcare practice there are two choices available:
 - products intended for re-use after they have been decontaminated and subjected to any

necessary reprocessing (for example terminal sterilization);

- single-use products, that is, to be discarded after use.

1.20 Products that are intended to be re-used should be decontaminated, in accordance with the manufacturer's instructions, in one of the following ways:

- manual cleaning followed by disinfection and sterilization;
- machine cleaning followed by disinfection and sterilization;
- automated machine decontamination incorporating cleaning.

1.21 See MHRA's Device Bulletin MDA DB 2006(04) – 'Single use medical devices: implications and consequence of reuse'.

Note

HSC 1999(178) – 'Variant Creutzfeldt-Jakob disease (vCJD): minimising the risk of transmission' states that "medical devices designated as single-use should never be re-used".

Choice of washer-disinfector

1.22 The choice of washer-disinfector should be determined by the nature of the loads to be cleaned and disinfected.

1.23 Purchasers should be aware that items suitable for a particular type of washer-disinfector might still require different operating cycles, which need to be specified before purchase.

1.24 Guidance on the modification of operating cycles to suit different loads is given in [Chapter 4, 'Operational management'](#). Advice on individual cases should be sought, if necessary, from the AE(D) before any decision is made. Use should be made of the 'Particular specification' for washer-disinfectors given in [Appendix A](#).

1.25 Once the type of washer-disinfector has been decided, preliminary enquiries should be made with a number of manufacturers. The use of the 'Particular specification' in Appendix A will enable data provided by the tenderer on technical points as well as financial data to be compared. Not only will this enable the purchaser to confirm the acceptability of current services, spatial requirements and portage but it will enable a like-

for-like tender analysis to be made. Tender analysis will be best achieved by formalizing tender comparison with respect to performance and cost in all key areas. Qualifying statements by the tenderer should be taken into account and their effect on tender content or eligibility should be made prior to a choice being made. See Appendix A for the 'Particular specification' for washer-disinfectors.

1.26 The dispensed volume of the chemical additives (e.g. detergents), including the accuracy and reproducibility of the dosing system(s), should be specified.

1.27 For washer-disinfectors employing jet-washing systems, the pump pressure and water flow are also critical variables. Water quality will also be a critical variable where water is used for cleaning or disinfection. Where drying is performed air temperature and air-flow rate are critical variables.

1.28 For ultrasonic cleaners, the frequency, amplitude and power are critical variables. For cleaning efficacy the following parameters may also be included: water temperature, detergent type and detergent concentration.

1.29 In all cases, the duration of each process stage should be determined with sufficient accuracy to ensure that consecutive cycles have the same efficacy.

Assessment of workload and throughput

1.30 Advice on equipment numbers and instrument throughput is included within Health Building Note 13 – 'Sterile services department'.

Rigid endoscopes

1.31 Many rigid endoscopes and most of the reusable surgical accessories used for minimal access therapy (MAT) can withstand steam sterilization and may be processed through washer-disinfectors employing a thermal disinfection stage to make them safe to handle during packing etc.

1.32 The washer-disinfector should be designed or adapted to ensure that during the flushing, cleaning, disinfecting and drying stages, process fluids (including drying) flow through the lumen(s) of the device.

1.33 Many washer-disinfectors can be equipped with dedicated load carriers to process rigid endoscopes, gas cannulae etc. There are also a number of

dedicated endoscope cleaners including ultrasonic cleaners.

- 1.34 The capacity of the washer-disinfector should be assessed on the number of items of each type that can be processed in a single load.

Powered devices

- 1.35 Some powered devices rely on the flow of a pressurised fluid or compressed air for operation (for example, dental handpieces, phaco handpieces, orthopaedic saws and drills). If such devices are to be processed through a washer-disinfector and the internal surfaces are to be flushed, they should be connected via dedicated connectors to ensure that process fluids flow through the lumen(s) of the device during the flushing, cleaning and disinfection stages. Flow and pressure through lumens should be in line with the manufacturer's recommendations for the device. Specialist filtration of process fluids passed through fine lumens and powered devices may be necessary to avoid adverse effects on their performance.

Specification and contract

Introduction

- 1.36 This section discusses general specifications for washer-disinfectors and the steps to be taken in inviting tenders and issuing a contract.

CE marking

- 1.37 Washer-disinfectors are covered by a number of European Directives and are thus required to be in conformance. Relevant Directives include but are not restricted to: Medical Device Directive (93/42/EEC and 2007/47/EEC) Class II a, Electromagnetic Compatibility Directive (2004/108/EEC), Low-voltage Directive (2006/95/EC), Pressure Equipment Directive (97/23/EEC) and the Machinery Directive (98/37/EEC and 2006/42/EEC).

Preparing a specification

- 1.38 It is essential that the preparation of procurement specifications be carried out by a qualified and competent person. The purchaser should employ the services of an AE(D) for this purpose.
- 1.39 Purchasers should refer to BS EN ISO 15883 Parts 1 and 2, CFPP 01-01 Part B and the 'Particular specification' for washer-disinfectors included in

Appendix A when preparing a specification for a washer-disinfector.

Construction materials

- 1.40 Materials of construction should comply with the requirements of BS EN ISO 15883.

Integral air compressors

- 1.41 Washer-disinfectors might require a supply of compressed air for either the operation of valves and powered door systems and/or during the drying stage of the cycle.
- 1.42 When compressed air is intended to come into contact with the washed and disinfected product the compressed air supplied should be "medical grade", i.e. it should be oil and particulate free (see Health Technical Memorandum 02-01 – 'Medical gas pipeline systems').
- 1.43 Built-in air compressors should be suitable for the duty imposed upon them.

Integral calorifiers and tanks

- 1.44 All integral calorifiers should conform to BS 853-1, and should be designed and constructed to allow thermal disinfection to be achieved throughout the calorifier and associated pipework before water/steam can be supplied to the washer-disinfector during the thermal disinfection and subsequent stages.
- 1.45 Water tanks within the washer-disinfector should be self-draining and located so that they are cleanable by the Operator and fitted with a drain down system which either works automatically when the machine is switched off or which is accessible to the User.
- 1.46 All tanks should be fitted with an overflow (see the Water Supply (Water Fittings) Regulations 1999).
- 1.47 When water is to be heated the heat source should be controlled by a thermostat and it should employ a heating medium as specified by the purchaser. The heat sources should be removable for replacement or maintenance purposes.

Dosing systems

- 1.48 The washer-disinfector should be fitted with not less than two systems for controlling the admission of chemicals (detergents, additives etc) and should be provided with the facility for at least one additional dosing system to be fitted.

- 1.49 Each dosing system should be provided with means to adjust the volume admitted. Access to the means of adjustment should require the use of a key, code or tool. The means of adjustment can be manual or automatic.
- 1.50 The stage(s) in the process cycle at which each dosing system admits chemical to the washer-disinfector should be under the control of the automatic controller.
- 1.51 Each dosing system should be provided with means to determine the volume admitted and the time within the operational cycle when the admission occurred. This data should be available to the operator.
- 1.52 Failure to admit the specified minimum volume should cause a fault to be indicated and a failure indication raised.
- 1.53 The accuracy and reproducibility of chemical dosing systems should meet the requirements of BS EN ISO 15883-2.
- 1.54 The washer-disinfector should be fitted with a system that will indicate when there is insufficient chemical(s) available for the next cycle and that cycle should not be allowed to start.
- successful operating cycle, i.e. without showing a fault.
- 1.58 If a fault develops, it should only be possible to open the loading door.
- 1.59 It should not be possible for an Operator to open or close a door at the opposite end of the washer-disinfector or for more than one door to be open at one time.
- 1.60 A visual display should be provided at both ends of the washer-disinfector to indicate when the cycle is in progress.
- 1.61 The indication “cycle complete”, or an equivalent indication, should be cancelled when the unloading door is unlocked, and the loading door should remain locked until the unloading door has been locked again.

Loading systems

Door controls

Control of manually operated doors

- 1.55 An explanation of the manual action required to lock the door should be provided for the Operator. In addition, if the unlocking procedure is not the reverse of the locking procedure, there should be an indication to the Operator of the manual action required to unlock the door. The indication should be clearly displayed either on the door or on its handle or handwheel. Explicit instructions should be displayed on the facing panel adjacent to the door or on the Operator's control panel.
- 1.56 The door mechanism should be such that the force to be applied by an Operator in order to either lock or unlock the door does not exceed 250 N at the intended point of grip.
- 1.62 The washer-disinfector should be provided with carriers to locate the load during the washing and disinfection process. When interchangeable load carriers and baskets are provided, each load carrier should be capable of being fitted and removed without the use of tools. Washer-disinfector loading systems should be designed with regard to the Manual Handling Operations Regulations 1992 (as amended).
- 1.63 When the washer-disinfector is supplied with a system for supporting the load/transferring the load into and out of the chamber, the following should apply:

- the load should be wholly supported and retained within the usable chamber space for the duration of the operating cycle;
- the force required by the Operator, either directly or by the application of a mechanical device supplied with the equipment, to remove the whole, or part, of the load from the chamber should not exceed 250 N when loaded and operated in accordance with the manufacturer's instructions;
- the load carrier should either be retained in the chamber by a mechanism which is only released when the transfer system is in place, or remain stable when withdrawn for a distance equal to two-thirds of the chamber length, and be fitted with a retaining device, which has to be released if the load is to be withdrawn further.

Control of doors in a double-door washer-disinfector

- 1.57 In double-door washer-disinfectors, the control initiating the automatic cycle should be at one end only. When the loading door is closed and locked, it should not be possible to open the unloading door until the washer-disinfector has completed a

- 1.64 Means should be provided such that the transfer of the load into and out of the chamber does not cause damage and wear to the chamber.
- 1.65 Systems which cause high levels of local stress, e.g. point loadings, might also initiate corrosion in stainless steel materials.
- 1.66 The system used to support the load should be constructed from durable, corrosion-resistant materials and should withstand, without damage, the environment within the chamber.
- 1.67 The system used should neither prevent the attainment of the pre-set cycle variables nor the free drainage of water from the load and the penetration of water/steam into the load. The load carrier(s) should be designed so that they cannot be mis-positioned in a manner that will prevent such attainment.
- 1.68 Any accessory used for handling the load which can be used outside the washer-disinfector (e.g. a trolley) should remain stable when it is supporting its maximum design load and a force of 250 N is applied horizontally in any direction to the highest point of the load or accessory.
- 1.69 The trolley should be designed to allow the operator to align the trolley with the washer-disinfector for ease of loading and unloading.
- 1.70 The trolley should be provided with means to collect liquid residues from the load to prevent these from dripping onto the floor. The means provided should be detachable for cleaning and disinfection in a washer-disinfector or, where necessary sterilization at 134–137°C in a porous-load sterilizer.
- 1.71 The trolley should be provided with swivel wheels to facilitate manoeuvring.
- 1.72 The trolley should be provided with a parking brake.
- 1.73 The trolley should be designed to secure the load carriers on the trolley during loading and unloading, and while traversing a gradient at a slope of up to 1 in 20.
- 1.74 Trolleys intended for use with single door machines should be designed and constructed to facilitate cleaning and disinfection of the trolley between use for dirty and clean loads.
- 1.75 Load conveyors outside the washer-disinfector that are intended to, or might reasonably be expected to, come into contact with soiled/contaminated

goods should be designed and constructed to be easy to clean and disinfect.

Cleaning the washer-disinfector

- 1.76 The design, construction and operation of the washer-disinfector should ensure that during the process the surfaces of the chamber and the load carrier presented to the operator are cleaned and disinfected.
- 1.77 For manually filled and emptied cleaning machines with no disinfection cycle, for example stand-alone ultrasonic cleaners, the manufacturer should advise on the cleaning/disinfection method.

Steam

- 1.78 Steam may be used in a heat exchanger as a source of indirect heating for water or air to be used in the cleaning, disinfection or drying stages of a washer-disinfector operating cycle.
- 1.79 Steam may also be used to heat process water directly, or to heat the load directly during the thermal disinfection stage, and for this purpose may be supplied either with an integral steam generator or from an external (mains) supply.

Steam for indirect heating

- 1.80 Steam heat exchangers used for heating water or air may be of the shell, tube or plate design. In all cases, the steam supplied should be substantially free from non-condensable gases and free from oil since these contaminants seriously impair the efficiency of the heating process. The effect of thin films of air on the surface of the heat exchanger may increase heating costs by 25% or more.
- 1.81 The steam service should be designed to meet the maximum demand of the washer-disinfector, or bank of washer-disinfectors, while keeping the fall in pressure before the final pressure reducing system to no more than 10%.
- 1.82 Except for vertical rises between floors or at intermediate points on long runs, the pipework should have a continuous fall such that any condensate flows by gravity in the same direction as the steam. Air vents and steam traps should be fitted at each vertical rise.
- 1.83 The condensate discharge system should be sized to ensure that the high volume of condensate found during the initial stages of heating can be discharged without waterlogging the heat exchanger.

- 1.84 When the steam supply pressure at the inlet to the washer-disinfector exceeds the maximum value specified by the manufacturer, a pressure reducing valve should be fitted to the supply pipe at least 3 m from the washer-disinfector.
- 1.85 Careful attention should be paid to the siting of all pressure relief valves to ensure that the washer-disinfector is properly protected.
- 1.86 Relief valves and their discharge pipes should be large enough to prevent the pressure in the supply pipes rising to more than 10% above the design pressure for the heat exchanger. The discharge pipe should terminate in a safe position outside the building.
- 1.87 Steel and copper piping traditionally used for steam supply are acceptable for this application.
- 1.88 Excessive moisture in the steam supply will impair the heating efficiency of the heat exchanger and so should be avoided.

Integral steam generators

- 1.89 Some washer-disinfectors are equipped with small electrically heated steam generators to raise steam to heat the load directly for thermal disinfection.
- 1.90 They may be of the open-boiler type, which are so designed and constructed that they are unable to generate an internal pressure above atmospheric pressure. The design should ensure that, under a single fault condition, for example obstruction of the steam discharge port, the boiler cannot become pressurised.
- 1.91 Integral steam generators which are pressure vessels should be in accordance with PD 5500.

Condensate recovery

- 1.92 Condensate from steam heating systems (calorifiers, dryers) and steam traps on the pipeline is suitable for recovery and should be returned to the steam generating plant when recovery is economically justifiable. It should not be possible to contaminate instruments if a failure of the heating coil occurs.

Compressed air

- 1.93 A compressed air supply might be required for the operation of controls and air-drying. If the washer-disinfector does not contain an integral air compressor, the air may be supplied from a piped service (mains supply) or from a local compressor.

Mains supply

- 1.94 If air is supplied by pipeline from a central air-compressor system, a Bourdon-type pressure gauge conformant to BS EN 837 should be fitted on the supply line to the washer-disinfector via an isolation valve.
- 1.95 A reducing valve, or other automatic device, should be fitted to reduce the pressure of air delivered to the washer-disinfector to no more than the maximum supply pressure specified by the manufacturer. A pressure relief valve will normally be required.

Local compressor

- 1.96 When it is not practical to obtain compressed air from a mains supply, a dedicated compressed air system should be installed to supply the washer-disinfectors.
- 1.97 The compressors might be too noisy to install with the washer-disinfector and might need to be located in a dedicated location away from noise sensitive areas.
- 1.98 Components of the compressed air system that require servicing and maintenance, such as dryers and filters, should be located where they are readily accessible for service or exchange.

Air quality

- 1.99 The quality of air can be critical for some applications and some washer-disinfectors will incorporate appropriate filters. When the purchaser is to be responsible for the provision of filtered air the CP(D) should ensure that the quality of air available meets the washer-disinfector manufacturer's specification or the requirements given below.
- 1.100 Air that could come into direct contact with the load, such as air used for drying the load or testing the free-passage of lumens, should be:
- oil-free (i.e. should have no more than 0.5 mg of oil per cubic metre of free air measured at 1013 mbar and 20°C; see ISO 554);
 - filtered to an efficiency of at least 95% when tested in accordance with BS 3928; and
 - free of bacteria (see Health Technical Memorandum 02-01).
- 1.101 Air for control purposes should be free of liquid water, filtered to 25 µm (5 µm for precision

controls) and lubricated with micro-fog oil particles of 2 µm or less.

Drainage

- 1.102 All effluent from a washer-disinfector is potentially contaminated and should be disposed of to the main drain. All washer-disinfector and associated equipment should be connected to the drain in a manner that provides backflow protection and complies with local regulations.
- 1.103 Effluent can originate from each of the stages of the process, which may include:
- flushing to remove gross contamination;
 - washing with detergent and/or enzymatic cleaners;
 - rinsing, with or without the addition of a neutralizer, rinse aid or instrument lubricant;
 - chemical disinfection or thermal disinfection;
 - post-disinfection rinsing;
 - drying particles.
- 1.104 Effluent from early stages (a) and (b) of the process might contain significant concentrations of organic contaminants and potentially infectious microorganisms. Effluent from the middle stages (b), (c) and (d) can contain some organic contaminants and potentially infectious microorganisms and high concentrations of process chemicals. Effluent from the latter stages (d), (e) and (f) can be at high temperatures (90°C–100°C).
- 1.105 Effluent from washer-disinfectors should pass via an air break into a tundish or tank before being discharged to drain. The air break should be preserved at all times to prevent the washer-disinfector and its associated pipework being contaminated by reverse flow from the drainage system.
- 1.106 Where a tank supplies water to a pump on the washer-disinfector, the overflow discharge from the tank should also include an air break.
- 1.107 The drainage system from the installation should be trapped and designed to pass the flow rate of water, air and condensed steam specified by the manufacturer, with account taken of the peak output during the operating cycle.
- 1.108 The drainage system should be designed to pass and maintain in suspension the solids removed

from the load during the flushing process. The minimum diameter of the drainage system should be greater than the maximum diameter of the most restricted section of the discharge from the washer-disinfector chamber.

- 1.109 Means should be provided to prevent, as far as possible, flash steam being liberated into the atmosphere or causing condensation on electrical equipment.
- 1.110 The discharge temperature from a washer-disinfector may be as high as 95°C.
- 1.111 The materials used for the construction of the discharge system should be chosen to withstand temperatures up to 100°C.

Hazardous effluents

- 1.112 The discharge of soil from washer-disinfectors should be regarded as being no more, but no less, hazardous than the discharge from any other sanitary appliance, e.g. a WC.
- 1.113 The discharge of process chemicals, including detergents and microbicides, may require special attention. The local water undertaking should be consulted before such chemicals are discharged into the drainage system, as it may be necessary to neutralize or inactivate them before discharge.
- 1.114 A sealed and vented drain should be used for the discharge of chemicals with a significant vapour pressure, determined at the maximum attainable temperature of effluent in the drain, which may be hazardous to health or a nuisance. Possible backflow from the drain should be prevented by the inclusion of a check valve and a vacuum breaker.

Machine ventilation

- 1.115 Washer-disinfectors are often run under a slight negative pressure to minimise the potential for the discharge of aerosols into the environment.
- 1.116 Washer-disinfectors not equipped with an air extract system may require siting under an extraction hood. The capture velocity in the vicinity of the process is a design issue based on a local assessment via design consultation. The advice of a specialist should be sought in designing the ventilation system.
- 1.117 Additional guidance is given in Part B of Health Technical Memorandum 03-01 – ‘Specialised ventilation systems in healthcare’.

Chemical additives: storage

- 1.118 Safe storage provision is needed for containers of chemical additives used in the washer-disinfector. These chemicals are frequently corrosive, irritant and toxic and provision should be made in, or adjacent to, the storage area for an emergency eye wash station and a source of running water to dilute any spillage. Reference should be made to local COSHH (Control of Substances Hazardous to Health) provisions.
- 1.119 In large installations with two or more machines, as might often be required in SSDs, bulk storage tanks for chemical additives required for the process might be preferred with a piped distribution system to each washer-disinfector.
- 1.120 For each chemical additive to be used, there should be two storage tanks in parallel (one of which might be a small reserve tank):
- to permit cleaning and maintenance of the of the large tank without interrupting the use of the washer-disinfectors;
 - to facilitate segregation between separate batches of chemical additive; and
 - to allow for an orderly change to a different formulation if required.
- 1.121 The liquid concentrates are often viscous and chemically aggressive. The pipework, valves etc used for the distribution of these chemicals will need to withstand the corrosive effects of these materials. Advice should be sought from the manufacturer of the chemical additives on suitable materials, construction and pumping systems for the distribution system.

Chemical additives

Introduction

- 1.122 Chemical additives are not necessary for all applications; while they can enhance the removal efficacy it is then necessary to remove them during the rinsing stage. For applications in the laboratory and in the preparation of components and equipment used in manufacturing medical devices and medicinal products, chemical additives should be avoided if their use is not essential.
- 1.123 In choosing the various chemical additives for effective cleaning and disinfection, the formulation of each chemical additive should be compatible with:

- the materials of construction of the washer-disinfector;
- the process being operated in the washer-disinfector;
- the quality of water available;
- the items to be processed and their intended use;
- any other additives to be used in the washer-disinfector process;
- any intended subsequent decontamination process (e.g. sterilization).

1.124 The required concentration should be accurately and reproducibly generated by the dosing system(s) on the washer-disinfector.

1.125 It is not sufficient to determine only the compatibility of the principal active constituents, as the precise formulation of the chemical additive will affect its compatibility.

Compatibility with the materials of construction of the washer-disinfector

1.126 The cleaning formulation should have no long-term effects on the components of the washer-disinfector.

1.127 Chemical additives that can be absorbed into, or adsorbed onto, surfaces of the washer-disinfector (e.g. plastic pipework) may be carried over into subsequent stages of the process (see also next section and [paragraph 1.134, 'Compatibility with subsequent decontamination processes'](#)).

Compatibility with the process

1.128 The performance of the additive should be matched to the physical characteristics of the operating cycle (e.g. jet-washing action systems require low foam detergents), if the washing action is not to be impaired.

Compatibility with the items to be processed

1.129 The chemical additives used should be compatible with the materials of which the load items are constructed and should not cause chemical or physical damage, e.g. phenolic compounds used in detergents and disinfectants might cause material changes in rubber and plastics, while the anodic coating on the surface of anodized aluminium is removed by strongly acid or strongly alkaline

compounds. Use of a rinse-free cycle for laryngeal masks may be considered.

- 1.130 The chemical additives used should be readily removed from the load items by rinsing with water and should be biologically compatible with the intended use of the load items. Chemical additives that are intended to persist on the surface of items processed through the washer-disinfectant (e.g. lubricants) should be biologically compatible with the intended use of the load items.

Compatibility with the quality of water

- 1.131 The activity of many detergents and disinfectants are seriously impaired by hard water.
- 1.132 Detergent formulations intended for use only with soft water might give rise to precipitation if used with hard water, particularly at elevated temperatures. Once this precipitation has occurred on the surfaces of the washer-disinfectant or the load it is particularly difficult to remove (see paragraph 1.129, 'Compatibility with the items to be processed').

Compatibility with other chemical additives

- 1.133 The additives used should be both compatible with other chemicals used in the same process stage and, as far as may be practicable, with those used in preceding and subsequent stages to minimize the adverse effect of any carryover.

Compatibility with subsequent decontamination processes

- 1.134 Chemical additives that might persist on the surface of items processed through the washer-disinfectant should be compatible with any subsequent decontamination process that might be required, such as terminal sterilization. An in-process instrument lubricant that deposits a lubricant film on all surfaces of the instrument should only be used if it has been demonstrated to be compatible with any subsequent sterilization process.
- 1.135 The specified concentration of chemical additives should be attained in order to ensure effective processing. Too little will impair the process, while too much is wasteful, might impair the process and might contribute to unacceptably high residual levels.
- 1.136 Suppliers of chemical additives should provide product data sheets and material safety data sheets

for the products supplied. These should include details of biocompatibility studies. Reference should be made to local COSHH provisions.

- 1.137 Suppliers of chemical additives normally provide details of the analytical methods that can be used to detect residual concentrations of product. The sensitivity of the method should be sufficient to determine the presence of the compound below the level at which any adverse biological reaction would be determined.

Detergents

- 1.138 For most applications, mild alkaline detergents in the pH range 8.0–11.0 are preferred. Alkalinity improves the efficacy of detergents both by enhancing their inherent cleaning capabilities, neutralizing and helping to remove acid soils, emulsifying oils and fats and peptidising proteins, and by synergistic action with other detergent compounds. Many surfactants work better in the presence of alkaline “builders” such as sodium tripolyphosphate (STPP).
- 1.139 Cleaning agents for use in washer-disinfectors should be:
- liquid (to facilitate accurate dispensing);
 - non-abrasive;
 - low foaming;
 - free rinsing;
 - biodegradable.
- 1.140 Detergents should not contain:
- artificial colouring agents;
 - optical brighteners;
 - perfumes;
 - halides at an in-use concentration greater than 120 mg/L;
 - fatty soaps, glycerine or lanolin;
 - toxic residue.

Enzymatic cleaners

- 1.141 Enzymes are organic catalysts through which the normal metabolism of most living organisms takes place. Although produced by living organisms they are not themselves alive. They are large molecules whose steric configuration (shape) affords them the ability to catalyse many reactions in the living cell.

- 1.142 Enzymes are classified into groups depending on the nature of the chemical reaction that they catalyse. Generally the enzymes used in enzymatic cleaners are hydrolases, i.e. they promote the hydrolysis of the substrate with which they interact.
- 1.143 Enzymatic cleaners are themselves proteins and can be sensitising or allergenic agents.
- 1.144 A considerable proportion of the soiling found on medical items contains proteins that act as binding agents. Particulate dirt can be bound by the coagulation of these proteins on the surface.
- 1.145 If the binder proteins can be broken down into a simpler molecular form this binding action is destroyed and the bound soil, as well as the protein, can be released from the surface.
- 1.146 Formulations will often include buffering agents to maintain the pH within the preferred range. Most enzymes have an optimum pH at which their activity is greatest and a pH at which the enzyme itself is most resistant to thermal degradation, although these two values are not necessarily the same.
- 1.147 For example, the proteolytic hydrolase derived from *Bacillus subtilis*, subtilisin A, withstands temperatures up to 60°C and displays its greater stability at pH 9.4.
- 1.148 The importance of the enzymatic solution being at the correct temperature and pH, as well as being used for the specified contact time, cannot be too strongly emphasised.
- 1.149 Enzymes are not themselves cleansing agents. A properly balanced detergent might still be needed to remove the simpler molecular forms resulting from the enzymatic action.
- 1.150 Enzymatic formulations for cleaning solid surfaces are available in two forms:
- a pre-soak formulation which is used to digest proteinaceous soil and is followed by normal washing process using detergent;
 - a combined enzyme and detergent formulation.

Cleaning additives for ultrasonic cleaners

- 1.151 Only detergents specifically intended for use in ultrasonic cleaners should be used. The use of other detergents can impair rather than enhance the cleaning process.

Rinse aids

- 1.152 Rinse aids are generally formulated from surfactants and are designed to make the water more free rinsing.

Lubricants

- 1.153 The addition of oil-based compounds to the cleaning process is wrong in principle, as they cause contamination over the entire cleaned surface. Only the water-soluble type should be used. Mineral oils have poor biocompatibility and can inhibit the penetration of steam or sterilant gases on terminally sterilized product.
- 1.154 Lubrication should only be applied to those areas where it is required during the inspection/packing process after thorough cleaning of the instrument

Disinfectants

Choice of disinfection method

- 1.155 Thermal disinfection using moist heat is the preferred method and should be used whenever it is compatible with the load to be processed.
- 1.156 Temperatures in excess of 65°C and up to 95°C (or in some cases 100°C) can be used for disinfection; the lower the temperature the longer the exposure time in order to obtain the same reduction in microbial population. The thermal disinfection process is reliable, reproducible, free from toxic residues and capable of easy and economical physical monitoring and recording.
- 1.157 Chemical disinfection should only be used for products that cannot be treated using thermal disinfection methods.

Materials compatibility

- 1.158 The disinfectant should not cause damage to either load items or the washer-disinfector in which it is used. Damage that may occur with incompatible disinfectants includes corrosion, embrittlement or swelling of plastics, degradation of lens cement in optical systems, etc. The potential for electrolytic attack to occur as a result of different metals in the load and the washer-disinfector coming into contact, via a powerful electrolyte, should not be overlooked.
- 1.159 The material of construction of the washer-disinfector and of the items in the load should not inhibit the disinfectant.

Safety of disinfectants

1.160 Many of the compounds that are most effective as disinfectants are potentially human health hazards. Employers are required by law to do everything that is reasonably practicable to protect the health of their workers. The safe use of these compounds is covered by the Control of Substances Hazardous to Health (COSHH) Regulations (as amended). It should be noted that some disinfectants require dilution or activation before use.

Washer-disinfectors for surgical instruments and associated equipment

- 1.161 This section gives recommendations for washer-disinfectors used for cleaning and disinfecting surgical instruments and associated equipment including anaesthetic accessories, bowls, dishes, receivers, utensils and glassware.
- 1.162 The guidance given here assumes that the washer-disinfector is to be used during decontamination of medical devices and that the essential requirements of the EU directives are met. Further information and requirements can be found in BS EN ISO 15883 Parts 1 and 2.

Single-chamber machines

- 1.163 Single-chamber cabinet washers for surgical instruments and associated equipment may be designed to accept interchangeable load carriers, typically with rotating spray arms or other devices to ensure a uniform wash action with several layers of load items.
- 1.164 The spacing between layers should be designed to accommodate a number of wire mesh baskets full of instruments or should be more widely spaced to accept and correctly position large bowls, instrument trays, reusable rigid containers and similar items. Spacing should also allow for anaesthetic accessories to be located and processed, as well as specialist carriers with connections for particular instruments such as rigid endoscopes, MAT instruments, lumen devices and powered devices.
- 1.165 Since all stages of the cycle take place in the same chamber, it is not possible to get physical separation between the dirty and clean stages of the cycle. Assurance that the load will not be recontaminated is dependent upon the efficacy of the cleaning and disinfecting stages in

decontaminating the interior of the washer-disinfector as well as the load.

Multi-chamber machines

- 1.166 Continuous process washers, other than those designed as automatic ultrasonic cleaners only, are usually designed to accept interchangeable load carriers (see paragraph 1.163, 'Single-chamber machines'). Where the first chamber is used for first (that is, cold) wash only there may be a build-up of deposits, which will require additional cleaning at set intervals.
- 1.167 Compared with single-chamber machines they have a higher throughput and, for a similar process, achieve some decrease in overall cycle time.
- 1.168 Since the load is moved through the machine as the cleaning and disinfection cycles proceed it is possible to get excellent physical separation between dirty and clean load items.
- 1.169 There may be some loss of operational flexibility when this type of machine is used for several applications at a time; for example, if it is used to process anaesthetic accessories, the increased drying time, necessary for this application, will slow the passage of other loads passing through the washer-disinfector.
- 1.170 Washer-disinfectors of this type are large, expensive pieces of equipment and they should only be used in centralized production units.

Standard specifications

- 1.171 Washer-disinfectors for surgical instruments and associated equipment should conform to BS EN ISO 15883 Parts 1 and 2 and the safety specifications in BS EN 61010-1.
- 1.172 Use should be made of the 'Particular specification' for washer-disinfectors in [Appendix A](#).

Load-handling equipment

- 1.173 A number of different types of carrier should be used to accommodate the range of items to be processed. The range of carriers that might be necessary for an SSD include:
- a multi-layer carrier for instruments in wire mesh baskets (wire-mesh baskets to include a number with retaining systems for small instruments);

- a two layer carrier for small hollowware and instrument trays;
- a single layer carrier for large bowls, Edinburgh trays etc;
- a rigid endoscope/MAT instrument/lumen device carrier;
- an anaesthetic accessories carrier;
- sterilization containers;
- carrier with dedicated connectors for powered devices.

- 1.174 The load carriers should protect instruments from mechanical damage during the wash process and should also orientate the instruments to facilitate proper cleaning providing, when necessary, a direct connection between the water flow and the lumen of the load item.
- 1.175 The specification for load handling equipment should include the provision of appropriate tabling to permit sorting of instruments and loading of load carriers and, after processing, the unloading of load carriers.
- 1.176 When double-door washer-disinfectors are used, means should be provided to return load carriers from the unloading to the loading end. Where this passes through the wall between the packing room and decontamination room, there should be a pass-through hatch with interlocked doors.

Test connections

- 1.177 Test connections should be provided to permit the connection of thermocouples and pressure transducers to be used during validation and periodic testing.
- 1.178 When additional monitoring is provided, a separate test connection should be provided for each sensor to permit periodic verification of the installed system by comparison with a calibrated test sensor.

Ultrasonic cleaners

- 1.179 The guidance given here assumes that the washer-disinfector is to be used to clean and disinfect medical devices and that the essential requirements of the EU directives are to be met.
- 1.180 Ultrasonic cleaners may be Type 1 or Type 2 machines (ultrasonication may be integrated into a Type 2 washer-disinfector). They are mainly used

as a preliminary cleaning stage prior to processing through a washer-disinfector.

- 1.181 Ultrasonic cleaners work by exposing the items to be cleaned to high frequency sound waves in the liquid cleaning medium. The high frequency sound waves are generated within the liquid by the vibration of one or more surfaces of the bath, which is caused by one or more transducers bonded to the outer surface(s). The transducers convert electrical power into vibrations of the required frequency and amplitude. The highly effective cleaning action occurs as a result of the penetrative agitation caused by cavitation, the rapid formation and collapse of tiny bubbles within the liquid, which are generated by the high frequency sound waves.

Applications

- 1.182 Ultrasonic treatment is particularly suitable for cleaning instruments of high-grade steel. Delicate instruments such as micro-surgery instruments and dental instruments can be effectively cleaned with little risk of damage.
- 1.183 Ultrasonic treatment is also particularly effective for cleaning instruments that have deep interstices that may be contaminated with body tissues, for example reamers, drills and burrs.
- 1.184 When combined with appropriate connection to an irrigation or flushing system, ultrasonicators are also effective for cleaning internal and external surfaces of cannulated instruments.
- 1.185 Ultrasonic cleaners are less effective when used to clean plastic and similar readily compressible materials since they absorb much of the ultrasonic energy.

Standard specifications

- 1.186 Safety specifications for ultrasonic cleaners are included in BS EN 61010-1.

Additional specifications

- 1.187 The ultrasonic cleaner should be fitted with means to drain the tank with the cleaner in situ. The tank should be free draining so that no pools of water are left in the tank after draining.
- 1.188 The tank should be heated electrically and the heaters should be thermostatically controlled.
- 1.189 The ultrasonic cleaner should be fitted with a timer to control the duration of exposure.

- 1.190 The ultrasonic cleaner should have a lid; the lid should be interlocked with the operating system to prevent normal operation if the lid is open and should fit securely to prevent the emission of aerosols when the cleaner is in operation.
- 1.191 The lid interlock should ensure that no part of the Operator's body can be immersed in the ultrasonic cleaner during operation.
- 1.192 The ultrasonic cleaner should be effectively insulated to prevent high frequency sound transmission at a power that could cause a health hazard. The casing and lid should provide adequate sound proofing so that harmonic frequencies within the audible range are not obtrusive.
- 1.193 The manufacturer will normally recommend the chemical additives (detergents/enzymatic cleaners) that are compatible with the process. Low foaming detergents should be used. Liquid detergents used for washing dishes ("washing-up liquid") are not suitable.
- 1.194 The manufacturer is obliged to specify how the cleaner may be disinfected. This might be by the provision of a high-temperature (for example 80°C) cycle option, or by means of a suitable disinfectant solution. In the absence of guidance the Microbiologist should be asked to advise on a suitable procedure.
- 1.195 The manufacturer is also obliged to specify the degassing time(s) to be used on start-up and, when necessary, between each load of instruments processed.

Wash cycle

- 1.196 The ultrasonic frequency used should be within the range 35 kHz \pm 5 kHz and the energy input used may range from 5 W L⁻¹ to 20 W L⁻¹.
- 1.197 Ultrasonic cleaners can be designed to operate at a single frequency, across a frequency range, or with a feedback control system claimed to adjust the frequency in response to the loading conditions.
- 1.198 For medical applications aqueous solutions should be used. Although ultrasonic cleaners containing aqueous solutions can be effective at temperatures up to 90°C it is normal practice to operate those for medical applications at temperatures between ambient and 40°C. This minimizes the rate of coagulation of proteinaceous material in the soiling and is compatible with the use of enzymatic cleaners, many of which are rapidly destroyed at higher temperatures.

Type 1 and Type 2 ultrasonicators

Type 1 ultrasonicators

- 1.199 A mechanical lifting device should be used when the ultrasonicator is intended to process heavy sets of instruments.
- 1.200 The load container, usually a wire mesh basket, should be of appropriate size for the longest instrument to be processed.
- 1.201 When it is intended to process micro-surgical instruments or instruments with fine points the load handling equipment should provide means of retaining these in position so that the points are not blunted by the impacts resulting from fine mechanical shaking.

Type 2 (continuous process) ultrasonicators

- 1.202 Continuous process washer-disinfectors may incorporate an ultrasonic cleaning stage within the cycle programme.
- 1.203 Ultrasonic cleaners are also available in continuous process format with a thermal disinfection stage and with the option to provide a hot air drying stage.
- 1.204 Ultrasonic cleaners with a solvent drying stage are no longer commercially available since the solvents used were CFCs, which are now prohibited under the Montreal Protocol.
- 1.205 If complex tabling or conveyors are required these should be specified, and preferably illustrated with a sketch plan, when seeking tenders.
- 1.206 When it is intended to process micro-surgical instruments or instruments with fine points the load handling equipment should provide means of retaining these in position so that the points are not blunted by the impacts resulting from fine mechanical shaking.
- 1.207 The ultrasonic cleaner should be fitted with a temperature indicator; provision should be made for a recorder to be fitted if requested by the purchaser.
- 1.208 The ultrasonic cleaner should be fitted with an indicator to show the power consumption (in watts), or electrical demand (in amps) of the ultrasonic transducers; provision should be made for a recorder to be fitted if requested by the purchaser.

2 Validation and verification

Testing of washer-disinfectors

- 2.1 Washer-disinfectors are used to carry out the processes of cleaning and disinfection consecutively. These processes require validation to demonstrate they will consistently clean and disinfect all instrument scenarios likely to be presented including those most difficult to clean and disinfect. Choice of test loads and dryness of soils need careful consideration to ensure that this objective is met.
- 2.2 In some instances a visual inspection for residual contamination may be considered sufficient for monitoring the adequacy of the cleaning process before use. However, this is not true in all cases; for example, visual inspection will not detect soiling on the internal surfaces of instruments with lumens and will not detect low, but potentially significant, concentrations of soiling (for example proteins) or residual chemical additives from the WD remaining on load items.
- 2.3 There is no simple method to verify by inspection or test the efficacy of the disinfection process on product prior to use.
- 2.4 In consequence, cleaning and disinfection processes should be validated before use, the performance of

the process should be monitored during routine use, the calibration of controls and instrumentation should be verified, and the equipment should be subjected to a suitable maintenance programme.

- 2.5 The control protocols recommended in this section provide the means for ensuring that the washer-disinfector is fit for its intended purpose and includes tests and checks carried out during manufacture, after delivery, during validation and periodically thereafter. Tests are also recommended before a washer-disinfector is returned to service after repairs that affect one or more components which influence the attainment of critical process control variables or after modification.

Interrelation of test programmes

- 2.6 The tests recommended in this section are intended for use in type tests, works tests, commissioning (installation and operational tests), PQ (thermometric tests, microbiological tests, cleaning efficacy tests and load dryness tests) and routine periodic tests.
- 2.7 The interrelationship of the various test programmes, the place where they would usually be conducted and the responsibility for conducting the tests are shown in Figure 2.

Figure 2 Interrelationship of test programmes

Location	Production of WD		Responsibility
	Serial Production	'One-off' WD	
Factory	Type Test → Works Test	Type Test	Manufacturer
On-Site	VALIDATION Commissioning Checks & Tests	INSTALLATION OPERATIONAL	
		PERFORMANCE QUALIFICATION	
		PERIODIC ROUTINE TESTS ANNUAL REVALIDATION TESTS	User

2.8 The programmes of tests should be applied to all washer-disinfectors where relevant. Details are given under the test schedules for particular types of washer-disinfector.

Schedule of type tests and works tests

2.9 The manufacturer will carry out type tests on representative samples of washer-disinfectors in serial production to demonstrate compliance of the washer-disinfector design with BS EN ISO 15883 Parts 1 and 2.

2.10 The manufacturer will carry out works tests on each washer-disinfector before it leaves the manufacturing site to ensure that each washer-disinfector meets the specification. These tests should be as given in BS EN ISO 15883 Parts 1 and 2 including any additional tests required by the User and AE(D).

2.11 The manufacturer should make the results of type tests and works tests available to the purchaser on or before delivery of the washer-disinfector.

2.12 It may be necessary for the purchaser, or their representative, to visit the manufacturer's works to witness works testing. The advice of the AE(D) should be sought.

Schedule of installation tests

Checks on ancillary equipment

2.13 Ancillary equipment should, whenever practicable, be installed and commissioned before validation of the equipment begins.

2.14 When the checks on ancillary equipment require the washer-disinfector to be in operation, the CP(D) should carry them out in co-operation with the contractor for the washer-disinfector.

2.15 The contractor for the equipment is not responsible for the correct functioning of services and ancillary equipment unless this was agreed in the purchase contract.

Engineering services

2.16 Checks should be made for the following services:

- a. the engineering services should be installed correctly, should be adequate to meet the demands of the decontamination equipment, should not leak and all necessary isolating valves/switches and test points should be installed;

- b. the drains should remove effluent effectively when all plant in the vicinity, including the decontamination equipment, is connected and operating;
- c. the water treatment plant (if fitted) should operate correctly and the quality of water supplied for each stage of the process should be in accordance with the specification;
- d. the water economy system (if fitted) should operate correctly;
- e. the provision for storage, handling and connection to the washer-disinfector for all process chemicals should meet the requirements for safe handling of potentially hazardous chemicals;
- f. the exhaust ventilation and/or condenser unit fitted to the washer-disinfector should be adequate to remove the hot, humid air evolved from the washing, thermal disinfection and drying and unloading processes;
- g. for washer-disinfectors employing volatile process chemicals, the exhaust ventilation should maintain the environmental concentration below any limit specified for occupational exposure and the discharge should be to a safe place.

Checks on washer-disinfectors

Preliminary checks

- 2.17 Check that the electrical equipment on the equipment is correctly connected to the electrical service. Carry out the following electrical tests:
- a. insulation resistance;
 - b. phase sequence (for three-phase installations);
 - c. polarity;
 - d. bonding and earth continuity;
 - e. emergency stop.
- 2.18 After the equipment has been installed, check that the following recommendations are met:
- a. the manufacturer has supplied all the documents specified in the contract;
 - b. the equipment has been supplied and installed in accordance with the contract;

- c. calibration verification certificates for the measuring instruments and controller(s) on the equipment have been supplied;
- d. no defects are apparent from a visual inspection of the equipment;
- e. all supports, bases and fixings are secure and without imposed strain from service connections;
- f. thermal insulation is in good condition and securely attached;
- g. security and settings of door safety switches are in compliance with data supplied by the manufacturer;
- h. keys, codes or tools required to operate locked controls and control over-rides have been supplied, operate correctly and only operate the control for which it is intended; and cannot unlock controls on other machines in the vicinity;
- j. loading conveyors and trolleys, load carriers and load baskets are effective and safe in use;
- k. All appropriate connections for irrigation are available.
- e. Throughout the cycle there are no leaks of water, steam aerosols, toxic chemicals, air, gas or effluent.
- f. There is no evidence of interference to or from other equipment connected to the same services.
- g. There is no evidence of electromagnetic interference to or from other equipment.
- h. Operation and reading of all instruments appear to be satisfactory.
- j. The temperature of surfaces routinely handled by the operator does not exceed those specified in Chapter 1, 'Design and pre-purchase considerations'.
- k. The effluent temperature does not exceed that specified in Chapter 1.

2.20 At the end of the cycle check that the following recommendations are followed.

- a. The door opening system cannot be opened until the cycle has been completed without causing the cycle to abort and a fault/ incomplete cycle indication produced, that is, the automatic controller has operated in accordance with its specification.
- b. For systems incorporating one or more cycle stages at pressures 200 mbar above or below atmospheric pressure:
 - (i) the door opening system cannot be operated until the chamber has been vented to atmosphere and the chamber pressure is within 200 mbar of atmospheric pressure;
 - (ii) the door retaining parts cannot be released until the seal between the door and chamber has been broken, and the chamber is effectively vented to atmospheric pressure.
- c. The door interlock system is either fail-safe or is fitted with at least two independent interlocks. Failure of one interlock, or any one service, does not allow the door to be opened when conditions within the chamber would cause a hazard (for example, pressure in excess of 200 mbar or unacceptable level of chemical vapours).
- d. The automatic controller has operated in accordance with the parameter values determined at validation.

Functional checks

2.19 During an operating cycle, with an empty chamber, check that the following recommendations are followed (several cycles may be necessary to complete all the checks):

- a. The selection of automatic or manual control is by key code or tool. The selection of one control mode inactivates the other control mode.
- b. Under automatic control, water, steam, compressed air or chemical additives cannot be admitted into the chamber, and the operating cycle cannot start until the door is closed (locked and sealed).
- c. Under manual control the operator can advance the cycle only sequentially through each stage. Any stages designed to remove chemical additives from the chamber and load cannot be circumvented.
- d. Throughout the cycle the indicated and recorded values of cycle variables are within the limits specified by the manufacturer.

Response to external faults

- 2.21 The decontamination equipment should be checked to ensure it reacts correctly and safely when exposed to a number of external fault conditions; that is, a safety hazard is not created and a false indication of satisfactory completion of a cycle is not obtained.
- 2.22 During each stage of an operating cycle, check the response of the decontamination equipment to the following simulated faults (as appropriate to the type of machine):
 - a. operation of the emergency stop button;
 - b. power failure;
 - c. water pressure too low;
 - d. water pressure too high;
 - e. steam pressure too low;
 - f. steam pressure too high;
 - g. compressed air pressure too low;
 - h. compressed air pressure too high;
 - j. failure of chemical additive supply (detergent, sterilant, disinfectant etc);
 - k. failure of extract ventilation.

Schedule of operational tests

- 2.23 Full testing protocols for washer-disinfectors for surgical instruments, anaesthetic equipment, bowls, dishes, receivers, utensils and glassware can be found in BS EN ISO 15883 Parts 1 and 2.
- 2.24 Protocols for washer-disinfectors are shown in Table 1.

Table 1 Testing protocols for washer-disinfectors

Installation tests – contractor
1. Verification of calibration of WD instruments
2. Automatic control test
3. Water quality tests – hardness
4. Water supply temperature
5. Water supply pressure
Operational tests – CP(D)
1. Weekly safety checks
2. Automatic control test
3. Verification of calibration of WD instruments

4. Water system: <ul style="list-style-type: none"> – overflow test – volume of water used per stage – chemical purity – bacterial endotoxins
5. Drainage: <ul style="list-style-type: none"> – drain seal integrity – free draining – efficacy of discharge – pipework residual volume
6. Venting system: <ul style="list-style-type: none"> – steam venting – load contamination from ductwork
7. Doors and door interlocks <ul style="list-style-type: none"> – In-cycle interlock – Double-door WDs – On sensor failure – Door opening force – Failed cycle interlock – Fault indication on sensor failure
8. Water vapour discharge test
9. Chemical additive dosing tests: <ul style="list-style-type: none"> – reproducibility of volume admitted – low level detection¹
10. Load carriers
11. Air quality
12. Cleaning efficacy test: <ul style="list-style-type: none"> – test soil – reference load
13. Chamber wall temperature test ²
14. Load carrier temperature test ²
15. Over-temperature cut-out test ²
16. Thermometric test for disinfection reference load ²
17. Load dryness test reference load ³
18. Test for ultrasonic activity
19. Process residue test

¹ Only required for WDs with automatic chemical dosing

² Only required for WDs with a thermal disinfection stage

³ Only required for WDs with a drying stage

NOTE: Checks are also required on the independent monitoring system, including calibration, limits and trip points/alarm actions

Performance Qualification tests

Surrogate devices

- 2.25 Many of the devices that constitute the most difficult loads to process in a washer-disinfectors, which therefore require PQ, are difficult to monitor

either thermometrically or microbiologically, and are in short supply and are extremely expensive; examples include fibre-optic endoscopes, videoscopes etc.

- 2.26 A surrogate device is a test piece designed and constructed to emulate the characteristics of a device to facilitate appropriate monitoring of the cleaning and disinfecting processes.
- 2.27 An example of a surrogate device might be a rigid endoscope emulated by a similar length of stainless steel tube of appropriate diameter and bore. The surrogate device can be constructed to incorporate the appropriate temperature sensors and so that it can be separated into sections to facilitate the evaluation of residual test soil or survivors from a microbial challenge.
- 2.28 The surrogate device should have similar geometry and thermal mass and, as far as practicable, should be constructed of the same materials and with the same surface finishes as the device it is designed to emulate.
- 2.29 When an instrument presents particular problems in validation the manufacturer of the instrument should be requested to provide details of the method by which they recommend that PQ studies should be performed.

Cleaning efficacy tests

Representative soiling

- 2.30 Cleaning efficacy tests are intended to demonstrate the ability of the washer-disinfector to remove or reduce, soiling and contamination that occurs during normal use of reusable items to acceptable levels.
- 2.31 Naturally occurring contamination shows considerable variation both in the nature and proportion of constituents and in the extent of soiling which can occur during use.
- 2.32 Test methods based on the detection of naturally occurring soiling are difficult to standardise and show poor reproducibility due to:
- variation in the composition of the soiling, which can affect the ease with which soiling is removed;
 - changes in sensitivity of detection that can occur due to variation in composition of the soiling;
 - variation in the extent of soiling.
- 2.33 A number of methods exist for estimating the residual levels of proteins on surgical instruments.
- 2.34 Common practice in the past has been to rely solely upon visual inspection to detect unacceptable levels of residual soiling. This method should not be used, as it has poor sensitivity, is very subjective and can be greatly influenced by a number of factors including the intensity and nature of the illumination in the inspection area.

Note

The development path for this guidance recognises the current work in the piloting of high sensitivity post-decontamination protein quantification tests. As these become validated and available, this guidance will be amended to incorporate their use.

Test soils

- 2.35 Artificial test soils are designed to simulate the nature of native soiling and to be equally or more difficult to remove.
- 2.36 Appropriate marker substances should be incorporated to provide improved sensitivity of detection.
- 2.37 Soil tests may be used to aid determination of loading levels.
- 2.38 Test soils avoid any hazard that might be associated with native soiling (for example blood-borne viruses), which can be of particular concern with the more extensive handling necessary for test work.
- 2.39 Worldwide, many different test soils have been specified for testing washer-disinfectors but they generally fail to meet the key criteria necessary for a test soil. These criteria include:
- a chemically defined formulation (the traditional soils use substances such as flour, wallpaper paste, fresh egg yolk, horse blood etc which introduce a significant variability);
 - a quantitative method of applying the test soil to the surfaces of all types of item to be processed;
 - a quantitative method of detection of soiling remaining after the washing–disinfection process;

- validation with a known relationship to native soiling for ease of removal, relevant residual levels etc;
- safe to handle, easy and economical to use.

Standard test soils

2.40 Test soils are defined in ISO/TS 15883-5 as follows:

Instrument type	ISO/TS 15883-5 Annex reference
Surgical instruments	Annex N – see description below
Anaesthetic accessories	Annex O

2.41 The procedures within the UK are described as follows:

Typical test soil

Constituents

2.42 The following ingredients are required:

- fresh egg yolk 100 mL;
- defibrinated blood, 10 mL (horse or sheep);
- dehydrated hog mucin 2 g.

Preparation and storage

2.43 Mix all the constituents together and agitate in a stomacher to give a liquid of uniform consistency.

2.44 Use immediately or store in an air-tight container at 2°C to 5°C for not more than one week.

Application and use

2.45 If the soil has been stored, allow it to equilibrate to room temperature before use.

2.46 The following apparatus is required:

- paintbrush, 25 mm in width, soft;
- disposable gloves, for example latex;
- drainage tray.

Method of application

2.47 The method is:

- Don the protective gloves. Apply the soil to the test pieces by fully immersing the items in the soil or, for larger items, applying an even coat of soil using the paint brush.
- Allow excess soil to drain from the items and allow them to dry at room temperature (15°C to 25°C) for not less than 30 minutes and not more than 2 hours.

Method of detection

2.48 Visual examination. There are alternatives available. The advice of the AE(D) should be sought.

Process residues

2.49 The nature and level of process residues that are of concern depend on the chemical additives and quality of water used during the process and the intended use of the washed and disinfected product.

2.50 The water used for the process can give rise to a number of chemical residues on processed items. The most obvious of these is limescale from hard water.

2.51 The water used for the process might also give rise to contaminants of microbial origin. Bacterial endotoxins, primarily derived from the cell wall of Gram-negative bacteria, may give rise to adverse (pyrogenic) reactions when introduced into the body.

2.52 Items intended for surgically invasive use or for the preparation or administration of parenteral fluids should have suitably low levels of bacterial endotoxins.

2.53 The chemical additives used during the process (detergents, rinse aids etc) might not be completely removed by the rinsing process. The residual level that may be tolerated depends upon the nature of the chemical and the intended use of the product. The supplier of any chemical agent used will normally provide data on the chemical composition of the chemical agent and the biocompatibility of the components of the chemical agent. The supplier will also normally provide details of the method of detection, which can be used to determine whether processed items are free from residuals at the specified levels.

Disinfection

Thermometric tests

2.54 Thermometric tests should be carried out for both thermal disinfection processes and chemical disinfection processes.

2.55 The requirements for thermal disinfection are defined in BS EN ISO 15883 parts 1 and 2.

2.56 Normally, microbiological testing is not recommended for thermal disinfection processes.

Load dryness tests

- 2.57 The presence of residual water on cleaned and disinfected items should be avoided as it can interfere with the correct functioning of the item, promote re-contamination and microbial growth, or prevent attainment of sterilizing conditions. In many cases these data will already be available from the published literature.
- 2.58 The ability of the washer-disinfector to dry the load can be evaluated either visually, when appropriate, or by drying to constant weight and determining the mass of residual water present at the end of the washer-disinfector process cycle.
- 2.59 Performance qualification tests are shown in Table 2.

Table 2 Performance qualification tests for washer-disinfectors (carried out by the CP(D))

Cleaning efficacy test for a full load of particular items not represented adequately by the reference load: – test soil
Load dryness test for a full load of particular items not represented adequately by the reference load
Process residues: chemical additives
For temperature of internal surfaces of processed devices (e.g. anaesthetic and respiratory tubing, lumen devices and powered devices), see BS EN ISO 15883-2 clauses 4.4 and 6.3.3
For cleaning efficacy test for internal surfaces of processed devices (e.g. anaesthetic and respiratory tubing, lumen devices and powered devices), see BS EN ISO 15883-2 clause 4.4 and 6.2
For verification of flow through lumen and powered devices, see BS EN ISO 15883-2 clauses 4.4, 5.1.1, 5.1.2 and 6.3.3

Schedule of periodic tests

Introduction

- 2.60 Periodic tests are carried out at daily, weekly, quarterly and yearly intervals. Whilst protocols are defined in BS EN ISO 15883 Parts 1 and 2, it is recommended that periodic testing is performed as defined in Table 3.
- 2.61 The yearly test schedule is identical to that required for revalidation. It contains the tests recommended for re-commissioning and for re-qualification of the performance of the equipment.
- 2.62 Tests should only be undertaken after completion of the planned maintenance tasks described in Chapter 4, 'Operational management'.

- 2.63 Each test is cross-referenced to a detailed description of the test procedure in paragraph 2.72ff, 'Test methods'. Unless otherwise specified in 'Test methods', the tests should be carried out with the machine at normal working temperature, which might require a warm-up run to be carried out before commencement of testing.
- 2.64 A number of the tests recommended can be carried out concurrently on the same operating cycle and this is also indicated in 'Test methods'.
- 2.65 The results of periodic tests should be filed securely, for example in the plant history file.

Weekly safety checks

- 2.66 The CP(D) should make the following safety checks before starting the sequence of weekly tests:
- examine the door seal(s);
 - check the security and performance of door safety devices.
- 2.67 For equipment that includes a pressure vessel or pressure system (for example steam or compressed air) the following checks should be made:
- that safety valves or other pressure limiting devices are free to operate;
 - any other checks required by the CP(PS) in connection with the written scheme of examination for the pressure vessel.

Yearly safety checks

- 2.68 In order to ensure the continued safe functioning of decontamination equipment the CP(D) should conduct a series of safety checks before starting the yearly tests.
- 2.69 The AP(D) should draw up a documented programme of the yearly safety checks necessary for a particular installation with input, if required, from the AE(D). Reference should be made to the manufacturer's maintenance schedules.
- 2.70 The original installation checks and tests can be used as a basis for the yearly safety checks paying particular attention to those factors that affect safety and especially to those that might have changed since the previous annual safety check (or installation test).
- 2.71 The adequacy and safe connection of all engineering services should be verified.

Table 3 Schedule of periodic tests

Daily tests – User
1. Check spray arm rotation for free movement
2. Check spray nozzles for blockage (paying particular attention to those fitted to carriages for cannulated instruments)
3. Remove and clean strainers and filters etc
4. Ensure sufficient additives available and that dosing system is functioning
Weekly tests – User or CP(D)
1. Weekly safety checks
2. Carry out daily tests
3. Water hardness (all process stages)
4. Water conductivity (final rinse stage)
5. Automatic control test
6. Cleaning efficacy test by residual soil detection
Quarterly tests – CP(D)
1. Weekly safety checks
2. Automatic control test
3. Verification of calibration of WD instruments
4. Thermometric test for thermal disinfection
5. Cleaning efficacy test: – reference load ¹ general instruments endoscopic/MAT instruments – test soil
Yearly and revalidation tests – CP(D)
1. Yearly safety checks
2. Automatic control test
3. Verification of calibration of WD instruments
4. Water system: – chemical purity – bacterial endotoxins
5. Drainage: – free draining – efficacy of discharge
6. Doors and door interlocks: – Cycle start interlock – In-cycle interlock – Failed cycle interlock
7. Fault indication on sensor failure
8. Water vapour discharge test
9. Chemical additive dosing tests: – reproducibility of volume admitted – low level detection
10. Load carriers
11. Air quality

12. Cleaning efficacy test: – test soil – reference load ¹ general instruments endoscopic/MAT instruments
13. Over temperature cut-out test
14. Thermometric test for thermal disinfection – reference load ¹
15. Load dryness test – reference load ¹
16. Process residue test

Notes

¹ Additional test loads and alternative test soils may be required for WDs that are also intended for use with hollowware and/or anaesthetic accessories. The additional testing should also include tests on the load carriers that will be used with these additional loads.

Calibration, limits and function (including fault/alarm) of independent monitoring system should be checked during quarterly and yearly tests.

Test methods

2.72 Test methods, equipment and requirements are detailed in BS EN ISO 15883 parts 1 and 2. The following information, where relevant, will apply to those washer-disinfectors not covered by this Standard. For thermal (moist heat) disinfection the disinfection conditions are specified either by an A0 value to be achieved throughout the load or by a disinfection temperature band, defined by a minimum acceptable temperature, known as the disinfection temperature, and a maximum allowable temperature. Time/temperature bands meeting the requirements of an acceptable A0 of 600 are as follows:

Disinfection temperature (°C)	Maximum allowable temperature (°C)	Exposure time period
70	75	100 minutes
80	85	10 minutes
90	95	1 minute

2.73 Other time/temperature relationships may be used as long as achievement of the acceptable A0 has been validated to be achieved.

2.74 The higher the disinfection temperature, the shorter the holding time that should be used to achieve the same level of disinfection.

Drainage

Drain seal integrity

- 2.75 When it is impractical to vent the washer-disinfector externally a condenser can be used to allow venting into the workspace without discharging hot, humid air. The restricted flow associated with this system can produce a back pressure in the chamber. If the back pressure is excessive the water seal between the chamber and drain might be broken.
- 2.76 This test is designed to establish that the seal integrity is maintained under normal operating conditions. The test is intended for use both as a type test (and as such is a requirement of BS EN ISO 15883-1) and as an operational or installation test.

Apparatus

- 2.77 The following equipment should be used:
- a full load of the type the washer-disinfector is designed to process;
 - a test trap, of the same type and dimensions as normally fitted, but manufactured from a transparent material (type test only).

Method for type test

- 2.78 Fit the transparent trap in place of the normal trap and connect it to a suitable outlet. Pour sufficient water into the chamber to charge the trap to the normal level. Verify that there are no leaks.
- 2.79 Place a full load in the chamber, close the door and initiate an operating cycle. Without opening the door between cycles run a further four cycles. During each operating cycle observe the trap and establish whether the water seal in the trap has been broken.

Method for operational test

- 2.80 Carry out the test on the installed washer-disinfector with all services connected. Verify that the trap is charged with water to the normal working level.
- 2.81 Place a full load in the chamber, close the door and initiate the operating cycle. At the end of the cycle, remove the load and examine the water level in the trap. This may be done either visually or using a dipstick as appropriate. Without any delay (which would allow the chamber or load to cool), reload the machine with the same full load and run a further cycle. Repeat the same procedure until five

consecutive cycles have been run. After each operating cycle, observe the trap and establish whether the water seal in the trap has been broken.

Results

- 2.82 The water seal should not be broken during the test.

Free draining (tanks, chamber, load carriers, pipework)

- 2.83 Residual water that does not drain from the internal pipework of the washer-disinfector can provide an environment for microbial growth; these microorganisms might then be available to recontaminate the disinfected load.
- 2.84 The following checks should be carried out during type-testing, works testing and installation testing to verify that as designed, built and installed the washer-disinfector will effectively discharge all the water from the system.

Method

- 2.85 Test the free draining of chamber and load carriers by visual observation at end of the cycle.
- 2.86 Test the free draining of tanks by visual observation on draining the tanks.
- 2.87 Test the pipework flow to discharge point by visual observation including use, when necessary, of a spirit level.

Purging of the trap (type test only): efficacy of discharge through the trap

- 2.88 The test is intended to verify that the operating cycle is effective in purging the trap of all waste and soil.
- 2.89 The test can be carried out as part of the cleaning efficacy test during operational testing.

Apparatus

- 2.90 The following apparatus should be used:
- test soil appropriate to the type of washer-disinfector being tested (see [paragraph 2.268](#), 'Periodic tests');
 - sampling tube of sufficient length to reach the water trap in the drain of the washer-disinfector and a sampling pump (for example a pipette pump or syringe).

Method

- 2.91 On completion of an operating cycle to test the cleaning efficacy by processing a full load contaminated with an appropriate test soil, place the sampling tube into the water trap and remove a sample for examination.
- 2.92 Examine the water sample from the trap for residual test soil using the detection method appropriate to the test soil.

Results

- 2.93 The water in the trap should be free from residual soil at the same level of detection as that specified for the load items. Residual soil in the trap can present an infection or recontamination hazard.

Estimation of dead volume of pipework

- 2.94 Residual water that does not drain from the internal pipework of the washer-disinfector might provide an environment for microbial growth; these microorganisms might then be available to re-contaminate the disinfected load.
- 2.95 This test is intended primarily as a type test but might also be of value as an operational test or when investigating microbial contamination occurring in a washer-disinfector.
- 2.96 The test should only be carried out once the checks for free drainage recommended in paragraph 2.83, 'Free draining (tanks, chamber, load carriers, pipework)' have been satisfactorily completed.

Apparatus

- 2.97 Volumetric measuring vessels of appropriate size should be used.

Method

- 2.98 Ensure that the pipework of the washer-disinfector is dry (either following disassembly and re-assembly or purging with dry compressed air for not less than 30 min), then flush with a known volume of water (simulating the flow that would occur in normal use).
- 2.99 Measure the volume of water discharged and the dead volume, estimated as the volume retained, calculated from the difference between the two values.
- 2.100 When the washer-disinfector has two or more pipework systems that are entirely separate (e.g. for flushing water, wash water, rinse water or chemical

disinfectant solution) each system may be tested separately.

Results

- 2.101 The volume of retained water should be less than 1% of the volume of water used.

Venting system*Steam venting*

- 2.102 For washer-disinfectors in which the load is heated/thermally disinfected by steam heating the chamber should be protected against a rise in pressure above the designed working pressure of the chamber and any discharge of steam should be only from the vent provided.

Apparatus

- 2.103 A pressure gauge should be used (if the washer-disinfector chamber is not already equipped with a pressure gauge in a known state of calibration).

Method

- 2.104 Close and seal the chamber of the washer-disinfector in the manner specified by the manufacturer and start an operating cycle. Override the automatic controller to allow the continuous admission of steam to the chamber. Observe where steam is vented. Record the maximum value obtained on the pressure gauge.

Results

- 2.105 Steam discharge should occur solely through the vent. The recorded value for the pressure within the chamber should not exceed the specified working pressure or, for a chamber designed to work at atmospheric pressure, 200 mbar.

Load contamination from ductwork

- 2.106 The evolution of water vapour from the chamber during the washing stage, thermal disinfection stage and drying stage can result in condensation occurring in the ductwork and in the condenser (if fitted). The ducting is commonly arranged to allow this condensate to drain back into the chamber. Should this condensate become contaminated there is a risk that it could contaminate the load. The test is designed to establish that any condensate draining back into the chamber will not contact the load. The test is intended for use both as a type test (and as such is a requirement of BS EN ISO 15883-2) and as an installation or operational test.

Apparatus

- 2.107 The following apparatus should be used:
- vessel of not less than 500 mL capacity having a discharge port at its base connected to a flexible tube fitted with an on/off valve and a flow control valve;
 - stopwatch;
 - load carrier and full load for the washer-disinfector;
 - paper towels.

Method

- 2.108 Disconnect the external ducting to the washer-disinfector 1 m above the chamber.

Note

If it is not possible to disconnect the ducting at this position the ducting should be disconnected at the chamber and a spare 1 m length of ducting should be connected to the chamber.

- 2.109 Position the vessel approximately 1 m above the level of the chamber discharge to the vent.
- 2.110 With the on/off valve closed, fill the vessel with 200 mL \pm 20 mL of cold water. Open the valve and adjust the flow control valve so that the contents of the vessel are discharged in 1 min \pm 5 s.
- 2.111 Refill the vessel with 200 mL \pm 20 mL of cold water. Feed the flexible tube into the ducting so that the open end of the flexible tube is 600–800 mm above the top of the chamber.
- 2.112 Load the chamber with a full load of dry load items in accordance with the manufacturer's instructions. Close the chamber door and then open the on/off valve. Record the time required for the vessel to empty.
- 2.113 Within ONE minute of the vessel emptying open the chamber door and remove the load and any removable load containers. Place all the removed items on absorbent paper and examine all surfaces of the load and the absorbent paper for traces of water.
- 2.114 Repeat the above procedure for the full range of load carriers that the washer-disinfector is designed to process.

Results

- 2.115 There should be no visible water on the load or load carriers.

Overflow test

- 2.116 For washer-disinfectors that incorporate one or more water storage tanks within the washer-disinfector, the capacity of the overflow(s) to discharge all excess water, as intended, without spillage into the washer-disinfector or working area should be verified.

Type/works test method

- 2.117 Ensure that the washer-disinfector is connected to all necessary services and the water supply pressure adjusted to not less than 6 barg under the conditions of flow that prevail with the supply valve(s) fully open.
- 2.118 Fully open the supply valve(s).
- 2.119 Observe the level of water in each tank or cistern until this has been unchanged for not less than 2 min.

Installation test

- 2.120 The washer-disinfector should be connected to all necessary services.
- 2.121 Fully open the supply valve(s).
- 2.122 Observe the level of water in each tank or cistern until this has been unchanged for not less than 2 min.

Results

- 2.123 The washer-disinfector and installation should be regarded as satisfactory when equilibrium conditions have been attained within the tank(s) without discharge of water other than by the intended (piped) overflow.

Volume of water used per stage

- 2.124 During type-testing, the manufacturer should be required to determine the volume of water used during each stage of the cycle. These data are used in calculations of the service requirements.
- 2.125 In addition, during installation or operational testing, the volume of water used for each stage of should be verified. If the volume of water used is insufficient the efficacy of the cleaning and disinfection process might be adversely affected. If the volume is greater than that specified an

unexpectedly heavy demand might be placed upon the water supply.

- 2.126 There are three methods that may be used for determining the volume of water used. The method should be chosen on the basis of which is most convenient for the particular installation.

Apparatus

- 2.127 A water flow meter (or volumetric measuring equipment) should be used.

Method

- 2.128 Fit a water flow meter in each of the water supply pipes, consecutively or concurrently, and determine the volume used by comparison of the reading before and after each stage of the process cycle. Operate the washer-disinfector with the chamber empty, apart from the chamber furniture. Follow the water meter manufacturer's instructions for installation. Pay particular attention to the length of uninterrupted straight pipe required on either side of the meter.
- 2.129 When the washer-disinfector is supplied from a readily accessible tanked supply interrupt the make-up to the tank and mark the water level. Determine the volume of water required to restore the level after an operating cycle stage by the addition of a measured volume of water.
- 2.130 For those washer-disinfectors which discharge all the water from each stage at the end of each stage, obtain a suitable estimate of the volume used by volumetric measurement of the discharge from the drain.

Results

- 2.131 The volume of water used for each stage of the cycle should be within $\pm 5\%$ of the volume specified by the manufacturer.

Doors and door interlocks

Cycle start interlock

- 2.132 The interlock should prevent a cycle being started with the door open.

Method

- 2.133 Leave the doors open and unlocked. Ensure that all services are connected. Make an attempt to initiate an operating cycle.
- 2.134 Close and lock the doors and make a further attempt to initiate an operating cycle.

Results

- 2.135 It should not be possible to initiate a cycle with the door(s) left open. With the door(s) closed it should be possible to initiate an operating cycle.

In-cycle interlock

- 2.136 An interlock is required to ensure that the door(s) cannot be deliberately or inadvertently opened while the washer-disinfector is in operation.

Method

- 2.137 Close and lock the door(s) and start the operating cycle. While the operating cycle is in progress attempt to unlock each of the doors. Where practicable, visually inspect the interlocks to verify engagement before attempting to open the door.

Results

- 2.138 In these circumstances it should not be possible to unlock any of the doors.

Double-door washer-disinfectors

Method

- 2.139 Both during and between cycles, make attempts to open either or both the loading door and unloading door of the double-door washer-disinfector.

Results

- 2.140 It should not be possible to open the unloading door after initiation of a cycle until a cycle has been completed satisfactorily.
- 2.142 It should not be possible for both doors to be opened at the same time.
- 2.143 It should not be possible to open the loading door until a cycle has been satisfactorily completed and the unloading door has been opened and closed.

On sensor failure

Method

- 2.144 Disable each sensor in turn and attempt to open each of the door(s). Where practicable, avoid the undertaking of checks during an operating cycle.

Results

- 2.145 In each case it should not be possible to open the door(s).

Door opening force

2.146 The mechanism for opening the washer-disinfector door should not require the use of excessive force.

Apparatus

2.147 The following equipment should be used:

- a. spring balance calibrated in kilograms with a range including 0–250 kg and with an accuracy of ± 1 kg over the range 0–250 kg;
- b. non-extensible means of attachment of the spring balance to the door mechanism.

Method

2.148 Measure the force required to initiate and sustain the movement of the door opening mechanism by interposing a spring balance, aligned co-axially with the direction of movement of the door opening mechanism, between the operator and the mechanism.

2.149 Attach the spring balance to the door opening mechanism. Open the door, record the force required to initiate the movement and to sustain the movement.

Results

2.150 The indicated value required to initiate or sustain the movement of the door opening mechanism should not exceed 25 kg.

Failed cycle interlock

2.151 The interlock should prevent the Operator from removing a load in the normal manner at the end of a cycle that failed.

Method

2.152 During an operating cycle interrupt one, or more, of the services to the washer-disinfector sufficiently to cause a cycle failure.

Results

2.153 A fault should be indicated. It should not be possible to open the unloading door (if fitted); it should only be possible to open the loading and/or unloading door by means of a special key, code or tool.

Fault indication on sensor failure

2.154 A failure of any sensor used as part of the control system of the washer-disinfector should cause a fault to be indicated by the automatic controller.

Method

2.155 Disable each sensor providing information to the automatic controller in turn to establish that a fault is indicated.

2.156 Test each sensor as follows. Start an operating cycle. During, or before, the stage of the cycle at which the sensor is intended to provide data used to determine the control of the cycle, disable the sensor.

2.157 Test each sensor in both “open circuit” and “short circuit” failure modes.

Result

2.158 A fault should be indicated during or at the end of the cycle. It should not be possible to open the door on a single-ended washer-disinfector or the unloading door of a double-door washer-disinfector.

Chemical dosing

Reproducibility of volume admitted

2.159 This test is intended to verify the setting for the dispensed volume of chemical additive(s) and to ensure that it is reproducible within defined limits. The test should be carried out for each chemical dosing system on the washer-disinfector.

Apparatus

2.160 Two measuring cylinders that conform to BS EN ISO 4788 should be used. The size of measuring cylinder should be appropriate to the volume of chemical additive to be dispensed.

Note

When compatibility with the chemical additive to be measured has been established, conformity to BS 5404-2 is sufficient.

Method

2.161 Disconnect the supply line to the chamber as close as possible to its discharge point into the chamber or water circulation system in order to discharge into the measuring cylinder.

2.162 Place detergent in the first cylinder to a known mark. Place detergent to a known mark in the second cylinder.

2.163 Actuate a normal cycle and, at the end of the dosing stage, top up the first cylinder to the

original mark from the second cylinder. Calculate the detergent added from the second cylinder.

- 2.164 Repeat the test three more times; record the volume added on each test.
- 2.165 Care should be taken since many of the concentrates used are irritant or corrosive. Water might not be an acceptable substitute because, for many dosing systems, differences in viscosity can affect the dispensed volume.

Results

- 2.166 The result of the first test should be ignored.
- 2.167 The mean collected volume from the final three tests should be within $\pm 10\%$ of the nominal dispensed volume.

Indication of insufficient chemical additives

- 2.168 The correct volume of chemical additive(s) for the correct functioning of the washer-disinfector should be used. The washer-disinfector should be equipped with means to ensure that a cycle is not initiated when there is insufficient chemical additive remaining in the reservoir to complete a cycle.
- 2.169 The test should be carried out for each chemical dosing system on the washer-disinfector.

Method

- 2.170 Place a low level of additives in the dispenser reservoir and run repeated cycles.
- 2.171 Fill an otherwise empty container with sufficient chemical for more than two but less than four operational cycles. Run the washer-disinfector on three consecutive cycles. Estimate the volume remaining at the end of each cycle (pre-marked container, dipstick, or weight).

Results

- 2.172 The washer-disinfector should indicate at the beginning of the third or fourth cycle that there is insufficient chemical remaining to complete a cycle.

Water vapour emissions

- 2.173 Faulty or damaged door seals, or faulty condensers, can give rise to vapour emission into the working area and the leakage of potentially infectious material from the washer-disinfector.

- 2.174 Excessive and persistent leakage carries the risk of scalding the operator and causing deterioration of walls and their surface finishes.

Apparatus

- 2.175 The following equipment should be used:
- absorbent paper wipes (of a type which change colour density when damp);
 - one or more mirrors 50 mm × 50 mm or larger.

Method

- 2.176 Load the washer-disinfector, close the door and wipe the joints between the door and the door surround to remove any moisture. Carry out an operating cycle.
- 2.177 Throughout the operating cycle use the mirror(s) to check if water vapour escapes from the door seal or from the condenser (if fitted).
- 2.178 At the end of the operating cycle, with the door still closed, use the absorbent wipes to wipe the joints between the door and the door surround as close as possible to the door seal. Examine the wipes for dampness.
- 2.179 A further four operating cycles should be run with the checks described above being carried out on the final cycle.

Results

- 2.180 There should be no misting of the mirror(s), which would be evidence of vapour emission, and no dampness of the absorbent wipes, which would be evidence of vapour or liquid emission.

Instrumentation fitted to a washer-disinfector

Verification of calibration

- 2.181 The calibration of instrumentation and any independent monitor fitted to the washer-disinfector should be verified by comparison with calibrated test instruments during steady state conditions, for example the temperature during the disinfection hold period. Compliance to BS EN ISO 15883-1 clause 5.12–17 should be met.
- 2.182 This may be carried out concurrently with other testing, for example during the automatic control test during quarterly periodic testing (see [paragraph 2.304, 'Automatic control test'](#)).

Load carriers

2.183 Load carriers come in a variety of forms including trolleys, carriages and baskets. Their correct functioning is essential to the successful outcome of a washer-disinfector operating cycle. It is important that they cannot easily be misaligned, that they function correctly and that, when applicable, they make good connection with service supply points in the chamber and with load items (when necessary).

Method

2.184 Verify by visual observation the alignment of load carriers, their connection to water, air or chemical additive supply in the chamber (when applicable) and their connection to load items, for example cannulated instruments (when applicable).

2.185 Check load carriers with rotary spray arms to ensure that the spray arms are free to rotate, both when the load carrier is empty and when fully loaded.

Thermometric tests

2.186 Thermometric tests are carried out to verify the attainment of the specified conditions throughout the chamber and load during the operating cycle. Continuous process washer-disinfectors and multi-chamber washer-disinfectors in which the use of recorders with fixed sensors is impractical should be tested using single channel data loggers that can be processed through the washer-disinfector. Biological indicators should not be used as a substitute for thermometric testing.

Chamber wall temperature

Apparatus

2.187 A temperature recorder, according to the recommendations given in the 'Decontamination equipment: test equipment and materials' section of CFPP 01-01 Part B, with not less than 12 sensors should be used.

Method

2.188 Locate thermocouples as follows: one in each corner of the chamber, one in the centre of the two side walls, one in the centre of the roof of the chamber and one adjacent to the temperature sensor used as the reference sensor for chamber temperature.

2.189 Measure the temperature attained throughout four operating cycles, the first of which should be at

least 60 min since the machine was last used (a cold start) and the final three with not more than a 15 min interval between cycles (a hot start).

2.190 The washer-disinfector should be operated empty except for chamber furniture (for example load carriers).

2.191 Multi-chamber washer-disinfectors may be tested with each chamber tested consecutively or concurrently. In the latter case eight sensors should be used for each chamber.

Results

2.192 The results should be as follows.

- a. The temperatures recorded on the surface of the chamber should be within the range 0–5°C of the disinfection temperature throughout the holding period for the disinfection stage.
- b. The temperatures recorded on the surface of the chamber should be within ±5°C of the set temperature for the relevant stage throughout the holding period for each of the other stages.
- c. The temperature indicated/recorded by the washer-disinfector instruments should be within ±2°C of that recorded by the test instrument from the sensor adjacent to the reference sensor throughout the holding period for the disinfection stage.
- d. The temperature profile obtained for the operating cycle should be consistent within ±2°C for the last three test cycles.

Load carrier temperature

Apparatus

2.193 A temperature recorder (according to the recommendations given in the 'Decontamination equipment: test equipment and materials' section of CFPP 01-01 Part B) with not fewer than four sensors should be used.

2.194 Three independent data loggers and a temperature recorder having at least one sensor may be used as an alternative.

Method

2.195 Locate temperature sensors at two diagonally opposite corners of the load carrier, in the approximate geometric centre of the load carrier and adjacent to the temperature sensor used as the reference sensor for chamber temperature.

- 2.196 Measure the temperature attained throughout four operating cycles, the first of which should be at least 60 minutes since the machine was last used (a cold start) and the final three with not more than a 15-minute interval between cycles (a hot start). Ensure that the washer-disinfector is empty except for chamber furniture (for example load carriers).
- 2.197 Replace the load carrier between cycles with a load carrier at ambient temperature.
- 2.198 Test each chamber of multi-chamber washer-disinfectors consecutively using independent data loggers to record the temperature of the load carrier. Use a temperature recorder with fixed sensors to record the temperature adjacent to the reference sensor.
- 2.199 Keep the washer-disinfector in continuous operation, with not more than 15 minutes between cycles, when the length of cycle/the number of dataloggers available precludes re-use of the dataloggers, so that when the second and subsequent tests are initiated not more than 15 minutes has elapsed since the first chamber completed a cycle.
- 2.200 This test may be run simultaneously with the chamber wall temperature test.

Results

- 2.201 The results should be as follows.
- The temperatures recorded on the surface of the load carrier should be within the range 0–5°C of the disinfection temperature throughout the holding period for the disinfection stage.
 - The temperatures recorded on the surface of the load carrier should be within $\pm 5^\circ\text{C}$ of the set temperature for the relevant stage throughout the holding period for each of the other stages.
 - The temperature indicated/recorded by the washer-disinfector instruments should be within $\pm 2^\circ\text{C}$ of that recorded by the test instrument from the sensor adjacent to the reference sensor throughout the holding period for the disinfection stage.
 - The temperature profile obtained for the operating cycle should be consistent within $\pm 2^\circ\text{C}$ for the last three test cycles.

Over-temperature cut-out

- 2.202 The washer-disinfector is fitted with an over-temperature cut-out to ensure that, in the event of the automatic control failing to control the temperature in the washer-disinfector, the temperature will not rise to a level that would damage the load in the washer-disinfector.

Apparatus

- 2.203 A temperature recorder, according to the recommendations given in the 'Decontamination equipment: test equipment and materials' section of CFPP 01-01 Part B, with not less than four sensors should be used.
- 2.204 Three independent data loggers and a temperature recorder having at least one sensor may be used as an alternative.

Method

- 2.205 Locate temperature sensors at two diagonally opposite corners of the load carrier, in the approximate geometric centre of the load carrier and adjacent to the temperature sensor used as the reference sensor for chamber temperature.
- 2.206 Operate the washer-disinfector on a normal operating cycle, empty except for the load carrier. For multi-cycle machines test the two cycles that have the highest and lowest operating temperatures.
- 2.207 During the stage of the cycle when the maximum temperature is attained, disable the temperature control system.

Results

- 2.208 The over-temperature cut-out should operate at a temperature not more than 5°C higher than that provided by any temperature control or temperature limiting device.

Load dryness

- 2.209 If the washer-disinfector includes a drying stage, the drying efficacy should be tested on a test load as follows:
- Within five minutes of the end of a normal cycle initiated from cold, place the load on a sheet of coloured crepe paper.
 - Observe any water emanating from the load and carriage, and examine the crepe paper from any residual water from the load staining it.

- Lumen instruments should be examined by blowing dry compressed air through the lumen onto a mirror surface.

Results

2.210 No residual water should be observed from the load or carriage, on the crepe paper or, where relevant, on the mirror surface.

Residual chemical additives

- 2.211 The nature of the residues and the level of such residues that might be of concern depend on the chemical additives used during the process and the intended use of the washed and disinfected product.
- 2.212 The chemical additives used during the process (detergents, rinse aids etc) might not be completely removed by the rinsing process.
- 2.213 The residual level that should be tolerated depends on the nature of the chemical and the intended use of the product. The supplier of any chemical agent used will normally provide data on the chemical composition of the chemical agent and the biocompatibility of the components of the chemical agent. The supplier will also normally provide details of the method of detection that may be used to determine that processed items are free from residuals at the specified levels.
- 2.214 The sampling method and analytical method should be capable of determining the presence of the chemical additive at concentrations below that specified as potentially harmful, i.e. as the maximum acceptable level.

Method

2.215 Test the efficacy of the rinse process by using twice the normal dose of the chemical additive on a normal operating cycle using a test load of the simulated product. Analyse the final rinse water and the simulated product using the method recommended by the manufacturer.

Results

2.216 The concentration on the simulated product should be lower than the specified maximum acceptable level.

Air quality

2.217 Many washer-disinfectors are fitted with air filters to remove particulate material from the air supplied to the drying stage. These filters are often

HEPA filters (for example, EU 12/13) of the type used to remove bacterial contamination from the air supply. When they are used as general particulate filters, performance tests for the filter or the filter housing are not necessary, except when the intention is to provide air free from bacterial contamination when the load is intended for use without further processing (for example sterilization).

2.218 Microbial sampling is not necessary for either system unless otherwise specified.

Method

2.219 The complete installation should be tested using the method specified in BS EN ISO 14644-1. Additional reference should be made to BS EN ISO 14644-1 for test procedures. A challenge aerosol of inert particles of the type produced by a dispersed oil particle generator should be introduced into the air upstream of the filter. The downstream face of the filter and its housing should then be scanned for leakage using a photometer.

Results

2.220 The reading on the photometer should be steady and repeatable and should not exceed 0.01% of the upstream reading.

Sound pressure test

2.221 The British Standard requires the manufacturer to carry out a sound power test as a type test for decontamination equipment. This test measures the total sound power radiated from the machine and should be performed in a specially designed and equipped test room. It is neither necessary nor practicable to repeat the test on an installed machine.

2.222 The perceived level of noise in the immediate vicinity of the equipment during operation is, however, of concern. The perceived noise level depends not only upon the sound power level of the equipment but also on the acoustic properties of the environment and other sources of noise. It should necessarily be determined with the washer-disinfector installed and working normally.

2.223 A failure of the sound pressure test need not be an indication that the machine is faulty. The problem might lie in the acoustic properties of the room in which the machine is installed.

2.224 The sound pressure test should be carried out in accordance with BS EN ISO 3746 by suitably trained and experienced personnel. The guidelines given in the following two paragraphs are intended only for additional guidance and are not the complete test method.

Method

2.225 Use the procedure specified in BS EN ISO 3746 for both the loading and unloading area if these are not common (and plantroom if present), to determine the following:

- a. the mean A-weighted surface sound pressure level;
- b. the peak A-weighted surface sound pressure level.

Results

2.226 The test should be considered satisfactory if the following recommendations are followed:

- a. the mean A-weighted surface sound pressure level does not exceed:
 - (i) 55 dBA for decontamination equipment installed in an operating suite, ward, treatment room or other noise sensitive area;
 - (ii) 70 dBA for decontamination equipment installed in a sterile services department, laboratory or pharmacy production area;
- b. in both the loading and unloading area the peak A-weighted surface sound pressure does not exceed the mean A-weighted surface sound pressure level by more than 15 dBA.

Electromagnetic compatibility

2.227 The British Standard for washer-disinfectors specifies that:

- when tested by one of the methods in BS EN 61000-4-3, the functioning of the automatic controller and the instrumentation should be unaffected by electromagnetic interference of severity level 3 as specified in BS EN 61000-4-3;
- when tested in accordance with BS EN 55014-1, any RF interference generated by the washer-disinfector should not exceed the limits specified in BS EN 55014-1.

2.228 The current proposal within the CEN TC102 committee preparing standards for washer-

disinfectors is to specify compliance with BS EN 61000-6-3 and BS EN 61000-6-1.

2.229 These tests are applicable as type tests or works tests only.

2.230 Since June 1993, washer-disinfectors classified as medical devices have been required to carry the CE marking. Under the Medical Devices Directive, any washer-disinfector which bears CE marking solely under the Electromagnetic Compatibility Directive should have the CE marking identified as being applied under the EMC Directive.

Cleaning efficacy tests

2.231 Cleaning efficacy tests are used to demonstrate the ability of the washer-disinfector to remove or reduce to acceptable levels, soiling and contamination that occurs during normal use of reusable items.

2.232 Test soils are used to simulate naturally occurring contamination since the latter shows considerable variation and is therefore more difficult to use for standardized testing.

Type tests

2.233 The cleaning efficacy should be determined using the relevant test soil applied to a reference load or simulated product of demonstrated relevance.

2.234 The manufacturer will normally establish worst-case conditions of temperature, detergent concentration, water hardness and water pressure/flow rate for use during testing.

2.235 By analysing the fraction of soil removed during the cleaning process when operated for various time periods shorter than those that will normally be used a quantitative comparison of cleaning efficacy can be made.

2.236 The recommended minimum operating conditions given by the manufacturer should be based on these data that should be made available to the User.

Operational tests

2.237 During operational tests of cleaning efficacy with test soils the thermal disinfection stage should be disabled.

2.238 During thermal disinfection the action of hot water/steam can also reduce the concentration of residual test soil.

- 2.239 The drying stage may also be disabled if this is necessary to facilitate the detection of residual soil.

Test soil

- 2.240 The choice of test soil to be used should be based on the intended use of the washer-disinfector and should be formulated to simulate the soiling which will be encountered in practice and which would be most difficult to remove.

Test loads

- 2.241 The test load should consist of items of similar size, mass and materials of construction to the range of products the washer-disinfector is intended to process. Care is needed if loads are mixed or lacking in uniformity.

Method for chamber walls and load carriers

- 2.242 Contaminate the chamber walls and load carrier with the test soil in accordance with the manufacturer's instructions for the test soil including the specified quantities to be used and any drying stage.
- 2.243 Run a normal operating cycle.
- 2.244 After completion of the wash cycle, and before the disinfection stage (except where this is combined with the rinse stage) abort the cycle.
- 2.245 For operational tests, carry out the test in duplicate for each type of operating cycle available on the washer-disinfector.
- 2.246 When used as a periodic test, carry out the test only once for each type of operating cycle available.

Method for reference loads

- 2.247 This test should be run only after satisfactory completion of the test for the efficacy of soil removal from chamber walls and load carriers.
- 2.248 Contaminate the test load with the test soil in accordance with the manufacturer's instructions for the test soil.
- 2.249 The specified quantities should be used and any drying stage should be carried out in strict accordance with the instructions.
- 2.250 Run a normal operating cycle for the load type under test.
- 2.251 Abort the cycle after completion of the wash cycle, and before the disinfection stage (except where this is combined with the rinse stage). Examine the test

load, chamber walls and load carrier for the presence of residual soil.

- 2.252 For operational tests, carry out the test in duplicate for each type of operating cycle available.
- 2.253 When used as a periodic test, carry out the test only once for each type of operating cycle available.

Results

- 2.254 The chamber walls and load carrier should be free from the test soil to the extent specified for the test soil employed.
- 2.255 The test load should be free from the test soil to the extent specified for the test soil employed and no test soil should have been transferred to the chamber walls or load carrier.

PQ tests

- 2.256 PQ tests of cleaning efficacy are necessary only when some of the items, or some of the loads, to be processed are more difficult to clean than the reference load.

Method

- 2.257 Repeat the tests described above for reference loads with actual loads to be processed specified by the user as being representative of the items or loads intended to be processed.

Note

The items to be processed might need to be replaced by surrogate devices when the design of the actual item makes subsequent examination for residual soil impractical.

- 2.258 Contaminate the test load with the test soil in accordance with the instructions for the test soil.
- 2.259 The specified quantities should be used and any drying stage should be carried out in strict accordance with the instructions.
- 2.260 Run a normal operating cycle for the load type under test.
- 2.261 Abort the cycle after completion of the wash cycle, and before the disinfection stage (except where this is combined with the rinse stage). Remove the test load and examine for the presence of residual soil. Run these tests in duplicate.
- 2.262 On satisfactory completion of this part of the test, run a further three cycles with actual loads, of the

type intended to be processed, contaminated with natural soiling in-use.

- 2.263 Run a normal operating cycle for the load type under test.
- 2.264 After completion of the complete operating cycle remove the test load and examine for cleanliness.
- 2.265 Assess the cleanliness of the processed items visually or by such other means as will be routinely used for acceptance testing (see paragraph 2.268, 'Tests for residual soil').

Results

- 2.266 The test loads should be free from the test soil to the extent specified for the test soil employed and no test soil should have been transferred to the chamber walls or load carrier.
- 2.267 The cleanliness of the processed items should be acceptable by the means that will be used routinely for acceptance testing as per paragraph 2.268, 'Tests for residual soil'.

Periodic tests

Tests for residual soil

Note

The development path for this guidance recognises the current work in the piloting of high sensitivity post-decontamination protein quantification tests. As these become validated and available, this guidance will be amended to incorporate their use.

Apparatus

2.268 The following apparatus should be used:

- sterile purified water;
 - ninhydrin reagent (0.30 g ninhydrin in 100 mL 70% v/v isopropanol);
 - cotton swabs (plastic handles not wooden); (Alternatively a small quantity of cotton wool held in the jaws of clean forceps may be used.)
 - glass microscope slide 25 mm × 75 mm;
 - oven set at 110°C ± 2°C;
 - 0.5 g/L arginine solution.
- 2.269 Alternative systems/kits may be commercially available. Use of these should demonstrate resolution and accuracy at a similar level to that of the technique described above.

- 2.270 The swabs used for sampling should not have a wooden handle as the wood can interfere with the reaction giving a false positive.

Ninhydrin method

- 2.271 The ninhydrin method uses the reaction of amino acids, peptides and proteins, with 1,2,3-indantrione monohydrate (Chemical Abstracts Service reference number 485-47-2).
- 2.272 This method detects a broad spectrum of substances that can occur as residuals from body fluids and from microbiological studies and is thus suitable for detecting residuals in many of the applications for which washer-disinfectors are employed. Because only the cotton swab moistened with sterile distilled water comes into contact with the item being tested the item may be used after testing if necessary.
- 2.273 Swab a defined area of the surface to be examined with a sterile cotton swab moistened with sterile water. The swab should be compatible with this procedure as many swabs contain sufficient protein for the result to be a positive result regardless of the condition of the instrument.
- 2.274 Air-dry the swab and then transfer it to the laboratory for examination or if suitable facilities exist on-site examine on site.
- 2.275 Test the swab for the presence of residual protein by adding a drop of ninhydrin reagent to the swab and heating the swab at 110°C for 30 min. After heating examine the swab for a purple colour.
- 2.276 Instruments with a lumen can be tested by using a pipe-cleaner style cleaning brush instead of the swab.
- 2.277 Each test, or series of tests, should be accompanied by a negative control. Moisten a cotton swab of the same batch with sterile water of the same batch and then immediately test using the ninhydrin reagent. The swab should be compatible with this test procedure.
- 2.278 A positive control can also be carried out to establish the sensitivity of the test method. Prepare a microscope slide by thorough cleaning, test using the method and demonstrate to be negative. Inoculate the surface with four 25 µL drops of a 0.5 g/L arginine solution, air dry and use as a positive control test surface.

Results

- 2.279 A purple discoloration on the swab indicates the presence of amino acid, peptide or protein residues.
- 2.280 The limit of detection for the method is equivalent to 2 mg m⁻² determined as arginine, or better.
- 2.281 There should be no purple discoloration on clean items.

Biuret test method

- 2.282 Commonly used as a test for the presence of proteins and detects the presence of two or more adjacent peptide bonds.
- 2.283 Flush the reagent solution through the channels of instruments with a long narrow lumen to detect residuals within the instrument. After testing the item clean to remove any residual reagent before it is used.

Results

- 2.284 A violet colour appears when a protein or tripeptide is treated with sodium hydroxide and copper (II) sulphate solutions.

Orthophthalaldehyde (OPA) method

- 2.285 Orthophthalaldehyde can be used, in the presence of a thiol compound, to detect a and e terminal amino groups of human blood proteins by formation of a fluorescent compound detectable at 360 nm in picomole quantities.
- 2.286 The method is only applicable where there is a suitably equipped laboratory available. The sensitivity may be greater than that provided by the ninhydrin method.

Visual inspection

- 2.287 High levels of residual soil, process chemical residues and residues from poor quality water can be detected by visual inspection of the items which have been cleaned.

Method

- 2.288 This should be done by unaided normal, or corrected, vision under suitable illumination (not less than 400 lux at the work surface).

Disinfection efficacy tests

- 2.289 Thermometric tests should be used for both thermal disinfection processes and chemical disinfection processes where temperature is a critical parameter. For thermal disinfection

processes the time temperature relationships are defined in BS EN ISO 15883 Parts 1 and 2.

Thermal disinfection test

- 2.290 During thermometric tests for thermal disinfection, the washing stages should be disabled or the controlled temperature reduced to ambient (20°C ± 5°C) in order to avoid pre-heating the load. Reducing the wash temperature to 20°C creates the worst-case conditions with which the disinfection stage might be expected to cope and ensures that disinfection conditions will be attained in the event of a failure of the washing stage.
- 2.291 Temperature monitoring of the load should be used to determine the attainment of the required time-temperature conditions.

Thermometric test for disinfection

- 2.292 This test is suitable for all washer-disinfectors and should be used to establish the adequacy of temperature control during chemical disinfection, as well as for verifying attainment of thermal disinfection conditions.
- 2.293 The load under test will consist of a reference load or a PQ load of discrete items of the type that the washer-disinfector under test is intended to process, or of surrogate devices used to simulate such load items.

Apparatus

- 2.294 The following equipment should be used:
- temperature recorder (see the 'Decontamination equipment: test equipment and materials' section in CFPP 01-01 Part B);
 - thermocouples or self-contained data loggers (see the 'Decontamination equipment: test equipment and materials' section in CFPP 01-01 Part B).
- 2.295 For type 1 machines and type 2 machines without physical separation of compartments (conveyor washer-disinfectors) sensors may be passed into the chamber through the thermocouple entry port into the chamber.
- 2.296 For type 2 (conveyor washer-disinfectors) the sensors to fixed positions may be passed into the chamber(s) through the thermocouple entry port; the sensors on the load might have to be fed in from one end and, when the load exits at the other end, detached and withdrawn from the back

through the washer-disinfector. Care should be taken to ensure that there is sufficient length of cable for this to happen.

- 2.297 When this is not possible the method used for type 2 (multiple cabinet) washer-disinfectors should be adopted.
- 2.298 For type 2 (multiple cabinet) washer-disinfectors the sensors to fixed positions may be passed into the chambers through the thermocouple entry port; the sensors on the load should be provided from self-contained data loggers within the load.

Method

2.299 Place temperature sensors in the following positions:

- at least one on an item at each level in the load carrier (up to a maximum of three if the load carrier accommodates load items on more than one level);
- one on an item in the region known to be slowest to attain the
- disinfection temperature*;
- one on an item in the region known to be fastest to attain the disinfection temperature*;
- one adjacent to the automatic control temperature sensor;
- one adjacent to the process recorder sensor (if fitted) in each chamber or compartment;
- three on the load carrier as follows: two at two diagonally opposite corners and one in the approximate geometric centre.

* These positions should be specified by the manufacturer and supported by data from type tests. If these data are not available from the manufacturer preliminary tests to map the temperature throughout the load will be necessary.

- 2.300 The sensors should be in good thermal contact with the item or installed sensor that they are monitoring and placed, if possible, in or on the part of the item that will be slowest to heat up.
- 2.301 The test should be performed in triplicate for PQ and commissioning tests but once for periodic testing.

Results

2.302 The test should be considered satisfactory if the following recommendations are followed:

- a. the recommendations of the automatic control test are followed;
- b. the holding time, as determined from the measured temperatures on the surface of the load items, is not less than that recommended to give an A_0 of 600 as defined in BS EN ISO 15883-1;
- c. during the holding time the measured temperatures are within the disinfection temperature band recommended for the operating cycle;
- d. the indicated and recorded chamber temperatures are within 2°C of the temperature measured at the automatic control sensor;
- e. the temperature measured on the surface of each load item does not fluctuate by more than $\pm 2^\circ\text{C}$ and does not differ from that in other load items by more than 4°C;
- f. at the end of the cycle: the temperature sensors have remained in position.

2.303 If having completed the commissioning tests based on a reference load the washer-disinfector fails to follow the above recommendations for the specific PQ load then it is possible that the washer-disinfector is not capable of processing loads of the type intended. Advice should be sought from the AE(D).

Automatic control test

- 2.304 The automatic control test is designed to show that the operating cycle functions correctly as shown by the values of the cycle variables indicated and recorded by the instruments fitted to the decontamination equipment.
- 2.305 It should be carried out once a week on most machines and is the main test for ensuring that the equipment continues to function correctly.
- 2.306 During the commissioning, yearly and quarterly test programmes the temperature (and for sterilizers, pressure) sensors for subsequent thermometric tests will be connected to the chamber during this test. If a sensor is placed adjacent to each of the sensors connected to the installed temperature measuring instruments the calibration of these instruments can be checked during periods of stable temperature in the automatic control test.

Method

- 2.307 Place the test load appropriate to the type of washer-disinfector, contained within any load furniture normally used, in the chamber.
- 2.308 For washer-disinfectors equipped with multiple cycle capability select the operating cycle to be tested. Start the cycle.
- 2.309 Ensure that a BPR is made by the recording instrument fitted to the machine. If the machine does not have a recorder observe and note the elapsed time indicated chamber temperatures and pressures at all significant points of the operating cycle, for example the beginning and ending of each stage or sub-stage, and the maximum values during the holding time.
- 2.310 Each stage should be independently timed and the indicated and recorded temperature during these stages logged.

Results

- 2.311 The test should be considered satisfactory if the following recommendations are followed:
- a. a visual display indicates “cycle complete”;
 - b. during the whole of the operational cycle the values of the cycle variables, as indicated by the instruments on the washer-disinfector and any independent monitor or shown on the batch process record, are within the limits established as giving satisfactory results either by the manufacturer or during performance qualification;
 - c. during the disinfection hold period determined from the indicated and/or recorded chamber temperature:
 - (i) the indicated, recorded and any independent monitor chamber temperatures are within the disinfection temperature requirement defined in BS EN ISO 15883 Parts 1 and 2;
 - (ii) the time for which the disinfection temperature is maintained is not less than that previously established, by either the manufacturer or performance qualification tests, as necessary to ensure that the load is maintained at temperatures within the disinfection temperature requirement defined in BS EN ISO 15883 Parts 1 and 2;

- d. the door cannot be opened until the cycle is complete;
- e. the person conducting the test does not observe any mechanical or other anomaly.

- 2.312 Where an independent monitoring system is used which has the necessary data-processing capability, process variability may be monitored automatically through presentation of suitable control charts displaying critical process data. Under these conditions, the need for automatic control tests may be restricted to quarterly, annual and revalidation testing.

Validation and periodic tests for ultrasonic cleaners**Introduction**

- 2.313 Ultrasonic cleaners are of the stand-alone ultrasonic bath type, or included as a separate dedicated chamber in a multi-chamber machine. Many ultrasonic cleaners do not incorporate a disinfection stage and are intended for use as an initial cleaning process before cleaning and disinfection in a washer-disinfector for surgical instruments.
- 2.314 Some ultrasonic cleaners are equipped with means to irrigate hollow instruments such as endoscopes. These washer-disinfectors should be tested both with the general reference load and the endoscope/ MAT reference loads. A testing protocol is shown in [Table 4](#).

Test for ultrasonic activity

- 2.315 The activity of an ultrasonic cleaner can be investigated by the erosion pattern that is created on aluminium foil exposed in the bath for a short period. The activity is not uniform throughout the ultrasonic bath. Tests carried out during commissioning are intended to establish the variation in activity at different positions and levels within the bath and the time required to obtain a characteristic erosion pattern.
- 2.316 The exposure time should depend on the thickness of the foil, the hardness of the foil, the operating frequency, the watt density and the temperature of the ultrasonic bath.

Table 4 Schedule of testing for ultrasonic cleaners

Test	IQ	OQ	PQ	Periodic
Automatic control test	X	X	X	W Q Y
Chamber wall temperature		X	X	Y
Chemical additive(s): low level detection		X	X	
Chemical additive(s): process residue			X	
Chemical: reproducibility		X	X	Y
Cleaning efficacy by residual soil		X	X	W
Cleaning efficacy with test soil		X	X	Q Y
Doors: in-cycle interlock		X	X	Y
Doors: cycle start interlock		X	X	Y
Doors: door-opening force		X	X	Y
Drainage: free drainage		X	X	Y
Fault interlock		X	X	Y
Load carrier temperature test	X	X		
Load carriers				Y
Load dryness test		X	X	Y
Over-temperature cut out test		X	X	Y
Remove and clean strainers or filters				D W
Weekly safety checks		X	X	W Q
Sound pressure test		X	X	Y
Test for ultrasonic activity		X	X	Y
Thermometric test for disinfection		X	X	Y
Verification of calibration	X	X	X	W Q Y
Water: hardness	X			Y
Water: overflow test		X	X	Y
Water supply temperature	X			Y
Volume of water used per stage		X	X	Y

IQ = installation qualification
 OQ = operational qualification
 PQ = performance qualification
 W = weekly
 Q = quarterly
 Y = yearly

Apparatus

2.317 The following apparatus should be used:

- aluminium foil, nominal thickness 0.015–0.025 mm;
- autoclave indicator tape;
- stopwatch, graduated in 0.2 s and with an accuracy over a period of 15 min of ± 0.5 s, or better;
- ruler/tape measure graduated in mm.

Method

2.318 Measure the depth of the bath from the level of the lid to the bottom of the bath. Let the depth be D mm. Cut strips of aluminium foil, 15–20 mm wide and (D + 120) mm. Carry out the manufacturer's recommended start-up procedure.

Note

This will normally include a period of operation to eliminate dissolved gases from the solution in the bath (the de-gassing procedure).

2.319 Ensure that the water in the tank is at the required level, that the required amount of any chemical additive specified by the manufacturer has been added and that the water in the tank is at the specified operating temperature.

2.320 Using strips of autoclave indicator tape across the top of the bath suspend nine strips of the prepared foil in the bath in a 3 × 3 grid.

2.321 The rolled end of each foil strip acts as a sinker weight to maintain the foil in an approximately vertical position. The sinker weight should be not more than 10 mm above, but not touching, the bottom of the bath.

2.322 Operate the bath for the predetermined exposure time. This can vary typically between 30 s for a watt density of 20 W dm⁻³ and 10 min for a watt density of 5 W dm⁻³.

2.323 Remove the strips from the bath, blot dry and examine.

2.324 File the strips by sticking them to an A4 sheet of plain paper using a transparent adhesive tape or by lamination in a clear pocket.

2.325 Drain the bath and clean to remove debris of eroded aluminium foil.

Results

- 2.326 The zones of maximum erosion should be at similar positions on all nine foils and each should be eroded to a similar extent (by visual inspection). For precise evaluation the foils should be weighed before and after exposure to ultrasonication and the loss in weight recorded. The variation in loss of weight should be such that the weight of any one foil is within $\pm 20\%$ of the mean loss of weight.
- 2.327 On re-testing, the extent of erosion and the erosion pattern should have remained consistent with those originally determined during commissioning.

Reference test loads

- 2.328 The test load should contain the following general equipment:
- 3 cuscoe speculae;
 - 3 artery forceps (Crile, Kelly or Spencer Wells) with box joints;
 - 3 No. 3 Scalpel handles;
 - 3 Yankauers or Pooles suction tubes;
 - sufficient additional instruments to make up a full load.
- 2.329 The test load should contain the following endoscope/MAT instruments:
- 2 Trochar and Cannulae;
 - 2 MAT forceps;
 - 2 surrogate endoscopes (see next paragraph);
 - sufficient additional instruments to make up a full load.

- 2.330 The surrogate endoscope should be constructed from 6 mm OD/4 mm ID stainless steel tubing. The overall length should be 450 mm. At the midpoint of the tube should be a 50 mm length of tubing connected to the tubing on either side with compression fittings.
- 2.331 The 50 mm demountable length can be used to provide a more readily visible section for determination of cleaning efficacy.

PQ tests*Load items*

- 2.332 Difficult to clean laboratory items, other than those of the type included in the reference load, or other items (for example hollowware) for which standard reference loads have been defined, should be reviewed to determine how well they are represented by the items of which the reference loads are composed. If the reference loads do not adequately represent the loads to be used further tests should be carried out using loads composed of items which will be in normal production loads.

Nature of soiling

- 2.333 Ultrasonic cleaners are often used for items that are contaminated with soiling which is difficult to remove by other cleaning processes.
- 2.334 The test soil for operational testing should be chosen to represent biological fluids that might be present. If other types of soiling will be encountered (for example orthopaedic bone cement) tests should be conducted using items soiled in the manner that occurs for normal production loads.

3 Water supply

Introduction

- 3.1 All the organisations responsible for water supply have the statutory power to make and enforce by-laws to prevent waste, excessive consumption, misuse or contamination of the water supply. Washer-disinfectors should be designed, constructed, installed, operated and maintained in accordance with the requirements of the relevant by-laws.
- 3.2 The number, nature and quality of water supplies required are dependent on the size and type of washer-disinfectors.
- 3.3 Washer-disinfectors can be supplied with both hot and cold water. When hot water is required as part of the operating cycle, it is generally advantageous to supply hot water to the washer-disinfectors rather than heat cold water to the required temperature within the washer-disinfectors.
- 3.4 The quality of water used at all stages in the decontamination process is critical to the successful outcome of the process.
- 3.5 At each stage the water quality should be compatible with:
- the materials of construction of the washer-disinfectors;
 - the load items to be processed;
 - the chemical additive used;
 - the process requirements of that particular stage.
- 3.6 The key factors are:
- hardness;
 - temperature;
 - ionic contaminants (for example heavy metals, halides, phosphates and silicates);
 - microbial population;
 - bacterial endotoxins.

Water hardness

- 3.7 Hard water is caused by the presence of dissolved salts of the alkaline earth metals (calcium, magnesium and strontium), which come out of solution and deposit as hard mineral layers (limescale) when water is heated or evaporated.
- 3.8 The fouling of electrical heating elements or heat exchange components by hard water dramatically reduces the heat-transfer efficiency and can quickly lead to an increase in heating costs of 50–100%.
- 3.9 The deposition of limescale within pipes and around the edges of spray nozzles will seriously impair the performance of a washer-disinfectors. Hard water will cause scaling on the edges of spray nozzles even when fed with only cold water.
- 3.10 The presence of hardness in water seriously impairs the efficiency of most detergents and disinfectants. If the use of hard water is unavoidable process chemicals containing sequestering agents should be used. This adds considerably to the cost of the process.
- 3.11 Using hard water in the thermal disinfection and final rinse stages of the washer-disinfectors cycle is one of the major causes of white powdery deposits on load items. These are not only unsightly and an unwelcome contaminant but act as a focus for soiling and recontamination of the item in use. In some applications (for example with optical systems) such deposits can seriously impair the utility of the item.
- 3.12 For the cleaning – that is, flushing and washing – stages most washer-disinfectors will operate with water hardness values up to 125 mg/L CaCO₃ but are more effective and cheaper to operate when the hardness of the water does not exceed 50 mg/L CaCO₃.
- 3.13 Some washer-disinfectors are fitted with integral water treatment systems.

3.14 The temperature at which water is supplied to each stage of the process has a major effect on the efficacy of the process:

- Water at too high a temperature during the flushing – i.e. the initial cleaning stage - stage can lead to the coagulation of proteins and thus serve to “fix” proteinaceous soil to the surface of the load items. BS EN ISO 15883-1 recommends that the initial temperature should not exceed 35°C. The flushing stage should be supplied with water from a cold supply.
- Water at too low a temperature during the washing stage – i.e. heated cleaning stage – of the cycle will often impair the ability of detergents used to remove soils composed largely of fats, oils or grease.

3.15 When enzymatic cleaners are used the water temperature should be maintained close to the optimum temperature specified by the manufacturer; too high a temperature will inactivate the enzymes.

3.16 When chemical disinfectants are used, the rate of activity generally increases with increased temperature. Too low a temperature will cause failure to attain the required microbial activation. However, too high a temperature with particular compounds can lead to degradation of the active components, evolution of toxic vapours or adverse reactions with the load items being processed.

3.17 The maximum temperature of rinsing water should be compatible with the items being processed; many items used in medical practice are temperature sensitive or might be damaged by thermal shock.

Ionic contaminants

3.18 To avoid the risk of corrosion, water used in the cleaning of stainless steel instruments should have a chloride concentration less than 120 mg/L and, for final rinse/disinfection, less than 10 mg/L Cl⁻. Chloride concentrations greater than 240 mg/L Cl⁻ cause pitting to occur.

3.19 Tarnishing of stainless steel instruments, shown by blue, brown or iridescent surface coloration, occurs when heavy metal ions, such as iron, are present in the process water. In hot water (over 75°C)

magnesium ions and silicates can cause similar discoloration.

Microbial population

3.20 The purpose of the cleaning and disinfection process is to remove soiling and reduce the microbial contamination to an acceptable level for the intended use of the items to be processed. The water used at each stage of the washer-disinfector process cycle should not increase the bioburden of the load items.

3.21 For terminal disinfection, where items are intended to be used without further decontamination processing, the nature and extent of the microbial population in the final rinse or disinfection water should not present a potential hazard to the patient, either through infection or by leading to a erroneous diagnosis. Appropriate treatment to control or reduce the microbial contamination in water might be required.

Concentration of bacterial endotoxins

3.22 Bacterial endotoxins (see [paragraph 3.194](#)) are thermostable compounds derived from the cell walls of bacteria which, when introduced into the human body, can cause a fever-like reaction and other adverse effects. They are not readily inactivated at the temperatures used for disinfection or sterilization.

3.23 Water used for the final stages of processing in a washer-disinfector, where there is a significant risk of residual water remaining on the load items, should not contain more than 0.25 EU mL⁻¹ when the washer-disinfector is being used to process surgically invasive items or those which are intended to come into contact with parenteral solutions.

Water treatment

3.24 Despite the cost involved in treating water from the public supply, to provide the optimum quality for use at each stage in the washer-disinfector process cycle, treatment is usually cost-effective. The use to which water of various qualities should be put is shown in [Table 5](#).

Table 5 Use to which water of various qualities may be put

Types of water	Application
Cold water	Flushing
Potable water	Flushing/washing
Soft potable water (<50mg/L / CaCO ₃)	Cleaning with detergent or enzymatic cleaners
Softened water-based exchange softener	Desirable in all water >50mg/L / CaCO ₃
	Essential in all water >125mg/L / CaCO ₃ for use in WDs
Deionised water	Final rinse water in all WDs
	Thermal disinfection in all WDs
Reverse osmosis	Final rinse water in all WDs
	Thermal disinfection in all WDs

Note: This table shows suitable applications for the various qualities of water commonly available. Although water of lower quality may be used, this will normally require additional chemical additives and might entail some impairment of the WD performance. The terms “softened”, “deionised” and “reverse osmosis” refer to defined purification processes.

Chemical purity

3.25 There are generally three methods of water treatment available for use on water supplies to be used in washer-disinfectors:

- water softeners;
- water deionisers;
- reverse osmosis (RO).

Water softeners

3.26 Water softeners, or “base-exchange” softeners, consist of an ion-exchange column containing a strong cation resin in the sodium form. Calcium and magnesium ions in the water are replaced by sodium ions. The column may be regenerated by treatment with a solution of common salt (sodium chloride). Base-exchange softeners can cause a significant increase in the microbial content of the water.

3.27 The concentration of total dissolved solids in the water is not reduced by this process. The sodium salts that remain do not readily form hard deposits to foul heat exchangers or spray nozzles but if used as the final rinse or disinfection will leave white deposits on the load items as they dry.

3.28 The process is simple to operate with an automated in-line system, will handle water with varying levels of hardness, and is simple and safe to regenerate. After regeneration, however, high levels of chloride ions might be present in the initial output from the softener, which should be run to waste.

3.29 In common with other water treatment systems, the base-exchange softener should run to a minimum volume of out-flow if the required water quality is to be achieved. This volume should be specified by the manufacturer of the treatment plant. The output from the softener should be to a water tank and the volume demanded each time additional water is fed to the tank should exceed the minimum flow.

Integral water softener

3.30 Some washer-disinfectors are available with built in base-exchange water softeners although these are generally laboratory washer-disinfectors.

3.31 Water softeners should be chosen based on the total demand of softened water in the unit, e.g. SSD, including when necessary provision for manual washing facilities and other plant.

Deionisers

3.32 Deionisation or demineralisation systems can remove virtually all the dissolved ionic material by ion-exchange using a combination of cation and anion exchange resins either in a single column (mixed bed) or in a separate column.

3.33 Operating costs of mixed bed deionisers are usually higher than for two-stage systems.

3.34 Routine maintenance (regeneration) of deionisation and demineralisation systems requires the use of a strong acid (hydrochloric acid) and a strong alkali (sodium hydroxide). For most types of installation, an exchange column service is available from the water treatment suppliers. The maintenance of these systems in line with manufacturers’ requirements is essential to safeguard output quality.

3.35 Deionised water may become contaminated with microorganisms and the resin column colonised. Deionised water should not be used for the final rinse of products intended for invasive use without further decontamination processing by heating, filtration etc. It is essential that your risk assessment in this area and related local policy establish safe

water supply appropriate for each stage of the process.

- 3.36 Systems are available in a range of sizes from small wall-mounted units in which ion-exchange resins are contained in disposable cartridges to large industrial units. Regeneration requires the use of strong acid (hydrochloric acid) and strong alkali (sodium hydroxide). For most types of installation an exchange column service is available from the water treatment suppliers.
- 3.37 For a given output volume, the initial cost of providing deionisation equipment will be lower than for reverse osmosis (RO). However, the inconvenience and cost of the regeneration process for deionisers, and the better microbial quality of the RO process, makes RO the preferred option.

Reverse osmosis (RO)

- 3.38 RO treatment plants remove almost all dissolved inorganic contaminants by passing the water, under pressure, through a semi-permeable membrane against an osmotic gradient. The process will also remove a high proportion of organic material, bacterial endotoxins and microorganisms. Some RO units are fitted with a final 0.2 µm filter to control bacterial numbers.
- 3.39 The initial capital cost of an RO plant is generally higher than for a deionisation system supplying a similar volume of water, but operational costs are lower. The water has a low microbial population. Measures are required to maintain the microbial quality of water during storage and distribution. The retention of this water quality requires a high level of understanding and maintenance.
- 3.40 The wastewater produced by properly RO plant may be designated as grey and reused appropriately.
- 3.41 Issues to be considered when installing RO systems include:
- An RO system removes bacteria, endotoxins and approximately 95% of chemical contaminants.
 - The system and associated pipework need to be sanitised regularly.
 - The system provides processed water over a long period with minimal maintenance.
 - Routine maintenance and membrane replacement are very important.
- The system usually requires a carbon filter to be fitted ahead of the RO unit to remove traces of chlorine from the water supply.
 - RO units use large volumes of water, much going to waste, which can be used as grey water.
 - Selecting an RO system will depend on the geographical area of the water source and its quality .
 - If the supply water is hard, a softening system will be required ahead of the RO unit.
- 3.42 Water storage is required, as RO units supply moderate volumes of water over a long period. Washer-disinfectors need large volumes of water quickly during various stages of the cycle.

Microbial purity

- 3.43 Potable water from the public supply has a low microbial content and should be free from pathogenic organisms, other than those that might cause opportunistic infections in immunologically compromised patients.
- 3.44 If stored in tanks or cisterns, the microbial content can increase considerably.
- 3.45 Attention is drawn to the requirement under the code of practice for control of legionella that water in intercepting tanks should be stored below 20°C or above 55°C.
- 3.46 The extent and nature of microbial contamination in the water supplied to a washer-disinfector will depend on the stage in the process cycle at which it is to be used and the intended use of the decontaminated load at the end of the process.
- 3.47 Water stored at 60°C or above may be assumed not to have a proliferating microbial population.
- 3.48 When water is treated by filtration, for example through a 0.22 mm filter to remove microbial contaminants, rigorous controls are needed to ensure that the system works effectively. These should include:
- either maintaining the pressure drop across the filter throughout its working life, a decrease in differential pressure being cause for rejection of the process cycle and a change of filter, or, a bubble point test before and after each process cycle (see BS 1752, ISO 4793);
 - a continuous recirculation system so that the filter is not left wet in static water;

- treatment of the circulating water to ensure that proliferation of microbial contamination is inhibited either by use of elevated temperature (for example $>60^{\circ}\text{C}$) or by the use of UV irradiation (wavelength $260\text{ nm} \pm 10\text{ nm}$; $>2\text{ J m}^{-2}$).

3.49 Verification of purification by filtration should be made by relevant TVC (total viable count) tests.

Pipework

- 3.50 The pipework used to supply the various grades of water should be appropriate to the quality of water carried
- 3.51 All pipework should be run with a continuous fall to the discharge point so that it is free draining. It should be free from dead ends and other areas where water can become stagnant.

Water supply by-laws

- 3.52 All the organisations responsible for water supply within the UK have the statutory power to make, and the duty to enforce, by-laws for the prevention of waste, undue consumption, misuse or contamination of the water supplied by them (see the Water Supply (Water Fittings) Regulations 1999).
- 3.53 By-laws 38 to 41 require storage cisterns to be fitted with warning pipes (and an overflow if in excess of 1000 L capacity).
- 3.54 The warning pipe and overflow should not comprise, or have connected to it, a flexible hose.
- 3.55 By-law 25 Schedule A gives examples of points of use or delivery of water where backflow is, or is likely to be, harmful to health due to a substance continuously or frequently present (By-law 25 (1) (a)).
- 3.56 This schedule also lists water softening treatment plant, bedpan washers, bottle washers, dishwashers and disinfection equipment and clearly applies to all washer-disinfectors.
- 3.57 The required protection is a Type A air gap at the point of use or an interposed cistern.
- 3.58 Water softeners, regenerated only by means of sodium chloride solutions, need only be protected by a Type B air gap.

Water system tests

- 3.59 A continuous supply of water of the specified chemical and microbial quality is essential to the correct functioning of all washer-disinfectors.
- 3.60 Water that is too hard or has too high a concentration of dissolved solids can impair the activity of detergents (or require the use of increased quantities of chemical additives) and cause deposits, scaling or corrosion of items being processed.
- 3.61 Water containing high numbers of microorganisms may recontaminate disinfected items. For all these tests the water should be sampled from the water supply pipe to each washer-disinfectors. Samples may need to be taken from additional points in the supply when trying to identify the cause of a non-conformity.

Water samples

- 3.62 Water samples should be obtained from draw-off points installed at convenient locations within the system, as close to the washer-disinfectors as possible.
- 3.63 The sampling procedure should be suitable for all the physical, chemical, and biological determinands of interest. It may be used for water samples throughout the water distribution system.
- 3.64 The sampling containers used should be specific for the determinands of interest. This should include, as appropriate:
- 250 mL sterile pyrogen-free single-use containers (for determination of bacterial endotoxin levels and/or total viable count);
 - 1 L acid-washed borosilicate bottles (for determination of cations);
 - 1 L polypropylene bottles (for determination of anions, total dissolved solids);
 - 100 mL high-density polyethylene bottles (for determination of pH, conductivity).
- 3.65 The first 50 mL of sample taken at each sampling point should be run to waste.
- 3.66 All samples should be taken in duplicate.
- 3.67 Samples should be stored at $2\text{--}5^{\circ}\text{C}$ and tested within 4 h of collection.

Water quality tests

3.68 The following clauses recommend analytical methods to determine the various biological, physical and chemical properties of water samples for the various qualities of feedwater to the washer-disinfector. A list of acceptable results for the analysis is provided in Table 6.

Table 6 Requirements for water quality: final rinse and process water

Determinant and unit	Maximum permitted values	
	Final rinse	Other stages
Appearance	Clear, colourless	–
Degree of acidity (pH)	5.5 to 8.0	–
Conductivity at 25°C (uS/cm)	30	–
Total dissolved solids (mg/100 mL)	4	–
Total hardness, CaCO ₃ (mg/L)	50	210
Chloride, Cl (mg/L)	10	120
Heavy metals, determined as Lead, Pb (mg/L)	10	–
Iron, Fe (mg/L)	2	–
Phosphate, P ₂ O ₅ (mg/L)	0.2	–
Silicate, SiO ₂ (mg/L)	0.2	2
Total viable count (TVC) at 22°C	100	–
at 37°C (cfu/100 mL)	100	–
Bacterial endotoxins (EU/mL)	0.25	–

3.69 The methods of analysis recommended to detect chemical contaminants at low concentrations with a high level of accuracy require the use of a laboratory with appropriate expertise, facilities and experience.

3.70 Other tests can be carried out on-site or with very simple laboratory facilities; these lack the precision and sensitivity of the laboratory tests but are sufficient for most purposes.

3.71 The following clauses contain detailed procedures for tests that may be carried out on-site or with very simple laboratory equipment at, or shortly after, the time of sampling.

3.72 The precision, accuracy, sensitivity and limits of detection of these methods are usually inferior to those of laboratory methods. They are useful,

however, in that they provide evidence of any gross failure and the results are available straightaway making them of diagnostic value during a fault-finding exercise.

3.73 They are generally economical compared with more sophisticated laboratory analysis and can be carried out by non-specialist personnel after thorough, but limited, training. The results should not, however, be used as evidence in cases of dispute.

Choice of method

3.74 For any given determinant there will usually be several methods that are suitable and cover the range of concentrations of interest. The methods recommended here are intended to be representative of those that might be suitable. They are chosen as examples of tests which can conveniently be carried out on site.

3.75 A number of test systems are available commercially. Before adopting one of these methods care should be taken to ensure that the test(s) provides results of sufficient accuracy and sensitivity.

3.76 It is not necessary to use experienced chemical analysts to undertake the on-site analysis of water samples described. It is, however, essential that personnel receive appropriate training before attempting to carry out this work. Recourse to more precise analysis might be needed in the event of a dispute between two parties.

3.77 Many contaminants can be detected by two or more of the determinations normally carried out for laboratory analysis. For example an increase in one or other of the ionic species present will cause an increase in electrical conductivity and an increase in the evaporative residue as well as showing an increase in the concentration of that particular ion.

3.78 Further guidance on appropriate test methods may be obtained from BS 1427.

3.79 Tests suitable for use on-site fall into three main categories:

- instrumental tests using portable instruments designed for on-site use, for example portable pH meters, ion selective electrodes etc;
- spectrophotometric tests based on measurement of the absorbance of a coloured reaction product; measurement can be visual or photometric and can be against a precalibrated

coloured disc or against standard reference solutions;

- c. titrimetric tests may be carried out using standard laboratory equipment or with commercially available apparatus designed for field use; the latter is usually much simpler to use.
- 3.80** For all the instrumental methods recommended there is commercially available equipment specifically intended for field use. All the variables for which instrumental methods are recommended are temperature dependent. The equipment used should be temperature compensated. Also the equipment should be allowed sufficient time on site, before it is put into use, to equilibrate to the local ambient temperature.
- 3.81** Commercially available test kits based on visual or photometric comparison with coloured discs have become an accepted standard for on site analysis. Manufacturers usually supply a complete test system, including kits of reagents. To ensure compatibility, and maintenance of the manufacturer's claimed sensitivity and accuracy for the method, the kit specified by the manufacturer should not be substituted.

Water supply temperature

- 3.82** The water supplied to the various stages of the washer-disinfector operating cycle should be at an appropriate temperature. If the temperature of the water supplied to the flushing stage is too high (>35°C) there is a risk of coagulating proteinaceous soiling, which inhibits the cleaning process. If the temperature of water supplied to the washing, rinsing and disinfection stages is too low, the washer-disinfector cycle can be greatly extended, with a significant reduction in throughput, while the water is heated to the required temperature within the washer-disinfector. Water supplied in the temperature range 25–40°C presents a serious risk of microbial contamination of the system.

Apparatus

- 3.83** An indicating or recording thermometer should be used.

Method

- 3.84** Measure the temperature of the water supply from a sampling point as close to the washer-disinfector as possible. Place the temperature sensor in the

middle of the flowing stream as close as practicable to the sampling point. Allow the water to flow for at least a minute before the temperature is read.

Alternative method (for periodic testing)

- 3.85** When it is not convenient, or practicable, to run the water to waste from a sampling point close to the washer-disinfector the water temperature can be estimated by measurement of the temperature of the outer surface of the supply pipe. For this method the correlation between the temperature of the water flowing out of the pipe and the surface temperature of the pipe at a particular point should be established during installation testing.
- 3.86** Measure the surface temperature using a sensor designed for the purpose and follow the manufacturer's instructions for ensuring good thermal contact with the surface. Record or note the temperature during a normal operating cycle not less than 30 s after the start of water flow through the pipe to the washer-disinfector.

Results

- 3.87** The noted value should be within the temperature range specified for the installation.

Water supply pressure

- 3.88** If the pressure of the water supply to the washer-disinfector is below the minimum pressure specified by the manufacturer, the performance and productivity of the washer-disinfector will be affected adversely.
- 3.89** If the pressure of the water supply to the washer-disinfector is above the maximum pressure specified by the manufacturer the capacity of overflow devices may be inadequate, the designed performance characteristics of valves etc may be exceeded and in extreme cases there may be the risk of damage to components of the washer-disinfector or to products being processed.

Note

For example, many flexible endoscopes are likely to be damaged if subjected to internal pressures greater than 35 kPa.

- 3.90** The test should be carried out as an installation and/or operational test. The test should be repeated when any change is made to the water services

supplying the washer-disinfector (including the connection or removal of additional machines).

Apparatus

- 3.91 A pressure indicator or recorder 0–10 barg should be used.

Method

- 3.92 Connect the pressure sensor to each of the water supply pipes to the washer-disinfector, as close to the washer-disinfector as practicable, on the supply side of the washer-disinfector isolating valve for that supply. Record or observe and note the static pressure when the valve is closed and the pressure indicated throughout a normal operating cycle. When the water service also supplies other equipment on the same supply line, run the test both with the other equipment operating throughout the test (or their operation simulated by an appropriate discharge to waste) and with no other equipment operating.

Results

- 3.93 The water pressure should remain within the supply pressure limits specified by the washer-disinfector manufacturer.

Appearance

- 3.94 All the water supplied to the washer-disinfector should be clean, colourless and free from visible particulate matter. The appearance of the sample should be assessed visually.

Apparatus

- 3.95 The following equipment should be used:
- a. clean, clear glass bottle and stopper;
 - b. filter paper (qualitative grade), filter funnel and holder.

Method

- 3.96 Transfer an aliquot to a clear colourless glass bottle, which should then be tightly stoppered. Shake the sample well and examine visually against a white background, preferably in a north light.
- 3.97 If the sample is turbid filter through a qualitative grade filter paper. Examine the filter paper and report a description of the retained material. Visually examine the filtrate as previously described.

- 3.98 Report the appearance in terms of both colour and the intensity of any colour. If the sample is coloured, examine it carefully to see whether visible evidence of colloidal material is present.

Results

- 3.99 All the samples tested should be clear, bright and colourless.

pH

- 3.100 Two suitable methods for on-site measurement of pH are available; the colour disc comparator and portable pH meter.

Apparatus

- 3.101 The following equipment should be used:
- a. pH meter:
 - (i) Several commercially available small portable pH meters are available. A number of these include built-in temperature compensation. These should not be used to measure pH solutions of low ionic strength, although they provide suitable accuracy for most general applications. Their use for the determination of pH in water of high purity might give unstable or unreliable readings.
 - (ii) Only those pH meters specifically designed for the measurement of low ionic strength solutions should be used for determining the pH of DI or RO water.
 - b. Colour disc comparator:
 - (i) Colorimetric tests for pH are suitable for high purity, low conductivity, water samples of the type required to be tested.
 - (ii) Since colorimeter methods are being used for other field tests this may be the more appropriate method.
 - (iii) The accuracy is limited and discrimination may not be better than 0.2 pH units. This is, however, quite suitable for field tests.
 - (iv) Colour, turbidity or strong oxidants in the sample all interfere with the test.
 - (v) A narrow range indicator (or two for use on successive samples) should be chosen

to cover the required range of pH 4 to pH 1. Manufacturers of colorimeters usually provide indicators to cover a range of 2 or 3 pH units. Wide-range indicators should not be used because of their poor discrimination.

Method

- 3.102 Operate the test kit in accordance with the manufacturer's instructions. Pay particular attention to using accurate volumes of both sample and reagent and monitoring both temperature and reaction time.
- 3.103 Match the colour of the reacted sample against the calibrated colour disc viewed through a blank sample. Read off the value in pH units directly from the disc.
- 3.104 Verify the calibration using standard buffer solutions made up in advance and kept in capped bottles until required. The buffer solutions should be chosen to have a pH in the midpoint of range of the calibrated colour discs to be used in the determination.

Results

- 3.105 The indicated value should be in the range 5.5 to 8.0. Photometric apparatus with somewhat better discrimination is also commercially available.

Electrical conductivity

Equipment and materials

- 3.106 There is a wide variety of portable conductivity meters available. The unit chosen should meet the performance criteria given in this clause, which will allow its use for measuring conductivity of very pure water through to boiler water containing in excess of 5000 ppm TDS.
- 3.107 Meters are available to cover the range 0.1–20,000 mS m⁻¹ at 25°C. Conductivity meters may also be calibrated in µS cm⁻¹. The meter, or meters, used should cover the following ranges:

Range	Resolution	Accuracy
0–199 µS cm ⁻¹	0.1 µS cm ⁻¹	±1% full scale
10–1990 µS cm ⁻¹	1 µS cm ⁻¹	±1.5% full scale

- 3.108 They should also be temperature-compensated over the range 0–40°C.

Method

- 3.109 Use the following method for electrical conductivity calibration.
- 3.110 Verify the calibration of the meter against 0.001 molar and 0.0005 molar solutions of potassium chloride and pure water as working standards. These give conductivities at 25°C of 14.7 mS m⁻¹ and 84 µS cm⁻¹ and <0.06 µS cm⁻¹.
- 3.111 Prepare the potassium chloride solutions by dilution of a 0.1 molar solution.
- 3.112 The working standards are stable for up to one week when stored in cool conditions, in a well-stoppered container.
- 3.113 A comprehensive range of standard conductivity reference solutions, including pure water standards (also known as absolute water) are available commercially, standardized at 25°C and traceable to national standard reference materials.
- 3.114 After calibration rinse the sample cup or immersion probe thoroughly with pure water.
- 3.115 Collect the sample in a high-density polyethylene bottle and test as soon as practicable.
- 3.116 Pour an aliquot of the sample into the sample cup of the conductivity meter or, for meters with an immersion probe, into the clean beaker. Follow the meter manufacturer's instructions for making the measurement; this will usually require a short stabilization period before noting the reading.

Results

- 3.117 The conductivity at 25°C should not exceed:
- 30 µS cm⁻¹ for reverse osmosis water;
 - 30 µS cm⁻¹ for deionised water;
 - 300 µS cm⁻¹ for softened or mains water.

Note

Conductivity levels in excess of this value are indicative of a high concentration of dissolved solids.

Total dissolved solids

- 3.118 The laboratory test for the determination of dissolved solids is a gravimetric method. This involves determining the weight of the residue obtained by evaporating a known sample volume to dryness.

- 3.119 When a water sample contains predominantly ionisable solids in solution, and the composition of the various constituents is reasonably constant, a good estimate of the total dissolved solids can be obtained from the electrical conductivity of the sample which can be used to determine concentrations up to 10,000 mg/L total dissolved solids.

Conductivity method

- 3.120 The following apparatus should be used:

- conductivity meter (see paragraph 3.106, 'Electrical conductivity');
- phenolphthalein indicator;
- 5% w/w acetic acid solution;
- 5% w/w sodium hydroxide solution.

- 3.121 Neutralize the test sample, using phenolphthalein as the indicator, by drop-wise addition of 5% w/w sodium hydroxide solution or 5% w/w acetic acid solution. Measure the conductivity of the sample (see 'Electrical conductivity') and multiply by the conversion factor to give an estimate of the TDS in mg/L.

- 3.122 The conversion factor can be derived experimentally for waters of consistent ionic composition by making direct comparison of the measured mass of total dissolved solids and the electrical conductivity.

- 3.123 Alternatively, an arbitrary factor can be used. The one most commonly chosen is based on sodium sulphate being the ionic species.

- 3.124 For conductivity at 25°C measured in mS m^{-1} :
- $$\text{TDS mg/L} = \text{electrical conductivity mS m}^{-1} \times 1.6.$$

- 3.125 Conductivity meters calibrated in TDS mg/L are also available. Care should be taken to ensure that the conversion factor used is appropriate.

- 3.126 Standard salt solutions are available commercially as ready-to-use standard solutions traceable to NIST standard reference materials. A TDS standard solution such as 1382 ppm NaCl, and a tenfold dilution of it, can be used to verify the calibration.

Results

- 3.127 The estimate of total dissolved solids should not exceed 4 mg 100 mL⁻¹ for purified water (RO or DI).

Evaporative residue method

- 3.128 The following apparatus should be used:

- silica or borosilicate dish or beaker of >150 mL capacity;
- oven set to 110°C ± 2°C;
- boiling water bath or heating mantle set to 100°C ± 2°C;
- 1 L polypropylene bottle;
- balance weighing to 0.1 mg;
- 100 mL pipette or measuring cylinder.

- 3.129 Collect a 1 L sample.

- 3.130 Take the silica dish (or equivalent), dried for 2 h in the oven set to 100°C ± 2°C and then cooled to ambient temperature, and weigh it to the nearest 0.1 mg.

- 3.131 Dispense 100 mL of the sample into the weighed dish and evaporate it over the boiling water bath until visibly dry. Evaporate two further 100 mL aliquots of the sample in the same dish in the same manner.

- 3.132 Dry the dish in the oven to constant weight to an accuracy of 0.1 mg.

- 3.133 Calculate the mass of residue in the dish and hence calculate the mass of residue per 100 mL of water.

Results

- 3.134 The evaporative residue should not exceed 4 mg 100 mL⁻¹ for purified water (RO or DI).

Hardness (as CaCO₃)

- 3.135 Hardness of water is due to the presence of dissolved salts of the alkaline earth metals (calcium, magnesium and strontium). Their presence causes limescale formation from heated or evaporated water, can inactivate detergents and disinfectants and causes scaling on load items.

Ion-selective electrodes (ISE) method

- 3.136 Ion-selective electrodes are available for calcium and also for divalent cations (total hardness). Ion-selective electrodes are not specific for a particular ion but have a relative selectivity for a particular ion or group of ions. They are sensors that provide a potentiometric response to the activity of the ions in solution. The activity is proportional to the

concentration for determinations carried out in solutions of the same ionic strength.

- 3.137 Adjust both analyte and calibration standard solutions to the same ionic strength. Use a high impedance millivoltmeter to measure the potential between the ion selective electrode and a suitable reference electrode. The measured potential is proportional to the logarithm of the concentration of the ion(s) in solution.
- 3.138 The optimum working pH range is 4 to 9 and the ionic strength of the sample should be adjusted for ionic strength. An adjustment buffer consisting of 4M KCl solution is often used. Phosphate buffers should not be used since the calcium activity will be lowered by complexation or precipitation.
- 3.139 The electrodes are free from any major interference except zinc ions. They are however poisoned by a number of biological fluids.
- 3.140 The calcium electrode requires a single junction reference electrode. Calibration is made against two or more standard solutions. These are commercially available.
- 3.141 The calcium selective electrodes that are available have a Nernstian response for concentrations from 1 M down to about 5×10^{-6} M and a selectivity ratio of better than 2000 against magnesium. This range is suitable for analysis of softened water and purified water (RO or DI).

Titrimetric method

- 3.142 Commercially available kits for the titrimetric determination of both total hardness and calcium hardness are available. They are based on the same reaction in which divalent cations are complexed with the disodium salt of ethylenediaminetetraacetic acid (EDTA). When the reaction is carried out, at pH 10 to 11, with eriochrome black as the complexometric indicator, all the calcium and magnesium ions are chelated by the EDTA and the absence of free calcium and magnesium ions causes a colour change in the indicator.
- 3.143 At pH values above 12 magnesium ions are precipitated as the hydroxide and do not react with the EDTA. Calcium hardness can be determined using Patton and Reeders indicator powder as a complexometric indicator.
- 3.144 The commercially available kits often use novel titration methods instead of burettes. The test

reagents are specific to each kit. The manufacturer's instructions should be followed.

- 3.145 Range: determinations within the range 5–400 mg/L can be made. The method is not applicable to purified water or condensate from clean or pure steam, which should have calcium concentrations well below the range for accurate determinations.

Results

- 3.146 The hardness expressed as mg/L CaCO_3 should not exceed 50 mg/L for softened water.
- 3.147 Water with values >210 mg/L should be regarded as unsuitable for use in washer-disinfectors without treatment.

Chloride

- 3.148 The presence of significant levels of chloride ions in water supplied to washer-disinfectors may cause pitting and corrosion in metallic items in the load (including stainless steel). Significant levels of chloride are present in untreated mains water supplies to which it is added for its anti-microbial activity. High chloride concentrations can be associated with breakthrough from a defective, or incorrectly operated, water softener or deioniser.

Ion selective electrode (ISE) method

- 3.149 The commercially available chloride selective electrodes have a working range from 1 M to 10^{-5} M. They work over the pH range 3–10; adjust the sample for ionic strength using an adjustment buffer consisting of 5 M NaNO_3 solution.
- 3.150 The electrodes show poor selectivity against other halides and cyanide ions. Sulphide ions should be absent.
- 3.151 The chloride electrode requires a double junction 0.1 M NaNO_3 reference electrode.
- 3.152 Conduct the calibration against two or more standard solutions.

Note

These are commercially available.

Silver nitrate titration method

- 3.153 Commercial titrimetric kits are available that are based on the method described in BS 6068-2.37:1990, ISO 9297:1989.
- 3.154 Titrate the sample at pH 5 to pH 9 with silver nitrate using a potassium chromate indicator solution.
- 3.155 The analytical range is 5–150 mg/L.
- 3.156 This method is not quantitative for purified water, which should have chloride concentrations well below the range for accurate determinations; it can be used, however, as a limit test. The BP limit test, based on comparison of the turbidity obtained from a known chloride concentration, can also be used.

Results

- 3.157 The chloride concentration in final rinse water for washer-disinfectors processing metal items should not exceed 10 mg/L.
- 3.158 The chloride concentration in other water supplies for washer-disinfectors processing metal items should not exceed 120 mg/L.

Heavy metals

- 3.159 Heavy metals are generally toxic in low concentrations and, as far as possible, should be absent from water used to process items that will be used invasively.

Method

- 3.160 Determine the total concentration of heavy metals using the limit test described in the British Pharmacopoeia (see also ISO 8288).

Results

- 3.161 The total concentration of heavy metals should not exceed 10 mg/L determined as lead.

Iron

- 3.162 The presence of significant concentrations of iron in water used to process stainless steel items promotes corrosion of those items and exacerbates the effect of any chloride ions that might be present.
- 3.163 One of the commercially available colour disc comparator kits should be used for this test. Typically these are based on the reference method

described in BS 6068-2.2:1983, ISO 6332-1982. The reaction of iron (II) with 1,10 phenanthroline in solution yields a red complex with peak absorption at around 510 nm. Most kits include methods and reagents for pretreatment to reduce any iron (III) compounds to the iron (II) form in which they can be analysed.

- 3.164 This method is generally suitable for determination of the concentration of iron in untreated water but is not sufficiently accurate for determination of the concentration specified for steam condensate which at ≤ 0.1 mg/L is at the limit of discrimination of most systems.

Apparatus

- 3.165 The following apparatus should be used:
- a. colour disc comparator kit;

Note

The analytical range depends on the calibrated colour disc supplied with the chosen test kit. A range of 0–5 mg/L is commercially available and provides adequate precision. Discs offering extended ranges are also available should not be used as the discrimination of intermediate concentrations becomes unacceptably poor.

- b. reagents;

Note

The prepackaged reagents available from the comparator manufacturer should be used.

- c. a standard 0.702 g/L iron (II) ammonium sulphate solution $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2$;
- d. a mercury in glass thermometer graduated in 0.5°C steps conforming to BS 1704:1985, ISO 1770:1981.

Method

- 3.166 Prepare a standard 0.702 g/L iron (II) ammonium sulphate solution $(\text{NH}_4)_2\text{Fe}(\text{SO}_4)_2$ which provides a standard solution of 100 mg/L iron. Prepare the solution as required and do not store. Prepare working standards spanning the usable range of the colour disc comparator by appropriate dilution.
- 3.167 Measure the sample temperature before commencing the analysis.

3.168 After the kit manufacturer's specified reaction time has elapsed use the colour intensity of the sample to estimate the concentration of iron in the sample.

3.169 For details of the colorimetric method, see the method for the determination of silicate.

Results

3.170 Untreated and softened water should have less than 2 mg/L iron present.

Phosphate

3.171 The test method measures only orthophosphate. Pre-treatment to convert other forms of phosphate to orthophosphate should be used if appropriate. Some other phosphates such as condensed phosphates and labile organic phosphates are slowly hydrolysed under the acidic conditions used for the test. Undue delay in reading the test can result in a gradual increase in phosphate concentration as hydrolysis proceeds.

3.172 The method depends on the reaction of phosphate in acidic solution with molybdate and antimony ions to form an antimony phosphomolybdate complex, which on reduction with ascorbic acid forms a blue coloured complex having maximum absorbance at 882 nm.

Note

Some commercially available comparators work at 700 nm but this is less sensitive and should not be used.

3.173 The presence of oxidising agents and sulphides will interfere with the reaction. Otherwise there are no particularly sensitive interferences.

Apparatus

3.174 The following apparatus should be used:

- a. colour disc comparator kit;

Note

The calibrated phosphate colour disc should be calibrated in P_2O_5 mg/L. A sensitivity range of 0–5 mg/L is commercially available and provides adequate precision. Discs offering extended ranges should not be used as the discrimination of intermediate concentrations becomes unacceptably poor. Commercially available kits are generally based on the reference method described in BS EN ISO 6878:2004, BS 6068-2.28:2004.

- b. reagents;

Note

The prepackaged reagents available from the manufacturer of the comparator should be used.

- c. sample container;

Note

Phosphate is readily absorbed on to many plastic surfaces. When polypropylene bottles are used as sample containers the sample for phosphate analysis should be transferred immediately to a glass container and assayed as soon as possible.

- d. glassware;

Note

This should be borosilicate which has been subject to acid hardening, that is, cleaned and allowed to stand overnight, filled with sulphuric acid (d201.84), then rinsed several times and stored filled with water, in the dark at 0–4°C until required for use. The glassware should not be allowed to come into contact with detergents or alkaline liquids.

- e. a standardised solution containing 100 mg/L potassium dihydrogen orthophosphate for preparation of calibration standards.

Method

3.175 React the sample in acidic solution with antimony and molybdate ions to form an antimony phosphomolybdate complex. Reduce with ascorbic acid to form a molybdenum blue complex.

3.176 Prepare a stock standard solution containing 100 mg/L potassium dihydrogen orthophosphate and dilute to provide suitable working standards for calibration verification. The concentrated stock solution is stable for several weeks.

3.177 Test the samples as soon as possible after sampling. If sampling will be delayed by more than 4 h store the sample(s) in suitable glass bottles at 2–5°C for up to 24 h.

3.178 Follow the manufacturer's instructions precisely.

3.179 For details of the colorimetric method, see the description given in the method for the determination of silicate.

3.180 The temperature has a significant effect on reaction time, at 20°C the reaction is typically

completed within 3–4 min. Before making the measurement ensure that the reaction is complete, but avoid excessive delays, which can cause errors from hydrolysis of other phosphates. Read the measurement at 10–15 min after the start of the reaction.

Results

3.181 The phosphate concentration of rinse water used for metal load items should not exceed 0.2 mg/L expressed as P_2O_5 .

Silicate

3.182 Silicate reacts with metal items, including stainless steel, causing corrosion and discoloration. This is accentuated at elevated temperatures.

3.183 This method is based on the use of one of the commercially available colour disc comparator kits. Typically these are based on the analytical method described in BS 2690-104, which is a recognised reference method. Reactive silica is reacted with ammonium molybdate under acidic conditions to form molybdosilicic acid which is then reduced to molybdenum blue.

3.184 The analytical range depends on the calibrated colour disc supplied with the chosen test kit. A range of 0–5 mg/L is commercially available and provides adequate precision. Discs offering extended ranges should not be used as the discrimination of intermediate concentrations becomes unacceptably poor.

3.185 The method is generally suitable for determination of SiO_2 level in softened and untreated water but is only sufficiently sensitive to act as a limit test for purified (RO or DI) water.

Apparatus

3.186 The following apparatus should be used:

- a. colour disc comparator kit;
- b. reagents;

Note

The prepackaged reagents available from the manufacturer of the comparator should be used.

- c. a standard 3.132 g/L disodium hexafluorosilicate solution (Na_2SiF_6);

- d. a mercury-in-glass thermometer graduated in 0.5°C steps conforming to BS 1704:1985, ISO 1770:1981.

Method

3.187 Prepare a standard 3.132 g/L disodium hexafluorosilicate solution (Na_2SiF_6), providing a stock standard solution of 1000 mg/L as SiO_2 . The solution is stable for several months after preparation stored in a sealed polyethylene bottle. Working standards spanning the usable range of the colour disc comparator can be prepared by appropriate dilution.

3.188 Measure the sample temperature before commencing the analysis using the mercury-in-glass thermometer.

Note

For most kits the temperature should be 15°C to ensure that the reaction will go to completion. If the sample temperature is below this, or the minimum temperature specified by the manufacturer, the sample should be warmed.

3.189 After the kit manufacturer's specified reaction time has elapsed use the colour intensity of the sample to estimate the concentration of silicate in the sample.

3.190 With the calibrated colour disc for silica in the comparator, an untreated water sample in the blank cuvette, and the reacted sample in the sample cuvette, placed in the comparator cell holder, visually match the colour density developed in the sample against the calibrated colour disc viewed through the untreated sample. Read off the displayed value of SiO_2 concentration from the calibrated disc.

3.191 Serial dilutions of the standard solution may be used to verify the calibration of the comparator disc.

Results

3.192 Untreated and softened water should have less than 2 mg/L silicate expressed as SiO_2 , determined as reactive silica, present.

3.193 Purified (DI or RO) water should have not more than 0.2 mg/L silicate expressed as SiO_2 , determined as reactive silica, present.

Bacterial endotoxins

- 3.194 When the intended use of the washer-disinfectors is for products that will be used invasively, for example surgical instruments, the water used for final rinsing should be tested for bacterial endotoxins (limulus amoebocyte lysate (LAL) test).
- 3.195 The method describes the detection of bacterial endotoxin by the limulus amoebocyte lysate (LAL) gel formation method.
- 3.196 Other LAL methods (chromogenic, turbidimetric or kinetic turbidimetric) are equally suitable.
- 3.197 The water sample is incubated, in a test tube, with an equal volume of lysate for 1 h at 37°C and examined for the formation of a solid clot that holds upon inversion of the test tube. The lysate, reconstituted from lyophilised LAL, is selected with the required level of sensitivity. Semi-quantitative results may be obtained by testing dilutions of the sample to be tested and by the use of lysates with different levels of sensitivity.

Apparatus

- 3.198 The following apparatus should be used:
- pyrogen-free water;
 - LAL reagent (reconstituted limulus amoebocyte lysate);
 - standard endotoxin;
 - reaction tubes;
 - dilution tubes;
 - hot air oven;
 - pipettes;
 - paraffin wax sealing film;
 - non-circulating water bath or dry block incubator capable of maintaining 37°C ± 1°C;
 - test tube racks for reaction tubes;
 - vortex mixer.

Method

- 3.199 Collect and preserve the sample: all apparatus used for collecting samples should be pyrogen-free; metal components may be depyrogenated by dry heat.
- 3.200 Samples should be collected in single-use sterile, apyrogenic polystyrene containers.

- 3.201 Test samples within 4 h of collection. When this is not practicable store samples at 2°C to 5°C for not more than 24 h before testing.
- 3.202 Verify each batch of containers as non-pyrogenic. This can be done by taking representative random samples, rinsing with pyrogen free water and testing the rinse water.
- 3.203 Pre-treat the sample: ensure that the pH of the reaction mixture (sample + lysate) is between 6 and 8.

Note

The presence of substances that denature protein, chelate cations, adsorb, absorb or bind endotoxin or alter the hydrophobic nature of the endotoxin can cause interference, which can be detected by the indication of significantly more or less endotoxin than expected when the sample is tested after the addition of a known amount of standard endotoxin.

- 3.204 Reconstitute the freeze-dried reagent: select the LAL reagent for the required sensitivity to endotoxin; 0.125 EU mL⁻¹ is used to test for compliance with a maximum permitted endotoxin level of 0.25 EU mL⁻¹.
- 3.205 After gently tapping the vial to ensure that all the white powder is at the bottom of the vial, remove the crimp cap and rubber stopper from the vial and discard them.
- 3.206 Add the required quantity of pyrogen-free water (usually 5 mL), allow several minutes for the lyophilized reagent to go into solution and then mix by gentle agitation and store on ice until used.
- 3.207 Re-seal the tube using paraffin wax film. Only the side in contact with the paper backing should be regarded as pyrogen-free.
- 3.208 If not all the reagent is to be used immediately, subdivide into suitable aliquots which may be frozen and stored at -20°C, or below, for up to three months.

Note

Frozen reagent should only be thawed once.

- 3.209 Before use, mix the contents of the vial gently to ensure homogeneity.

Note

Vigorous mixing may cause foaming and loss of sensitivity.

3.210 The following description of the test method is intended for guidance only. The directions given by the manufacturer should be strictly followed.

3.211 Reconstitute the control standard endotoxin: use a similar procedure to that described for the LAL reagent. Reconstitute a vial of freeze-dried control standard endotoxin (CSE). Mix by repeated vortexing, and prepare working standards by serial dilution using pyrogen-free water and pyrogen-free dilution tubes.

Analytical procedure

3.212 This procedure is as follows:

- a. Dispense 0.1 mL of LAL reagent into each of six reaction tubes.
- b. Add 0.1 mL of sample to each of two of these tubes (duplicate determination).
- c. Add 0.1 mL of pyrogen free water to each of two of these tubes (negative controls).
- d. Add 0.1 mL of the working dilution of control standard endotoxin to each of two of these tubes (positive controls).
- e. Incubate for 1 h at 37°C.
- f. After 1 h, carefully remove each tube and invert slowly. If a firm clot has formed which does not move, the sample contains at least 0.125 EU mL⁻¹ endotoxin.

Sources of error

3.213 These are:

- a. The various commercially available LAL reagents differ in their tolerance of interfering substances, their inhibition endpoints and in their buffer capacity. There might also be some variation between lots within supplies from any one manufacturer. Revalidation is therefore necessary for each change of lot/manufacturer and for water supply system.
- b. Turbidity can be confused with initial stages of gelation. It is essential that all apparent gel formation is verified by demonstrating a stable gel on inversion through 180 degrees.

- c. Vibration during the incubation period can prevent stable gel formation. An unstirred water bath should be used since the vibration from the stirrer motor may be sufficient to interfere with the reaction.
- d. Different surfaces have different affinities for endotoxin. New borosilicate glass test tubes have a relatively high affinity for endotoxin in aqueous solution and may give rise to artificially low readings if used to make dilutions of endotoxin. Strict adherence to the reagent manufacturer's recommendations for choice and preparation of test equipment is necessary.
- e. Standard deviation/precision: for a given batch of LAL of calibrated sensitivity EU, a positive result indicates the presence of endotoxin within the range
- f. Limit of detection: a limit of detection of 0.03 EU mL⁻¹ is achievable using the gel formation method (0.01 EU mL⁻¹ and 0.001 EU mL⁻¹ are achievable for chromogenic and kinetic turbidimetric methods respectively).

Safety information

- 3.214 No hazardous component has been identified in commercially produced limulus amoebocyte lysate. However, since it contains proteinaceous material (including human serum albumin) the possibility of allergic reaction, especially after prolonged or repeated exposure, should be considered.
- 3.215 Other reagents, although sterile and pyrogen free, are not intended for administration to animals or humans.
- 3.216 The lyophilized endotoxin that is used to prepare reference standards has a threshold limit value of 5 EU (kg⁻¹ h⁻¹). Over exposure can result in fever, nausea and shock. It is important to avoid inhalation of the powder and parenteral injection of the reconstituted solution.
- 3.217 Emergency first-aid procedures: if endotoxin solution has been injected treat the patient as for endotoxin shock, if the lyophilized endotoxin has been inhaled, remove the patient to fresh air and treat as for endotoxin shock.

Results

- 3.218 The endotoxin concentration in final rinse water for items intended for invasive use, parenteral

administration or the preparation of parenterals should not exceed 0.25 EU mL⁻¹.

Total viable count

- 3.219 When the operating cycle of the washer-disinfector requires that the product is rinsed after the disinfection stage the rinse water should be free from microbial contamination, which could compromise the intended use of the load. A total viable count should be made on the final rinse water.
- 3.220 The following test will only detect the presence of mesophilic aerobic bacteria that do not have specialized nutritional requirements. If particular microorganisms are of concern, other recovery conditions (growth medium, incubation temperature, etc.) should be used as appropriate. The advice of the microbiologist should be sought.
- 3.221 The following test should be carried out as an installation or operational test and as a weekly periodic test for endoscope washer-disinfectors.
- 3.222 For other washer-disinfectors the test should be carried out when requested by the User as an installation test and repeated annually thereafter.
- 3.223 The following test should be carried out by the Microbiologist.

Apparatus

- 3.224 The following apparatus should be used:
- 250 mL sterile single-use containers;
 - sterile filter membranes (47 mm diameter, ≤0.45 µm pore size);
 - suction filtration apparatus;
 - incubator set at 35°C ± 2°C;
 - any general-purpose non-selective agar such as tryptone soya;

Note

The medium should have been demonstrated as capable of recovering an inoculum of 10–100 cells of *Pseudomonas aeruginosa*.

- 70% alcohol wipes.

Method: sample collection

- 3.225 Wipe the discharge surfaces of the sampling point thoroughly with 70% isopropanol and allow to evaporate to dryness.
- 3.226 Run off 200 mL through the sampling point and discard.
- 3.227 Take the sample downstream of any filter or other device or equipment intended to remove or control microbial contamination in the water supply.
- 3.228 Using aseptic handling techniques collect not less than 200 mL of sample in the sterile container and close the lid securely. Label the container with details of the sampling point and the time and date the sample was collected.
- 3.229 Transfer the sample to the laboratory for testing within 4 h; if this is not possible the sample should be stored at 2°C to 5°C for not more than 24 h before testing.

Method: testing

- 3.230 Filter a 100 mL aliquot of the sample through a 0.45 µm filter. Aseptically transfer the filter to the surface of a agar plate and incubate at 35°C ± 2°C for 48–72 h. Carry out the test in duplicate.
- 3.231 Examine the filters at the end of the incubation and record the number of colony forming units that are visible.

Results

- 3.232 For washer-disinfectors in which the product is rinsed after the disinfection stage there should be no recovery of microorganisms from the rinse water.
- 3.233 All other water services supplied to washer-disinfectors should have less than 100 cfu/100 mL of water (determined as the mean of the duplicate tests).

4 Operational management

Introduction

- 4.1 This chapter covers the maintenance and operation of the various types of washer-disinfectors used in healthcare facilities.
- 4.2 Terminology used in decontamination has long been inconsistent and this has often led to ambiguities. European and International Standards adopt a common set of definitions for terms relating to decontamination. Reference should be made to these documents for definitions.
- 4.3 The testing, maintenance and reporting procedures recommended in this section are based upon good practice in both the United Kingdom and the rest of Europe, as formalized in new European Standards designed to support the new EU Directives. They are designed to prevent the possibility of gross failure and serious incident.
- 4.4 Good staff morale is important. Anomalous behaviour that may foreshadow a malfunction of a washer-disinfector is often first noticed by an alert operator or other relatively junior employee. It is vital that staff feel free to report such observations promptly, and that appropriate remedial action is taken.
- 4.5 It is important that all washer-disinfectors are effective in achieving the performance required to produce a clean and disinfected product, and that they are safe in operation.
- 4.6 Failure to achieve the required standard of cleanliness might also impair the capability of the process to achieve disinfection/subsequent sterilization.
- 4.7 The cleanliness and microbial safety of all products processed in a washer-disinfector ultimately depends upon the care taken by the personnel responsible for its design, manufacture, installation, operation, test and maintenance.
- 4.8 Cleaning and disinfection might appear to be relatively simple processes, but considerable care is needed to consistently achieve satisfactory results.

- 4.9 Responsibility for assurance on these points rests variously with the AE(D), AP(D), CP(D), the Microbiologist, the Control of Infection Officer and the User.
- 4.10 A permit to work system should be utilized as given in the 'General' section of CFPP 01-01 Part B.

Equipment damage

- 4.11 The User should ensure that chemical additives used in the decontamination process are compatible with the materials of which the washer-disinfector is constructed and also with the items to be processed.
- 4.12 Most washer-disinfectors are made partly or wholly of stainless steel; the water supplied to the chamber, and the detergent and other chemical additives used, should have a low chloride content to minimize the risk of corrosion.
- 4.13 Lubricants for squeeze tubes on peristaltic pumps and other dispensing devices should be chosen and used in accordance with the manufacturer's instructions.
- 4.14 Care should be taken to ensure that the walls of ultrasonic tanks are not scratched as this can cause serious cavitation erosion.
- 4.15 Operators should be instructed never to drop or rest load items directly on the bottom of an ultrasonic tank.

Maintenance

Introduction

- 4.16 Decontamination is a set of processes. The efficacy of which cannot be verified retrospectively by inspection or testing of the product before use. For this reason decontamination processes (cleaning and disinfection and/or sterilization) have to be validated before use, the performance of the process routinely monitored, and the equipment maintained.

4.17 Means of ensuring that a washer-disinfector is fit for its intended purpose will include the validation and testing programme recommended in 'Validation and verification', and also the programme of planned maintenance recommended in this section.

4.18 The philosophy of maintenance and testing embodies two main principles to ensure that the required standards of performance and safety are met and maintained:

- all washer-disinfectors are subject to a carefully planned programme of tests to monitor their performance;
- all washer-disinfectors are subjected to a planned programme of preventive maintenance.

4.19 Expertise on the maintenance of washer-disinfectors is available at three levels: the CP(D), AP(D) and the AE(D). The roles and responsibilities of the AE(D) are defined in the 'Staffing roles and responsibilities' section of CFPP 01-01 Part A, and those of the AP(D) and CP(D) are defined in the 'Responsibilities' section of CFPP 01-01 Part B, which also describes the permit-to-work system (see links below).

4.20 Recommendations for testing of washer-disinfectors are given in [Chapter 2, 'Validation and verification'](#).

Planned maintenance programme

4.21 The planned maintenance (PM) programme should be designed according to the following principles:

- a. All parts of the reprocessor vital to correct functioning or safety should be tested at weekly intervals. This is interpreted as follows:
 - (i) there is no need to test components individually in cases where any malfunction will be revealed by the periodic tests recommended in [Chapter 2, 'Validation and verification'](#) for weekly or more frequent intervals;
 - (ii) where the correct functioning of important components is not necessarily verified by the periodic tests recommended for the washer-disinfector, those components should be individually tested each week and reference to testing them should be included in the schedules of maintenance tasks. This applies, for example, to door interlocks which may

only be required to perform their safety function when presented with an abnormal condition;

- b. The maintenance programme should include, at appropriate intervals, those tasks such as lubrication and occasional dismantling of particular components (such as pumps) necessitated by normal good practice, manufacturer's advice and experience. Apart from these tasks, the maintenance programme should concentrate on verifying the condition of the washer-disinfector and its components by means of testing and examination without dismantling. Parts which are working correctly should be left alone and not disturbed unnecessarily.
- c. Maintenance should be carried out under a quality system such as BS EN ISO 9001. Spares fitted to washer-disinfectors constructed under a quality system should be sourced from a similarly approved quality system.

Design of a PM programme

4.22 The PM programme recommended by the manufacturer should be supplied and should be used. The maintenance programme may be modified subsequently to take account of equipment use, equipment history and local conditions after a suitable period of operational experience.

4.23 The manufacturer's schedule may be modified if necessary but only after discussion with the AE(D), AP(D) and User.

4.24 Although the manufacturer might carry out certain inspection and maintenance procedures under the terms of its guarantee, these might not constitute a full PM programme. The User should therefore ensure that the complete PM programme is carried out by the CP(D) (who may be an employee of the manufacturer) during the guarantee period. The User should also implement any reasonable instructions given by the manufacturer during this period. Failure to carry out maintenance tasks and periodic tests could affect safety. It could also allow a contractor to place some, if not all of his liability on to the management. Where maintenance is carried out under lump sum term contract (see [Chapter 1, 'Design and pre-purchase considerations'](#)) such failure is tantamount to breach of contract and can give the contractor cause to terminate the contract if it so wishes.

- 4.25 A set of procedures should be developed for each model of washer-disinfector, each containing full instructions for a particular maintenance task.
 - 4.26 The frequency with which each task will need to be carried out will depend, in part, on the usage level for the washer-disinfector and on the quality of the water/steam supplied to the washer-disinfector. It might be necessary to adjust the programme so that work is carried out more frequently on machines that are heavily used/supplied with hard water.
 - 4.27 It is important that maintenance is planned so that the washer-disinfector is out of service as little as possible. Maintenance should, where practicable, be scheduled to immediately precede the periodic tests as specified in **Chapter 2, 'Validation and verification'**.
 - 4.28 Systematic records should be kept of all maintenance work undertaken both to demonstrate that the work has been carried out and also to facilitate periodic review of the PM programme.
 - 4.29 Maintenance and facilities management software packages (e.g. WIMS) may be used to maintain a full technical and financial history of the equipment.
- d. any changes required to any maintenance procedure;
 - e. any additional training required by maintenance personnel.
- 4.33 Proposed changes to the PM programme should be made in consultation with the manufacturer whenever possible.

Modifications

- 4.34 Occasionally, modifications to the washer-disinfector might be recommended by the manufacturer or by the UK Health Departments for reasons of efficacy and safety. The User should arrange for such modifications to be carried out within a reasonable period, normally coinciding with a scheduled maintenance session.

Routine housekeeping

- 4.35 Certain maintenance tasks may be carried out by the User, or by the Operator under the User's supervision, and should be recorded in the washer-disinfector log. Examples of such tasks include:
- a. checking that the rotating spray arms are free to rotate;
 - b. checking that nozzles are not blocked;
 - c. removal, cleaning and replacement of strainers and filters;
 - d. checking that the supply of chemical additives is sufficient for the day's use and replenishing if necessary;
 - e. cleaning the inside of the chamber;
 - f. cleaning the external surfaces of the washer-disinfector;
 - g. washing of loading side conveyors and trolleys;
 - h. for washer-disinfectors with a built-in water softener, checking the level of salt in the regeneration tank and replenishing if necessary.

Warranty period

- 4.30 After the purchase of a new washer-disinfector the manufacturer might carry out certain inspection and maintenance procedures under the terms of the warranty. This might not be a full PM programme. The User should ensure that the complete PM programme is carried out by the CP(D) during the warranty period.
- 4.31 The User should comply with any reasonable instructions from the manufacturer during the warranty period.

Review of PM programme

- 4.32 The PM programme should be reviewed at least annually to ensure that the equipment is being fully maintained but without any unnecessary maintenance activity. The review should aim to identify:
- a. the adequacy of maintenance records and compliance with the PM programme;
 - b. any emerging defects;
 - c. any changes required to the PM programme;

Pressure Systems Safety Regulations

- 4.36 Requirements of the Pressure Systems Safety Regulations should be met following advice from the CP(PS).

Features requiring special attention

Steam generators

- 4.37 Steam generators, provided as an integral part of some washer-disinfectors, are steam boilers and they should be subject to a written scheme of examination for pressure vessels.
- 4.38 Steam generators constructed from stainless steel will be subject to the risk of stress corrosion cracking. To minimize the risk the manufacturer's guidance on feedwater quality should be followed.

Leak tightness of doors

- 4.39 The door(s) of the washer-disinfector are intended to prevent the escape of fluids into the surrounding environment, that is to ensure freedom from aerosols that might be potentially infectious.
- 4.40 Damaged door seals are the major potential source of leaks and should receive careful attention as advised by the manufacturer.
- 4.41 The working life of door seals can be prolonged by regular cleaning.
- 4.42 Door seals should be renewed with spares provided, or approved, by the manufacturer at recommended intervals or when there is any sign of damage/deterioration.

Door interlocks

- 4.43 The interlocks on door(s) of the washer-disinfector are intended to:
- prevent the operator gaining access to the load during processing;
 - prevent both the loading and unloading doors being open at the same time on "pass-through" washer-disinfectors;
 - prevent the operator gaining direct access to a load that has not been satisfactorily processed.
- 4.44 Maintenance and inspection of door safety devices and door interlocking and chamber sealing systems should be carried out in accordance with the manufacturer's written instructions.
- 4.45 Security and settings of door safety switches and interlocks should be checked at least monthly. The setting should be within the limits specified by the manufacturer.

Chemical dosing systems

- 4.46 The correct amount of chemical additive should be administered at the right time in the operating cycle to ensure the correct functioning of a washer-disinfector.
- 4.47 The chemical additive dosing system should be subjected to regular (at least daily) inspection, maintenance and test. This should include:
- visual inspection of all piping to ensure freedom from leaks;
 - visual inspection/testing to ensure that neither the delivery or pick-up piping is blocked by coagulated or hardened chemical additive (many of the chemical additives used are a viscous suspension), followed by cleaning or replacing piping as necessary;
 - lubrication of the pinch tubing on peristaltic pumps in accordance with the manufacturer's instructions;
 - ensuring that there is sufficient additives(s) available and that it (they) are being dosed.

Water sprays and jets

- 4.48 The correct flow and distribution of water and aqueous solutions throughout the chamber and load are essential to the correct functioning of a washer-disinfector. The spray system should be checked on a daily basis as part of the routine housekeeping tasks carried out by the User.
- 4.49 In addition maintenance staff should also check the system at least weekly; this should include:
- checking that the rotating spray arms, both installed within the chamber and located on load carriers, are free to rotate;
 - checking that nozzles are not blocked; clean and/or replace if necessary;
 - checking for wear in bearings of rotating parts; replace any worn parts as necessary;
 - checking the mating of any necessary connection between the load carrier and the water supply in the chamber.

Ultrasonic transducers

- 4.50 Many ultrasonic cleaners are not fitted with means to provide continuous monitoring of performance. Transducers can fail or become detached from the ultrasonic tank without being noticed by the

Operator (other than by the deterioration in the cleaning performance).

- 4.51 Periodic functional testing of ultrasonic cleaners is defined in this CFPP.
- 4.52 The tank of the ultrasonic cleaner should be cleaned with a suitable neutral detergent and soft brush at least weekly.

Ventilation plant

- 4.53 Correct operation of ventilation plant is essential to ensure:
- the safe operation of washer-disinfectors that include a chemical disinfection stage;
 - the efficient operation of the drying stage (where this is included);
 - the maintenance of a comfortable working environment.
- 4.54 All ventilation systems associated with a washer-disinfector should be inspected, serviced and maintained at least every six months. Guidance on maintenance is given in Health Technical Memorandum 03-01 Part B.
- 4.55 Washer-disinfectors that include a chemical disinfection stage should have the associated ventilation system examined and tested annually.
- 4.56 Before undertaking maintenance work on the machine covering/fascia, or its associated ventilation system, it may need to be decontaminated and the advice of the designated safety officer should be sought. A permit-to-work system should be in operation.

Returning a washer-disinfector to service

- 4.57 Whenever any work has been carried out on a washer-disinfector, whether or not this was part of the PM programme, the user should be satisfied that it is fit for use. Following major repairs, overhauls, etc. which might affect the performance of the washer-disinfector, the User and AP(D) with assistance from the AE(D), should draw up a schedule of checks and tests to be carried out before the washer-disinfector is returned to service. This should include some or all of the recommissioning (yearly) tests specified in **Chapter 2, 'Validation and verification'**. See also guidance on the permit to work system given in the 'General' section of CFPP 01-01 Part B.

Troubleshooting

- 4.58 A failure to clean all the items processed in a load through a washer-disinfector is the most frequently observed fault. The most common causes of this type of failure, and thus those that should be considered first in any investigation, are:
- incorrect loading:
 - items that are not correctly located in an appropriate load carrier will not be subjected to the intended washing process;
 - overloaded baskets and load carriers allow some load items to shield others from spray jets, etc.;
 - hinged instruments that are not opened prior to washing will not be effectively cleaned;
 - blocked spray jets, spray arms that are not free to rotate, or a blocked strainer in the chamber base;
 - soiled instruments which have been stored for prolonged periods before decontamination: blood and protein will coagulate if stored for more than 8 h making this hard to remove;
 - soiled instruments subjected to heat treatment before decontamination: this will lead to coagulation of blood and protein making this hard to remove;
 - incorrect choice or quantity of detergent:
 - the detergent chosen should be compatible with the loads to be processed, the soil to be removed and the quality of water supplied;
 - malfunction of the dosing system may cause the wrong quantity of chemical additive to be used: too little will not provide the detergency required but too much may also impair cleaning by causing excessive foaming, etc;
 - Inappropriate water quality:
 - initial flush with water that is too hot will lead to coagulation of blood and protein making this hard to remove;
 - the hardness of water used during washing should be compatible with the detergent chosen;
 - hard water used in the final rinse will leave deposits on the surface of instruments.

- 4.59 The choice of detergent should be based upon a number of factors. These include:
- a. the quality of water available;
 - b. the nature of the load;
 - c. the nature of the soiling to be removed;
 - d. the nature of the washing process.
- 4.60 Advice should be sought from both the washer-disinfector and detergent manufacturers.
- 4.61 The quality of water used for the final rinse stage is important in ensuring freedom from scaling, process residuals etc. Water purified by deionisation, distillation or reverse osmosis should be used for the final rinse stage since this gives the lowest levels of process residuals.

Appendix A: Particular specification for washer-disinfectors used for processing surgical instruments

This appendix contains a purchasing specification for use when when planning to buy those washer-disinfectors used for processing surgical instruments.

[Click here for the full specification.](#)

References

- BS EN ISO 15883-1.
BS EN ISO 15883-2.
MDA DB 2006(04).
HSC 1999(178).
Health Building Note 13 – ‘Sterile services department’.
Medical Device Directive (93/42/EEC).
Medical Device Directive (2007/47/EEC).
Electromagnetic Compatibility Directive (2004/108/EEC).
Low-voltage Directive (2006/95/EC).
Pressure Equipment Directive (97/23/EEC).
Machinery Directive (2006/42/EEC).
Health Technical Memorandum 02-01 – ‘Medical gas pipeline systems’.
BS 853-1.
Water Supply (Water Fittings) Regulations 1999.
Manual Handling Operations Regulations 1992.
PD 5500.
BS EN 837.
ISO 554.
BS 3928.
Health Technical Memorandum 03-01 Part B.
Control of Substances Hazardous to Health Regulations 2002.
- BS EN 61010-1.
ISO/TS 15883-5.
BS EN ISO 4788:2005.
BS 5404-2.
BS EN ISO 14644-1.
BS EN ISO 3746.
BS EN 61000-4-3.
BS EN 55014-1.
BS EN 61000-6-3.
BS EN 61000-6-1.
Chemical Abstracts Service.
BS 1752, ISO 4793.
BS 1427.
BS 6068-2.37:1990, ISO 9297:1989.
British Pharmacopoeia.
ISO 8288.
BS 6068-2.2:1983, ISO 6332-1982.
BS 1704:1985, ISO 1770:1981.
BS EN ISO 6878:2004, BS 6068-2.28:2004.
BS 2690-104.
BS EN ISO 9001.
Pressure Systems Safety Regulations 2000.

Particular specification for washer-disinfectors used for processing surgical instruments

Section 1

Refer to Choice Framework for local Policies and Procedures (CFPP) 01-01 Parts A, B and D also.

Name of Trust	Details
Purchaser	
Hospital site	
Department	
Name of SSD manager and contact details	
Name of estates contact and details	
Authorising Engineer (Decontamination)	

The machine(s) are to be supplied under the Trust Contract Conditions or the NHS Supply Chain framework agreement.

Note: Site visit(s) are required by the supplier to ensure that the machine(s) will fit correctly and that no problems will be encountered during the delivery process. All engineering systems and services should be surveyed during the visit(s).

Standards relevant to this equipment:

BS EN ISO 15883-1:2009. Washer-disinfectors. General requirements, terms and definitions and tests.

BS EN ISO 15883-2:2009. Washer-disinfectors. Requirements and tests for washer-disinfectors employing thermal disinfection for surgical instruments, anaesthetic equipment, bowls, dishes, receivers, utensils, glassware, etc.

Standards relevant to decontamination management

BS EN ISO 13485:2003. Medical devices. Quality managements systems. Requirements for regulatory purposes.

Standards relevant to safety requirements for decontamination equipment

BS EN 61010-2-040:2005. Safety requirements for electrical equipment for measurement, control and laboratory use. Particular requirements for sterilizers and washer-disinfectors used to treat medical materials.

BS EN ISO 13849-2:2008. Safety machinery. Safety-related parts of control systems. Validation.

Standards relevant to medical devices

BS EN 1041:2008. Information supplied by the manufacturer of medical devices.

BS EN ISO 17664:2004. Sterilization of medical devices. Information to be provided by the manufacturer for the processing of resterilizable medical devices.

BS EN ISO 14971:2009. Application of risk management to medical devices.

1 Washer-disinfector selection details

Total number of machines requiredas below

Washer-disinfector type	Type A Standard cabinet Pass-through double interlocked doors	Type B Cabinet single door	Type C Multi-chamber design
Numbers of machines			
Chamber capacity (nominal)			
Special requirements			
Door movement preferred (upward closing, drop-down door etc.			

Surgical specialties of department (optional details)

Types of surgery	Daily throughput	Loading system	General comments

The supplier should discuss details and installation with the customer prior to purchase.

Type C multi-chamber design and sections required:

Chamber function				

2 Washer-disinfector cycles requirements

Cycle	Required (yes or no)	Options and comments
General instrument cycle		
Laryngeal mask airway (LMA) – anaesthetic cycle (no rinse aid)		
Lubricants		
Descale cycle		
High-risk cycle		
Drying only		
Door interlock system (passing carriages back from the inspection, assembly and packing (IAP) room)		
Other		

3 Washer-disinfector process chemicals – options

Chemistry used within operating cycles	Required (yes or no)	Options and preferred choice
Wash detergent		
Rinse aid (to aid drying performance)		
Lubricants		
Descale		
Other(s)		

4 Details of delivery/installation requirements

4.1 Any comments on interim storage requirements or installation of the delivered equipment prior to final installation.

4.2 It is the responsibility of the supplier to establish the site access, route and requirements of delivery of the equipment to the final installation site.

Comments

5 Delivery details of packing methods

Select:

Standard packing for basic weather protection - A

Good weather covering to protect machines under delivery - B

Dust-proof packing and wrapping for further storage needs - C

Dust-proof packing and timber casing - D

--

6 Removal and disposal of existing plant, equipment and services

Details:
Washer-disinfectors
Plant
Services

7 Drawings

7.1 Layout drawings should be submitted to the User prior to tender to view the details of the installation.

7.2 Any drawings such as engineering services supplied by the supplier or required by the User should be clearly agreed and defined during the tender process.

7.3 All service(s) and connections should be agreed by the supplier and User (or representative) during the tender process. These connections will then be clearly illustrated on the drawings as submitted with the tender.

8 Documentation

8.1 Machine manuals should be supplied with the washer-disinfectors on site delivery.

8.2 Pressure vessel certificates are to be supplied on machine delivery, if required.

9 Air supply (if required)

9.1 Compressed air may be required for process control or for drying.

9.2 It should be agreed at the tender how the air will be supplied to the washer-disinfector(s).

Select one or more:

- A Individual machine compressors
- B Common supply
- C Air compressors paired up per two machines
- D None supplied with tender
- E Spare compressors supplied
- F Other
- G N/A

Comments and details

10 Heating medium

10.1 The washer-disinfectors may be steam- or electrically heated in the various stages. Medium choice will affect the cycle time. The supplier should discuss the options available and services required with the User.

Stages	Steam (yes or no)	Electric (yes or no)
Washing		
Thermal disinfection		
Drying		
Or ALL stages		

11 Steam

11.1 A pressure-reducing valve should be fitted in the steam line from the mains supply to the machine to protect any direct-acting valves on the machines from damage and to ensure a steady and safer supply pressure for the process.

11.2 Steam may be supplied indirectly to heat the water/air to the relevant control temperatures during the cycles through suitably manufactured heat exchangers.

Further details

12 Electrical supply

12.1 It should be agreed at the tender stage what the electrical loading/demand is for the installation.

12.2 Discussions should be held with the relevant estates department officers and the suppliers to determine the supplies in general, and if single- or three-phase is available or required.

Further details

13 Water supply

13.1 An assessment should be carried out on the supply water used in all phases to the washer-disinfectors prior to the procurement process.

13.2 A decision on what further treatment is required can be assessed in conjunction with the User and AE(D).

13.3 The supplier should provide advice on the minimum supply pressure(s) required at each stage of the process(s).

Further details

14 Washer-disinfector monitoring

14.1 It is a requirement that cycle independent monitoring is fitted to each washer-disinfector. This should retain at least 20 days of cycle data and should monitor all the parameters/items listed below.

14.2 Monitoring could be a built-in supervisor system, electronic independent system or data recorder as agreed with facilities services and the User.

14.3 It is a requirement that the instrumentation is connected to the hospital IT server and system.

Details:

Cycle details to be independently monitored

15 Consumables

15.1 At the time of delivery of the washer-disinfector(s), consumables such as printer roles and cartridges SHOULD be supplied to the unit for a minimum of a three-month operating period of constant use.

15.2 Detergents, rinse aids and other solutions should be supplied by the manufacturer at the time of installation for at least three months' operation.

15.3 Discussions should be made by the supplier with the SSD manager or representative prior to tender to discuss the solution options available with machines.

15.4 Rinse-aid selection or requirements for final disinfection stage should be discussed.

Comments if required

15.5 Consumables are required to be supplied with the machines.

Detergents(s) required (make/type)

Rinse aids

Descalers

Details of general consumables required by the user

16 Chamber furniture required

Numbers and types of loading trolleys
Numbers and types of loading carriages (internal)

Further comments for loading equipment

Compatibility of matching existing equipment (hatches, chamber floor heights)

17 Testing and validation

If witnessing of factory testing is required, this will be identified in the tender documents.

IQ, OQ and validation testing will be carried out by the manufacturer.

The AE(D) will be monitoring and auditing all test results.

The supplier will consult with the User and AE(D) for any technical advice required.

Further comments/requirements

Details of any special loads

Testing and maintenance contracts are to be quoted by the manufacturer during the tender for the costs to be analysed by the User for machine care after the warranty period. This should include 3 quarterly periodic tests plus maintenance, and 1 annual revalidation plus maintenance.

18 Service response times and costs

Details and User response time(s) requirements

Breakdown advice time required

Site attendance time required

Spares availability in delivery to site

19 Fascia and panelling

Details of panelling required

Access door (if required)

Extra panelling requirements

13-amp twin sockets to be provided in the front panel as close to the machines as possible for the test equipment.

Position of sockets:

20 Training requirements

20.1 Staff training is required before the machine(s) can be put into service.

20.2 The training will include the monitoring system and logging requirements.

20.3 Factory testing can be arranged by prior agreement with the manufacturer.

20.4 Training is required for both engineering staff and users.

20.5 Full operational training for SSD staff will cover all staff who will be required to work on the machines.

20.6 Estates staff training will be required to cover:

- Cycle control
- Machine controls and operating procedures
- Door operations
- Loading equipment
- Monitoring equipment
- Fault finding
- Repair/dismantle main components
- Cycle operation via the valves and operation components
- Basic cycle programming and fault analysis
- Demonstration of the maintenance manual

20.7 Numbers of staff required for training:

Operational staff	Estates staff
Details of shift times	Details of shift times

21 Warranty

Warranty details should be quoted and agreed with the User and the date from which it will commence.

Costs in section 2

- The agreement should be clear before the purchase is made
- Extended warranty options can be quoted and discussed with the User to cover both maintenance and testing as required
- Number of visits per year
- Cost of each visit

22 Contract testing/maintenance

Contracts can be built into the tender with full consultation with the User.

- Quarterly testing contracts
- Breakdown call-outs
- Response times
- Maintenance contracts as required
- Availability of spares

Details to be given in section 2 by the supplier

23 Gauges fitted for the washer-disinfectors

Gauge	Yes	No
Steam supply pressure – mains		
Reduced steam supply pressure at machine		
Water supply pressure		
Air supply pressure		
Door-seal pressure if active seal fitted		
Pump pressure		
Mains water supply temperature		
Tank water temperature		
Hot water service temperature		
Others		

24 Ventilation requirements of washer-disinfectors

24.1 Individual machines are or can be ventilated independently or linked into a common extract. The supplier should clearly indicate at the tendering stage what details are required to be supplied by them or as fitted by the hospital unit.

24.2 Drawings and air duties should be supplied with the tender documents.

24.3 The supplier should inform the purchaser if fan(s) are required as part of the machine as supplied.

24.4 The supplier should inform the purchaser if a complete system is required and where it will be terminated under the supply contract, i.e. for others to design and extend the system to a safe extract position outside the building.

24.5 The system should be constructed to allow correct drainage of possible condensate forming in the duct to drain. It can be constructed or manufactured of non-corrosive materials, and heat recovery systems can be considered.

Additional notes/requirements to be included in the tender

Section 2

This section is a guide for the type of information and energy duties that is required by the User for a good and effective installation.

INFORMATION TO BE COMPLETED BY SUPPLIER

Details of microprocessor control system

The following information should be provided by the supplier:

Details of independent body where complete program and software are lodged.

Details of interface and file protocol requirements for transfer of data in the storage device to an external computer.

Details of diagnostic checks incorporated in the system.

Details (including cost) of the data storage device.

Maximum ambient temperature within the protective case°C

with an ambient temperature of°C

Interim storage requirements

Suppliers are required to advise of the storage conditions required if different from final installed location.

If interim storage is needed – state storage conditions required

Details:

Warranty Details

Length of standard/free warranty period offered:

Number of included service visits during warranty period:

Conditions of Warranty

Projected mean time between failures:

Guaranteed up-time:

Please state definition of up-time:

.....
.....
.....
.....

Please state remedy available to purchaser if guaranteed up-time is not achieved:

.....
.....
.....

Extended warranty options for service and maintenance

Please complete the following schedule with regard to a planned preventative maintenance and emergency call out contract to cover all items shown in the individual site schedule and to commence 12 / 24 / 36 * months after acceptance if required by the purchaser:

Number of service visits per annum

Duration of service visits.....hours per machine

Normal working hours are 0800–1800 unless otherwise stated:

All emergency call-outs included: *YES / NO

Price for emergency call-out during normal working hours, if not included: £..... per hour

All out of hours working included: *YES / NO

*DELETE AS NECESSARY

Details continued

Price for Saturday working £.....per hour

Price for Sunday working £.....per hour

Price for evening working £.....per hour

Price for bank holiday working £.....per hour

Response time to emergency call-outs (engineer on site) hours

Latest time on a working day to guarantee engineer on site same day

Base of engineer to service this site

How many other sites does he/she service

Number of engineers available to service this site

All spare parts included *YES / NO

Please list any parts that are not included that appear on the following lists:

Ten most used commodities by volume

Description Part No Delivery lead time Price (exc. VAT)

1.....

2.....

3.....

Most used commodities by value:

Description Part No Delivery lead time Price (exc. VAT)

- 1.....
- 2.....
- 3.....
- 4.....
- 5.....
- 6.....

Location of spare parts

Delivery lead time for spare parts

Is remote maintenance and diagnosis via modem available: *YES / NO

Price for supply and installation: £.....

Software upgrades (during warranty or maintenance contract period):

Safety/defect upgrades *Free of charge / At cost

New Applications *Free of charge / At cost

*DELETE AS NECESSARY

Annual maintenance contract costs including validation

Contract price for one year £..... exc. VAT

Five-year maintenance contract £..... exc. VAT

Annual maintenance contract costs excluding validation:

Contract price for one year £..... exc. VAT

Five-year maintenance contract £..... exc. VAT

Contract price for five years paid annually (including warranty).

The maintenance contract will be at this price with no price increases. These costs are not to form part of the total costs, but are to be provided as options for consideration.

Service requirements

The following information should be provided by the supplier for each type of machine supplied (based on a standard cycle being processed).

SERVICE REQUIREMENTS

machine number.
water flow rate.
water supply pressure.....
water consumption per cycle.....
drain flow rate.
drain size.
drain type.....
drain vent size and type.
compressed air flow rate.
compressed air supply pressure.
compressed air consumption per cycle.
electricity voltage.
electricity current.
electricity maximum power kW.....
air filter (air removal) expected life.....
test procedure(s) for filter integrity

If steam heating is used:

steam flow rate – average
steam flow rate – maximum
steam consumption per cycle
steam supply pressure
safety valve outlet size
condensate flow rate
sound power per washer-disinfector
total sound power all specified washer-disinfector(s)
process chemical cost per cycle:
other.....
cost per cycle
total energy cost per cycle (please specify cycle type)
.....

Overall washer-disinfector dimensions

The following information should be provided by the Supplier.

m/c no

internal chamber dimensions (H x W x L) mm.....

max floor area

height

max floor

loading

force kN/m²

max fascia opening

portage details

weight

Total cost of processing one DIN 1/1 tray basket (480 mm x 250 mm x 50 mm) containing reference load validation & verification on cycle suitable for surgical instruments (including energy and process chemical costs):

Energy cost basis:

Mains cold water £ /m³

Hot water £ /m³

Electricity per kWh

Steam £ per 1000 kg

Other.....

Overall cycle(s) time(s)

Note: For continuous process machines this should be expressed as the time from the first basket in to the last basket exit for a production load.

The following information should be provided by the supplier

Load (including number of baskets) Maximum Cycle Time

The following information should be provided by the Supplier.

Details (including weight and dimensions):

Heat emission

The following information should be provided by the Supplier.

Heat emission during normal operation at ambient temperature of 25°C:

to fascia – with door closed W

to plant area W

Contract completion

The following information should be provided by the Supplier:

time required from receipt of order in works weeks

time required for installation and pre-commissioning on site..... weeks

time required for commissioning on site weeks

Detailed cost breakdown

The following information should be provided by the Supplier:

Item Washer-disinfector Type Model
Name/No
No. of

Unit Total Price Chamber furniture

Summary of tender

The following information should be provided by the Supplier:

£

supply [nos of]..... washer-disinfector(s) ex works
delivery, offloading & positioning of washer-disinfector(s)
installation of washer-disinfector(s)
supply and installation of services
supply and installation of fascia panelling
site commissioning, i.e. installation checks and tests
test equipment, test loads and materials (if required)
12-month periodic testing and service including 3 off quarterly visits plus 1 off annual
staff training, consisting of days
supply chamber furniture /racks/loading trollies/ etc. type
costs of consumables
costs of cleaning solutions – detergents etc
independent monitoring equipment
removal of existing equipment
supply set(s) of recommended service spares
contingency – to be set by Purchaser
SUB-TOTAL
..... VAT @ %

Hospital:

Site:

Department:

TOTAL £.....

Comments

Date of tender