


### Spring Stringer

<b>Brand Name of Product</b>	Spring Locking Instrument Stringers
<b>Generic Name of Product</b>	Instrument stringers
<b>Product Code Number(s)</b>	USPRS-10, USPRS-12, USPRS-14, USPRS-16, USPRS-6, USPRS-8
<b>Purpose of Product</b>	This product is used to expand instruments to better expose hard to reach areas during cleaning and sterilization.
<b>Range of Applications for Product</b>	
<b>Key specifications of product</b>	<ul style="list-style-type: none"> <li>Manufactured of stainless steel</li> </ul>

Instructions for Using Product	
<b>Description of Use (s)</b>	Keeps instruments maximally exposed for hard to reach areas during cleaning and sterilization.
<b>Preparation for Use</b>	
<b>Diagrams (drawings, pictures):</b>	
<b>Steps for Use of Product</b>	<ol style="list-style-type: none"> <li>Unlock the stringer by gently pulling the smaller center piece of the spring out from inside the larger piece.</li> <li>Once the latch is unlocked, extend both sides of the latch so they are completely extended. String ring handled instruments in set assembly order onto stringer by placing one ring on each of the stringer arms. Continue until all instruments are loaded onto stringer.</li> <li>Lock latch by closing the latch into a secure locking position. The smaller end of the spring will insert inside the larger spring.</li> </ol>
<b>Interpretation of Results</b>	
<b>Contraindications of Test Results</b>	
<b>Documentation</b>	
<b>Special Warnings and Cautions</b>	
<b>Disposal</b>	

Reprocessing Instructions	
<b>Point of use:</b>	
<b>Preparation for decontamination:</b>	
<b>Disassembly Instructions:</b>	
<b>Cleaning – Manual:</b>	<ul style="list-style-type: none"> <li>Before initial use of the stringer, cleaning is required. <ul style="list-style-type: none"> <li>Manually wash stringer(s) with a neutral ph detergent and warm water.</li> </ul> </li> </ul>
<b>Cleaning – Automated:</b>	
<b>Disinfection:</b>	
<b>Drying:</b>	
<b>Maintenance, inspection, and testing:</b>	
<b>Reassembly Instructions:</b>	
<b>Packaging:</b>	
<b>Sterilization:</b>	Stringers are materially compatible with steam and EO sterilization.