

TOSI® - FlexiCheck

**Routine monitoring & validation of cleaning efficiency
of flexible endoscope washers and detergents in accordance with EN ISO 15883**



Product Code: 20240-25 (Test plate) & 04000-01 (FlexiCheck holder)

The FlexiCheck monitors the effectiveness of flexible endoscope washers and detergents.

The FlexiCheck can be used to assess either system and can simulate the conditions that a flexible endoscope undergoes during the washing process.

The stainless steel plate has two test soils, enabling the user to check for blood or mucus residue after washing.

By carrying out this test on a regular basis the user is able to confirm that the washers are functioning at the optimum performance level.

The test soil has been assessed as the closest correlated to human blood and mucus without any of the associated hazards, with the ability to test well into the microgram range.

The TOSI FlexiCheck offers an exacting test that will complement existing cleaning protocols while offering a quick and easy to read visual check of the machine's operational effectiveness.

All data results can be collated and kept as evidence of compliance.

