

## HemoCheck-S

## A routine test method for detection of blood on surfaces. Complying with EN ISO 15883







Product Code:67067-12

While instruments may look clean there is still a high risk of cross contamination from material invisible to the human eye.

The HemoCheck-S is a quick and easy way to detect residues post decontamination and pre sterilisation.

This offers assured cleanliness before the risk of residue becoming baked-on during any subsequent sterilisation process.

Using an enzymatic reaction (which is also used in the forensic science field) the HemoCheck has the capability to detect blood to a very low level of 0.1 µg.

By simply rubbing the area in question with the swab and then placing the swab into the test vile, any detectable colour change will indicate the presence of blood.

No colour change means all is well and the item is clean.

HemoCheck offers an exacting test that will complement existing cleaning protocols while offering a quick and easy to read visual check of the decontamination procedure.

The data can then be collated and kept as evidence of compliance.



Web: www.minervastericlean.co.uk Email: info@minervastericlean.co.uk Tel: +44 (0) 1249 820178