


Famos Integrators

Integrators are easy to use indicator cards or self adhesive labels which can be attached or inserted into records, pouches or containers to monitor sterilisation parameters up to 134 deg Celsius from 3 to 18 minutes.



These self-adhesive integrators (Type 4) comply with **ISO-11140 Type 4** requirements and are available for Steam and Dry Heat, with future on-going development of EO, FO and Plasma labels

Steam Integrators (Type 5) integrators comply with **EN 867 Type D and EN ISO 11140-1 Type 5** and provide a simple single colour change reference.



The 3 Phase Integrator (Type 6) comply with **EN ISO 11140-1 Type 1**. 3 Phase Integrators are specially designed to check the time and temperature combination in three different steps & provide a more comprehensive assessment of sterilisation effectiveness (failure, safe & optimum).

The Prion Integrator complies with **EN ISO 11140 Type 5** and are suitable where higher sterilisation temperatures are reached over a longer period (18 mins).



Product Codes: 15.100.262,263 (Type4)
15.100.135,136,137 (Type 5)
15.100.065, 077 (Type 6)
15.100.064,078 (Prion)