

TOSI® - LumCheck

A cannulated cleaning efficacy test



Product Code: 20100-25 (Test plate) & 02000-01 (Holder)

Ensuring the cleanliness and efficiency of cannulated instruments has always been a challenge.

Utilising the test soil and stainless steel plate technology as used in the TOSI, placing the test in a cannulated stainless steel holder replicates the conditions of a surgical instrument.

We are able to check the efficiency of rigid scope cleaners by monitoring the washing cycle inside the cannulated instrument itself.

The test soil has been assessed to be the closest in correlation to human blood without posing any of the associated risks and hazards.

The test has the ability to detect residue well into the microgram range.

The TOSI LumCheck offers an exacting test that will complement existing cleaning protocols while offering a quick and easy to read visual check of the machine's operational effectiveness.

All data results can then be collated and kept as evidence of compliance.