

## The Aluminium Spray Gun

**For use with tap water and pressurised air or nitrogen**



**Product Code: 6000**

The standard set spray gun is the ideal instrument for cleaning laboratory glassware, pipettes, syringes, catheters, cuvettes, endoscopes, cystoscopes & cannulae etc.

The gun provides excellent dosage ensuring a highly satisfying cleaning result.

Attachments (available separately) are connected by seating them on the safety cone which is located at the tip of the gun, with confident engineering to ensure they remain in place during high pressure usage.

The gun is incredibly easy to use & can be used with either tap water, pressurised air or nitrogen.

**(They are not suitable for use with demineralised water)**

The gun is available as a stand alone item, or supplied with fabric reinforced PVC or stainless steel hosing in 1.5 or 3m lengths, with 1/2 or 3/4 inch fittings with the option of including a selection of accessories (spray nozzles, brushes, wall and table mounts).

Spare part items are also available.

Web: [www.minervastericlean.co.uk](http://www.minervastericlean.co.uk)  
Email: [info@minervastericlean.co.uk](mailto:info@minervastericlean.co.uk)  
Tel: +44 (0) 1249 820178