


Pyromol-E

**For rapid and confident testing for protein & blood traces
within cannulated instruments**



**Product Code: 77070- 06A, 06B, 06C, 06D
(1.7, 2.8, 3.8 & 5mm swab tips)**

The Pyromol test from Pereg is a simple protein test which when used with appropriately sized swabs can reliably and confidently be used to test the lumens of cannulated instruments such as flexible endoscopes, for the presence of protein, high levels of which can be found in blood.

The Pyromol-E test works on the practice of the test swabs passed entirely through the lumens of the instrument in question, followed by the tip of the swab being cut & placed into the appropriate test vial solution.

The presence of blood or protein is indicated by the appearance of a blue dot on the swab tip.

The test provides results that are easily interpreted & provide such accuracy as to ensure complete confidence when releasing instruments back into service.

The protein free spun polyester swabs are available in four sizes of tip (1.7, 2.8, 3.8, and 5.0 mm) & are attached to a long slender and ultra flexible & robust nylon shaft (approximately 2 meters in length) offering maximum control, manoeuvrability and penetration.

The Pyromol - E tests are supplied as a pack of six test vials & six swabs.

With a detection level as low as 1µg this self contained and pre-prepared test requires no further equipment and delivers rapid, reliable and objective test result.



Web: www.minervastericlean.co.uk
Email: info@minervastericlean.co.uk
Tel: +44 (0) 1249 820178